Here are some potential Q&A for this product on Amazon:

Q: How do I change between the Rapid mode and Basal mode?

A: We recommend watching the instructional video provided with the product that will demonstrate how to properly change between modes. The video will provide a visual demonstration.

- 1、Turn off
- 2. Press and hold the button for 6 seconds until the screen shows "Thermometer Icon"
- 3  $\sim$  Release the button, and the press the button 1 time, it is the rapid mode when the "Thermometer Icon" disappears
- 4 \ If you want to switch back to the BBT mode, following the steps until the screen shows "Thermometer Icon"

Q: Does this digital thermometer need any batteries?

A: Yes, this digital basal body thermometer requires 1 Lithium Metal battery to operate. The battery is included with the product.

Q: What are the two modes on this thermometer and how do they work?

A: This thermometer features two modes:

- 1. Rapid Mode: This mode takes a quick oral temperature reading in 8-12 seconds. It is suitable for general fever detection or temperature monitoring for adults and children.
- 2. Basal Body Mode: This mode takes 60-120 seconds to detect subtle changes in body temperature to track ovulation. It is recommended for women to detect ovulation and increase the chances of fertility. This mode detects changes as small as  $0.05\,^\circ$  C for accurate ovulation prediction.

Q: How accurate is the temperature reading on this thermometer?

A: This digital thermometer provides a high accuracy of  $\pm$  0.05  $^{\circ}$  C or  $\pm$  0.10  $^{\circ}$  F. It has a precision of 1/100th degree which allows for accurate ovulation and fertility temperature tracking. The accuracy helps in precise temperature-based ovulation prediction and family planning.

Q: Can I use this thermometer for my kids and babies? Is it safe?

A: Yes, this thermometer can be safely used for infants, babies, kids and children. It has a soft beep indicator and kid-friendly design. The rapid mode quickly provides temperature readings in 8-12 seconds, suitable for children who cannot keep the thermometer in mouth for long. The thermometer can also recall up to 30 previous temperature readings to track changes. It has a hygienic, washable tip for added safety. Always supervise kids during use and follow the included instructions.