



Manufactured by KBD, Inc.

FIJI SUN TANNING LAMP OWNER'S MANUAL

For proper usage, please read your Owner's Manual. Always wear protective eyewear during lamp operation.

KBD, Inc.
2550 American Court
Crescent Springs, KY 41017, USA
www.sperti.com

Phone 859-331-0800
Fax 859-331-0802
Toll Free 1-800-544-3757
Email info@sperti.com

The information provided in this booklet is for educational purposes only. The information is not medical advice and is not a substitute for medical care or professional legal advice. A medical practitioner or other expert should be consulted prior to any significant change in diet, sun exposure, exercise or any other lifestyle change. There shall be neither liability nor responsibility should the information provided in this booklet be used in any manner other than for the purposes of education.

CAUTION

PROTECT YOUR EYES DURING FIJI SUN OPERATION; wear protective eyewear. To avoid injury to eyes, protective eyewear (KBD, Inc. Cat. No. 31050A or FIN-EYECANDY) must be worn by all persons in room during lamp operation. Read the instruction booklet before use.

DANGER

Ultraviolet Radiation, follow instructions, avoid overexposure. As with natural sunlight, overexposure can cause skin injury and allergic reactions. Repeated exposure may cause premature aging of the skin and skin cancer. Wear protective eyewear, failure to do so may result in severe burns or long-term injury to the eyes. Medications or cosmetics may increase/decrease your sensitivity to ultraviolet radiation. Consult physician before using any ultraviolet lamp if you are using medications or have a history of skin problems or believe yourself especially sensitive to sunlight.

OPERATING INSTRUCTIONS FOR YOUR KBD-FIJI SUN LAMP

CAUTION: PROTECT YOUR EYES - WEAR PROTECTIVE EYEWEAR

To avoid injury to eyes, protective eyewear must be worn by all persons in the room during lamp operation.

- Put on protective eyewear
- Place unit 18" (inches) from skin
- Plug unit into 120 volt - 244 volt (50Hz/60Hz) power source
- Set timer to 15 minutes
- Turn power switch to on position
- Follow recommended exposure / maintenance schedule
- Unplug lamp for storage and service

DANGER: Ultraviolet Radiation, so follow instructions and avoid overexposure. As with natural sunlight, overexposure can cause skin injury and allergic reactions. Repeated exposure may cause premature aging of the skin and skin cancer. Wear protective eyewear, failure to do so may result in severe burns or long-term injury to the eyes. Medications or cosmetics may increase/decrease your sensitivity to ultraviolet radiation. Consult physician before using any Ultraviolet lamp if you are using medications or have a history of skin problems or believe yourself especially sensitive to sunlight.

RECOMMENDED EXPOSURE SCHEDULE

Different skin types require different exposure times and distances. The figures given for time and distances are safe for average skin, Type III. Do not exceed 3 exposures per week during tan build up. Initial exposure should not exceed 3 minutes for Type II skin or 4 minutes for type III, IV, or V skin at a distance of 18" per session. Exposure at distance less than 18" is not recommended. Distance is determined from metal guard, not light bulb.

Skin type I individuals, defined as one who "always burns easily: never tans (sensitive)" are advised not to tan. Less photosensitive skin types, while able to tolerate higher exposure doses without burning, also, by definition tan more easily and should not require greater exposure to achieve or maintain a tan. Consequently all other skin types (II-IV) use the same exposure schedule.

NOTE: Reddening or sunburn is usually not visible until 4-8 hours after exposure.

	SKIN TYPE	WEEK 1 3 EXPOSURES	WEEK 2 3 EXPOSURES	WEEK 3 3 EXPOSURES	MAINTENANCE 2 TIMES PER WEEK
Sensitive skin, very fair, red hair, dense freckles.	I	DO NOT USE THIS APPLIANCE			
Burns easily, Tans Minimally or lightly.	II	3 minutes	6 minutes	8 minutes	15 minutes
Often referred to as "average" complexion. Burns moderately, tans about average.	III	4 minutes	6 minutes	8 minutes	15 minutes
Burns minimally, tans easily and above average with each exposure.	IV	4 minutes	8 minutes	10 minutes	15 minutes
Rarely burns, tans easily and substantially.	V	4 minutes	8 minutes	10 minutes	15 minutes

0930
Rev. 1/2011



KBD, INC. HISTOLOGY

SPERTI SUNLAMP (manufactured/owned by KBD, Inc.), has been manufacturing tanning lamps since the 1930's. The original Sperti Sunlamp consisted of a filtered intermediate pressure mercury lamp, approved by the AMA in 1940 as a sunlamp. Today's KBD FIJI Sun Lamp is an updated version utilizing fluorescent technology, specifically formulated to produce ultraviolet energy in the specific SKIN TANNING range and eliminating any harmful UVC energy. In order to comply with modern safety guidance, the fifteen (15) minute timer provides an exposure, at 18" (inches) from the user's skin, of 1 (one) Standard Erythemal Dose (SED). The SED is an internationally standardized measure of UV radiation equivalency.

TUBE REPLACEMENT INSTRUCTIONS

(Recommended Annually)

Tools required: Phillips Head Screwdriver

Step 1 - Remove table stand. Position unit on a flat surface with wire guard facing ceiling.

Step 2 - Remove wire guard by compressing either side away from the metal housing.

Step 3 - Top of unit is the end WITHOUT the timer/switch. Remove the top front screw from both sides of the unit. Remove the only screw on the top of the unit.

Step 4 - Remove the top panel.

Step 5 - Rotate each bulb 1/4 turn either direction. Gently pull top of tube toward you. Lift out of base connector. Remove all tubes in this manner.

Step 6 - Align bottom pins of replacement tube front to back and insert into bottom lamp holder. Seat top set of pins by pushing into top lamp holder.

Step 7 - Rotate each tube 1/4 turn making sure the label/etching is facing you. You will notice the "lighter" part of the tube is in the back.

Step 8 - Replace top panel and wire guard.

After all tubes are replaced, your FIJI SUN lamp can be restarted.

FIJISUN & D/UV-4F VITAMIN D UNITS

SUBJECT: Fluorescent tubes do not ignite.

During the delivery process issues may occur preventing tube ignition.

PLEASE:

1. Check to be sure the tubes have not become "dislodged". If so, please reinstall the tubes.
2. Check plug for a snug fit into the appliance—thoroughly inserted into the unit.
3. Confirm electric outlet is energized.
4. Set timer to the maximum time, click on/off switch to the on position. Wait 3 seconds and click the switch to the off position. Repeat this procedure at least 8 to 10 times. Cycling the tubes with the on/off process allows any debris on the tube filaments to be cleared. This lighting process will not be required for future use.

Customer Service: www.sperti.com

[Sperti replacement bulbs are available at www.aviva.ca](http://www.aviva.ca)

AViVA
Natural Health Solutions

1224 St. James St. - Winnipeg MB
(204) 947-6789 or (866) 947-6789

www.aviva.ca

LIMITED WARRANTY KBD FIJI SUN TANNING LAMP
PROOF OF PURCHASE REQUIRED FOR ALL WARRANTY WORK

WHAT IS COVERED _____

All parts defective in materials or quality of workmanship.

FOR HOW LONG _____

Six (6) months from date of purchase.

WHAT WE WILL DO _____

We will decide either to replace or repair free of charge any defective part.

WHAT WE WILL NOT DO _____

Pay shipping or transportation costs.

WHAT YOU MUST DO _____

The warranty registration card must be completed within ten (10) days of the date of purchase. The completion and return of the warranty registration card to KBD, Inc.; 2550 American Court; Crescent Springs, KY 41017 within ten (10) days of purchase of this product is a condition of warranty coverage.

HOW TO GET WARRANTY SERVICES _____

For service, defective unit or defective part must be returned to the Factory Service Department: KBD, Inc., 2550 American Court, Crescent Springs, KY 41017.

OTHER CONDITIONS _____

This warranty covers the first purchaser only. This warranty does not cover commercial use. This warranty shall only apply to FIJI SUN Lamps that have not been tampered with or altered in any way or show evidence of misuse, neglect, or accident.

STATE LAW RIGHTS _____

Implied warranties are limited to the duration of six (6) months. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This warranty gives you specific legal rights and you may also have other rights which vary from state to state.

ABUSE OR MISUSE WILL VOID WARRANTY