

# seca robusta 813

The standing,  
taking everything lightly.



- Solid steel frame. For persons up to 440 lbs.
- Non-slip rubber base. For safe standing during measuring.
- Big LCD figures. For easy reading.
- Step-off function. For weighing without switching on scales.
- More than 165 years of experience. For maximum precision.
- From the global market leader for medical weighing and measuring.

## seca robusta 813:

The standing, making everything look a breeze.

**T**ake a look at the extra-wide weighing platform. Even very big people can stand here safely and securely. People of larger girth too can be accommodated. Our intention is to lighten the task for weighty people.



A scale of standing, making the sturdy look beautiful.

More and more people weigh more and more. Now there's a stable scale for them. Thanks to its solid steel frame, it can withstand 440 lbs and makes light of heavy use over many years. It's also no burden to step onto the scale, since its construction is extremely flat. That's a big advantage for people whose size makes it difficult for them to get about.

The standing, making the robust look elegant.

All rubber is not created equal. For this weighing platform seca has chosen a high-quality two-component nubby material which is soft like rubber and pleasant to the touch. But don't be fooled. This material is also tough. It fulfils two requirements at the same time: pleasant weighing and an easy-to-clean surface that stands up to all commercially available cleaning and disinfecting products.

The standing, paying attention to details.

Weight and girth are no problem on the seca robusta 813. You don't even have to bend over to switch it on, thanks to the unique Step-off function. An additional problem with many scales is the display. Not so here, because thanks to the 1.1 inch-high LCD figures, it is extremely easy to read the measurements. You see, a scale from seca makes life easier by paying attention to details.

*High-quality, two-component nubby material feels as soft as rubber.*



*Big LCD figures. For easy reading.*

## seca robusta 813

### Technical data

- Capacity: 440 lbs / 200 kg
- Division: 0.2 lbs / 100 g
- Dimensions scale (WxHxD): 17 x 1.8 x 14.7 inch
- Weight: 5.7 lbs / 2,6 kg
- Functions: Step-off, Auto-HOLD, automatic switch-off.
- Parts to be delivered: 4 AAA batteries, 4 extension feet for carpets