

Reynard Cleaning & disinfection range.

The ultimate clean team.

Stay one wipe ahead of germs.

Developed with convenience in mind, Reynard cleaning & disinfection range provides fast and effective solutions that save time and money. Suitable for a wide variety of cleaning purposes including surfaces and equipment. We also offer chemical-free cleaning as well as 100% biodegradable wipes.



Reynard is reassurance.

Cleaning & disinfection | Personal care | Dry wipes
Clinical skin cleansing | Hand hygiene | PPE



Premier detergent & disinfectant wipes

Our winning formula.

A non-ionic surfactant with non-streaking and superior wetting properties. Containing three active biocidal agents for maximum effectiveness and adaptability. Optimised for use in both clinical and non-clinical settings on surfaces and equipment. This single step solution kills 99.999% of pathogens and viruses after just 30 seconds. Tested and ARTG listed for use against SARS and COVID.

The all-in-one wipe that has the best in class performance for cleaning and disinfection.

Testing Data.

Effective against a wide range of pathogenic micro-organisms including fungi and bacteria, viruses such as SARS-CoV-1&2 and Mers-Cov.

Pathogen	EN Test	Time
Bacteria		
Kills at least 99.999% of pathogens. Effective from 60 seconds.		
Enterococcus Hiera	16615	1 min
Staphylococcus Aureus	16615	5 min
Fungi		
Aspergillus Niger	13624	5 min
Candida Albicans	13624	1 min
Candida Albicans	16615	5 min
Candida Auris	13624	1 min
Viruses		
Effective against enveloped viruses		
SARS-CoV-1&2 and Mers	14476	30 sec

Also tested to European standards EN1276, EN13727, EN14348, EN13704, EN1650, EN14675, EN14476, EN14651, EN1656M EN16615.



The cloth that is a cut above.

Made from a newly developed, soft apertured polyester substrate, that provides superb formula retention and extreme traction. This helps pick up dirt, bacteria, and pathogens where other, thinner cloths simply push it around the surface. Also available with a completely biodegradable, 100% bamboo cloth for an environmentally friendly alternative.

While germs aren't safe, your equipment is.

Developed with key clinical equipment providers to ensure an unparalleled range of device and material compatibility. Our formula was developed to ensure the wipes will effectively clean hard surfaces and equipment without causing harm.


BIODEGRADABLE




Product code	RHS210	RHS216	RHS218
Wipes per pack	50	50	100
Packs per carton	12	12	6
Wipe size	33 X 22cm	33 x 22cm	33 x 22cm



Product code	RHS200	RHS219	RHS205
Wipes per pack	200	200	280
Packs per carton	6	12	2
Wipe size	28 x 22cm	22 x 20cm	28 x 25cm



✓ Reynard Premier Disinfectant Wipe is a **28%** thicker cloth



Note: the above comparison is with the RHS216, RHS218 and RHS210 cloth. Our RHS219, RHS200 and RHS205 cloth is a thinner thermal bond cloth with similar thickness to the competitor cloth but still slightly larger.

Reynard is reassurance.



Neutral detergent wipes

Our expert formula works to lift soil and effectively clean surfaces without drying, spotting, or leaving residue. Created with the user in mind, this non-tainting, low-foaming formula is tough on germs but gentle on hands.

With unparalleled equipment compatibility, Reynard Neutral detergent wipes are approved for use by multiple local & global medical equipment manufacturers. Safe to use on most hard surfaces including mattress covers, tables & other hospital equipment as well as in the home. Our formula is alcohol and fragrance free and can be used on vinyl and leather.



Product code	RHS201	RHS203	RHS206
Wipes per pack	50	160	280
Packs per carton	12	12	2
Wipe size	33 X 20cm	22 x 20cm	28x 25cm

IPA Wipes

A hospital grade disinfection wipe impregnated with 70% isopropyl alcohol that will kill 99.9% of germs. An ideal solution for infection control or general cleaning, Reynard IPA wipes quickly and efficiently disinfect hard surfaces and equipment.



Product code	RHS402	RHS409
Wipes per pack	160	75
Packs per carton	12	12
Wipe size	33 x 20cm	42 x 14cm

Reynard is reassurance.



Absorb Cleaning Cloths


Initially designed for use in Central Sterilising Units and Operating Theatre environments, these versatile high-quality wipes are extremely durable and suitable for a wide variety of cleaning purposes.

The cloth is made from a unique non-woven, low lint material that is extremely strong and absorbent, making Reynard Absorb Cleaning Cloths ideal for general cleaning requirements within a hospital or healthcare environment.

		
Product code	RHS312	RHS313
Wipes per pack	100	50
Packs per carton	12	12
Wipe size	35 x 35cm	60 x 35cm

Extra strong General Purpose Wipes


The multitasker: highly absorbent and low-linting, they can be used wet or dry and are ideal for surface cleaning. Reynard Extra strong general purpose wipes are made from a strong, apertured, thermal bond polyester and viscose blend cloth ensuring they won't fall apart even when wet.

	
Product code	RHS316
Wipes per pack	100
Packs per carton	12
Wipe size	33 x 29cm

Disposable Microfibre range

What is a microfibre cloth?

Our cloths are a microfibre blend of polyimides, each fibre is finer than a single strand of silk. The material is split several times to make a finer thread which is then woven together. This allows the cloth to remove dirt more efficiently from surfaces due to the higher number fibres per cm² of cloth.



Due to their thinness and density, microfibre cloths can remove dirt using the mechanical force created by rubbing the cloth across a surface. The dirt particles are absorbed deep into the spaces created between the fibres. This allows the cloth to effectively Hoover up dirt, drawing it into the cloth rather than pushing dirt around the surface. Providing a superior clean without the use of any chemicals.

Benefits of using Microfibre

Cleaning performance

- Bacterial removal capability: tested and proven to remove 99.9% of Clostridium difficile spores and Pseudomonas aeruginosa spores
- Ideal for both wet cleaning and dry dusting
- Avoids streaking

Health and safety

- Chemical-free cleaning avoids pathogens developing a resistance
- Limits staff exposure to chemicals
- Less potential for slips on wet floors and other injuries
- No heavy water buckets to carry

Environment

- Reduces water usage
- Reduces plastic consumption
- Removes the need for laundering

Cost-effective

- Cheap enough to use once but durable enough to use again and again
- Disposable cloths remove costs associated with laundering cloths

Designed to deliver a high quality, chemical-free clean and reduce the potential for cross contamination.


Microfibre technology is proven to eliminate 99% of all dirt, dust, and pathogens. Our range can be used with or without steam, chemicals and other cleansers

Disposable microfibre cloths

		
Product Code	RHS315	RHS315-20
Wipes per pack	50	20
Packs per carton	12	28
Wipe size	40cm x 32cm	40cm x 32cm

Disposable microfibre mop heads

Reynard disposable microfibre mop heads are a fast and easy way to clean your floors with no chemicals required.

	
Product Code	RHS318
Mop heads per bag	25
Bags per carton	4
Mop head size	45cm x 13.5cm



Reynard is reassurance.



Biodegradable Multi-purpose wipes

Reynard Multi-purpose wipes are strong, large and extremely absorbent. The cloth is made with 100% viscous, making it completely biodegradable & environmentally friendly. Used wet or dry, these versatile high-quality wipes are extremely durable and suitable for a wide variety of cleaning purposes.

Our wipes are low linting and non-streaking, making them ideal for use on cleaning glass or other surfaces that require no residue to be left behind.



Product code	RHS601
Wipes per pack	50
Packs per carton	16
Wipe size	40 x 38cm



Australia
44 Hampden Road
Artarmon, NSW, 2064
enquiries@reynardhealth.com.au
1300 667 199

New Zealand
PO Box 8470
Havelock North
sales@reynardhealth.com
0508 494 737

United Kingdom
13 White Horse Business Park,
Stanford in the vale, SN7 8NY
uk@reynardhealth.com
+44 (0)845 263 7335

International
PO Box 8470
Havelock North
sales@reynardhealth.com
0508 494 737