

Reynard™ IV access and clinical skin cleansing. Our extensive range means you're prepped for anything.

The perfect solution for every
clinical skin cleansing situation.

Our broad range includes innovative packaging with 'easy-peel' and 'tear open' options. This, combined with unsurpassed value, ensure our prep pads and swabsticks are used daily throughout Australia. Providing every healthcare system provider in the country with the highest quality, TGA registered products.



Reynard™ is reassurance.



Reynard™ IV access & clinical skin cleansing.

The organisms may be micro,
but the impact can be big.

Surgical site infections (SSIs) are a major challenge for patients and health sector resources. Effective cleansing and disinfection using appropriate Reynard™ prep pads and swabsticks can help minimise the risk of surgical site infections and improve overall patient outcomes.

- SSIs negatively impact patient physical and mental health.
- Increased patient morbidity, mortality, and loss of earning during recovery are some of the indirect costs associated with infection.
- Intangible costs may also be incurred by the patient, such as pain and anxiety. In addition, patients may experience delayed wound healing and more susceptible to secondary complications, such a bacteraemia.
- Distress may also be caused to the patient and family member if the patient is absent from home and work for a prolonged period.
- Accordingly, prolonged hospitalisation and increased morbidity as a result of developing an SSI have been shown to negatively impact patient health-related quality of life.
- Increase in healthcare costs per patient and additional stress on healthcare systems and staff.

The perfect solution for every
problem.

Reynard™ provides the broadest range of prep pads and swabsticks in a variety of formulations and configurations to suit your specific clinical need and appropriate application.

- 0.5% chlorhexidine digluconate (CHG) and 70% isopropyl alcohol (IPA)
- 2% chlorhexidine digluconate (CHG) and 70% isopropyl alcohol (IPA)
- 70% isopropyl alcohol (IPA) – sterile
- 0.5% chlorhexidine digluconate (CHG)
- 2% chlorhexidine digluconate (CHG)



Reynard™ Health Supplies.

Reynard™ Health Supplies was established in New Zealand in 1999 and today is the number one supplier of surface cleaning & disinfection wipes, specialist patient care wipes and clinical skin cleansing products to hospitals in New Zealand and Australia.

Reynard™ is owned by a dedicated group of Australian and New Zealand shareholders who are all actively and directly involved in the day-to-day operation and development of the business globally.

Used extensively in the private sector and every public healthcare system in Australia.

Reynard™ is reassurance.

At Reynard™, we strive to provide products and support services that provide you with the reassurance you need to be able to focus on optimum patient outcomes.

This reassurance stems from:

- Quality & Range - a wide range of quality products to suit every need
- Cost Savings - cost effective and pricing stability provide real value
- Reliability - a reliable, proven supply chain that has demonstrated reliability and consistency throughout the global pandemic
- Support - Dedicated technical, educational and compliance support services
- Research & Development – focus drives continuous product development and a pipeline of innovative solutions
- Flexibility - the capability and desire to design new products and supply chain options



Reynard™ is reassurance.



Our extensive range means you're prepped for anything.

Size and type of product					
Formula	Prep pads				
	Substrate only				
	6 x 3cm	6 x 6cm	8.5 x 4.5cm	20 x 16cm	Easy Peel
70% IPA only	RHS451G	RHS417	RHS438		
2% CHG + 70% IPA	RHS461	RHS401	RHS437	RHS422	RHS422EP
0.5% CHG + 70% IPA					
2% CHG only			RHS421		
0.5% CHG only		RHS479			

Size and type of product								
Formula	Swabsticks							
	Cotton head				Foam Head			
	Regular	Triple Pack	Easy Peel	Tint	Regular	Maxi	Reg Tint	Maxi Tint
70% IPA only								
2% CHG + 70% IPA	RHS408	RHS408T	RHS408EP	RHS408TR	RHS473	RHS474		RHS472
0.5% CHG + 70% IPA					RHS470		RHS470TR	RHS471
2% CHG only					RHS475			
0.5% CHG only					RHS477	RHS478		

Skin cleansing considerations.

Product formulation	Recommended product application
0.5% CHG + 70% IPA	General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potentially allergies, skin reactivity and risk.
2% CHG + 70% IPA	General use where sensitivity etc. are not a concern. This solution provides broad capacity for most patients.
70% IPA Sterile	General use, also suitable for paediatrics and areas where use of CHG is not preferred due to potentially allergies, skin reactivity and risk
0.5% CHG only	General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potentially allergies, skin reactivity and risk. Alcohol (IPA) free, appropriate where use and application of IPA is not suitable.
2% CHG only	General use where CHG concentration sensitivity is not a concern. Alcohol (IPA) free, appropriate where use and application of IPA is not suitable.

All Reynard™ clinical skin cleansing products are suitable for a wide range of skin cleansing requirements including but not limited to:

- Blood culture collections
- Central venous catheters insertion
- General skin preparation, pre-procedure, injection, or incision
- Insertion of percutaneous devices
- Percutaneous device site maintenance
- Peripheral intravascular catheters insertion
- Site maintenance of central lines

Note: Only use on intact skin before any invasive procedures

Directions	Product should be applied with gentle friction using a repeated back-and-forth motion for manufacturer-recommended application time.
Dry time	Allow product to air dry completely; dry time can be affected by volume of product, body site, presence or absence of hair, humidity etc.

Reynard™ prep pad and Swabstick product application

Instructions for use:

1. Check allergy adverse reaction status with your patient.
2. Thoroughly clean and treat the chosen area with the prep pad using a repeated back-and-forth scrubbing motion for 20 to 30 seconds.
3. Use the necessary number of prep pads dependent on the clinical procedure.
4. Allow the skin area to fully air dry to be most effective against bacteria growth. Allow approximately 30 seconds for a dry site and up to 2 minutes for a moist site.
5. Discard of prep pad after single use, do not flush.
6. Dispose of prep pad in adherence with your facility infection control guidelines.

- Is this product flammable?**
- Refer to manufacturer's instructions for use, and any associated hazard warnings
 - Follow AORN guidelines for perioperative practice

Reynard™ is reassurance.



Prep pads

CHLORHEXIDINE & ALCOHOL



Product code	RHS401	RHS422	RHS422EP	RHS437
Swabs per box	200	100	50	100
Boxes per carton	25	10	20	10
Wipe size	6 x 6cm	20 x 16cm	20 x 16cm	8.5 x 4.5cm
Solution	70% IPA & 2% CHG	70% IPA & 2% CHG	70% IPA & 2% CHG	70% IPA & 2% CHG
Volume	0.7ml	4.5ml	5ml	1.6ml
GSM	45gsm	40gsm	40gsm	100gsm

CHLORHEXIDINE



Product code	RHS421	RHS479	RHS461
Swabs per box	100	200	200
Boxes per carton	10	25	25
Wipe size	8.5 x 4.5cm	6 x 6cm	6 x 3cm
Solution	2% CHG	0.5% CHG	70% IPA & 2% CHG
Volume	0.8ml	0.7ml	0.5ml
GSM	45gsm	45gsm	45gsm

ALCOHOL




Product code	RHS417	RHS438	RHS451G
Swabs per box	200	100	200
Boxes per carton	25	10	25
Wipe size	6 x 6cm	8.5 x 4.5cm	6 x 3cm
Solution	70% IPA	70% IPA	70% IPA
Volume	0.7ml	1.6ml	0.5ml
GSM	45gsm	100gsm	45gsm

Reynard™ is reassurance.



Skin surface area coverage for Reynard™ Prep pads*

22cm



22cm

Size	Product code	Coverage area
Small	RHS451G, RHS461	42cm ²
Medium	RHS401, RHS417, RHS479	100cm ²
Large	RHS421, RHS437, RHS438	225cm ²
X-Large	RHS422, RHS422EP	484cm ²

*Note: coverage area may vary in accordance with variations in atmospheric conditions and condition of the skin. Tests we conducted on hairless male skin at 21 degrees Celsius. The skin area was wetted sufficiently to ensure it was wet for a minimum 30 seconds before allowing to dry. In each test enough liquid was left in the product to allow up to a further 20% area to be covered. Diagram is not to scale.

Swabsticks

COTTON HEAD





Product code	RHS408	RHS408T	RHS408TR	RHS408EP
Sticks per sachet	1	3	1	1
Sachets per box	50	25	50	50
Boxes per carton	10	10	10	10
Stick type	Cotton	Cotton	Cotton	Cotton
Solution	70% IPA & 2% CHG	70% IPA & 2% CHG	70% IPA & 2% CHG	70% IPA & 2% CHG
Tint	No	No	Yes	No
Volume	1.6ml	4.5ml	1.6ml	1.6ml

FOAM HEAD




Product code	RHS470	RHS470TR	RHS473	RHS475	RHS477
Sticks per box	25	25	25	25	25
Boxes per carton	10	10	10	10	10
Stick type	Foam	Foam	Foam	Foam	Foam
Solution	70% IPA & 0.5% CHG	70% IPA & 0.5% CHG	70% IPA & 2% CHG	2% CHG	0.5% CHG
Tint	No	Yes	No	No	No
volume	2.5ml	2.5ml	2.5ml	2.5ml	2.5ml

LARGE FOAM HEAD



Product code	RHS471	RHS472	RHS474	RHS478
Stick per box	25	25	25	25
Boxes per carton	10	10	10	10
Stick type	Large Foam	Large Foam	Large Foam	Large Foam
Solution	70% IPA & 0.5% CHG	70% IPA & 2% CHG	70% IPA & 2% CHG	0.5% CHG
Tint	Yes	Yes	No	No
Volume	5ml	5ml	5ml	5ml

Skin surface area coverage for Reynard™ Swabsticks*



22cm

Size	Product code	Coverage area
Regular cotton swabsticks	RHS408, RHS408EP, RHS408T, RHS408TR	72cm ²
Regular foam swabsticks	RHS470, RHS470TR, RHS473, RHS475, RHS477	144cm ²
Large foam swabsticks	RHS471, RHS472, RHS474, RHS478	484cm ²

*Note: coverage area may vary in accordance with variations in atmospheric conditions and condition of the skin. Tests were conducted on hairless male skin at 21 degrees Celsius. The skin area was wetted sufficiently to ensure it was wet for a minimum 30 seconds before allowing to dry. In each test enough liquid was left in the product to allow up to a further 20% area to be covered. Diagram is not to scale.

Reynard™ focus on killing germs so you can focus on saving lives.

Germs.

They can proliferate in the cleanest of places. A constant bugbear, literally, to be eradicated, first time, every time. Reynard™ constantly adopt new technology, new research to develop new products in all our ranges. Working to protect patients, and everyone and everything they come into contact with. An ever-growing Reynard™ product range means you can choose a product that meets every operational need, whether cleaning, disinfecting, skin cleansing or patient care. Many applications developed in consultation with health professionals, in response to urgent clinical needs. So you can concentrate on doing what you do best. Caring for patients. Visit www.reynardhealth.com.au

Reynard is reassurance.



More from Reynard™ ...

Personal Care



Our diverse range of personal care products provide gentle and effective cleansing while minimising cross-infection risks.

Cleaning & disinfection



Developed with convenience in mind, this range provides fast and effective solutions for every cleaning need.

P.P.E.



Promotes infection control and works towards minimising the spread of bacteria, dust and pathogens.

Dry Wipes



Disposable, practical and absorbent, Reynard™ dry wipes are some of the softest and strongest wipes on the market.

Hand Hygiene



Killing 99.9% of germs while leaving hands soft and refreshed. Available with or without chlorhexidine.

Baby Care



Chemical and fragrance-free, our baby care range is designed to be ultra-soft and gentle on little ones.



Reynard™ is reassurance

Australia

44 Hampden Road
Artarmon, NSW, 2064
enquiries@reynardhealth.com.au
1300 667 199

New Zealand

PO Box 8470
Havelock North
sales@reynardhealth.com
0508 494 737

United Kingdom

Boston House, Downsview Rd,
Wantage, OX12 9FF
uk@reynardhealth.com
+44 (0)845 263 7335

International

PO Box 8470
Havelock North
sales@reynardhealth.com
+64 (0)6 650 0709