



medical  
sales &  
service

GYNAECOLOGY



**DTR Medical**<sup>®</sup>



YOUR STERILE SINGLE-USE SOLUTION  
WITHOUT COMPROMISE

# INNOVATION • QUALITY • SERVICE



## Matching reusable standards

### Save Time - Reprocessing

- First time sharpness
- Instrument damage/blockage is limited



## Suitable for high risk procedures

### Save Lives - Contamination

- Cross-infection is eliminated
- Risk is avoided



## Cost effective

### Save Money – Increase Productivity

- Theatre and clinic time IS REDUCED
- Diminished re-processing, repairs and replacement

YOUR STERILE SINGLE-USE SOLUTION

## Contents

<b>Aspirating Needle</b>	<b>4</b>
<b>Spackman Insufflation Cannula</b>	<b>5</b>
<b>Cervical Biopsy Punches</b> <i>Including Rotating &amp; Tischler</i>	<b>6-7</b>
<b>Bonney Uterine Polyp Forceps</b>	<b>8</b>
<b>Forceps</b> <i>Including Vulsellum Forceps &amp; Desjardins</i>	<b>9-10</b>
<b>Scissors</b> <i>Including American Pattern, Braun Stadler &amp; Busch Scissors</i>	<b>10-11</b>
<b>Obstetrics &amp; Maternity Packs</b> <i>Including Delivery &amp; Perineal Packs</i>	<b>12</b>
<b>Dental Syringe</b>	<b>13</b>
<b>Endocervical Specula</b> <i>Including Kogan &amp; Kurihara</i>	<b>13</b>
<b>Metal Cusco Speculum</b> <i>Including with and without smoke tube</i>	<b>14-15</b>
<b>Endometrial Sampling</b>	<b>16</b>
<b>LLETZ Loops &amp; Balls</b>	<b>17</b>

# Cannula and Biopsy Needles



- Stainless steel cannula with 18g needle for laparoscopic aspiration
- Single-use avoids difficult to clean cannula and needle with stepped joint
- Extra long 335mm length suitable for patients of all sizes
- Optimum first-time sharpness of needle for ease of use
- Available in 3.5 and 5mm diameter cannulas complete with Luer connection

## ASPIRATING NEEDLE

ASN2018	Aspirating Needle 335 x 3.5mm (18g) with Luer, double packed	20 units
ASN5018	Aspirating Needle 335 x 3.5mm (18g) with Luer, double packed	50 units
ASN2518	Aspirating Needle 335 x 5mm (18g) with Luer, double packed	20 units



## Spackman Insufflation Cannula

# Spackman Insufflation Cannula



- Single-use removes difficulty of cleaning occluded acorn end
- Adjustable vulsellum plate
- Moveable acorn option available



### Moveable Acorn



### SPACKMAN INSUFFLATION CANNULA

SPC1001	Spackman Insufflation Cannula, double packed	10 units
SPC1001-M	Spackman Insufflation Cannula with Movable Acorn, double packed	10 units



# Cervical Rotating Biopsy Punch



- Ideal for busy clinics
- Single-use eliminates cross contamination
- Metal jaw ensures a sharp, clean cut
- Rotating head provides easier positioning
- Useful for increased HPV testing demands
- Cost effective with no repair bills

“...slick, easy to use and sharp enough to take biopsy without problems”

*Burton Hospitals NHS Foundation Trust*




Did you know this instrument is renewably sourced?


*The handle contains 20% renewable material sourced from technical starch*




metal jaw ensures a sharp, clean cut



29mm tip length




5mm tip width



9mm bite when closed



rotating head secures visibility



ergonomic handle designed for comfort

## CERVICAL ROTATING BIOPSY PUNCH

CBR1001 Cervical Rotating Biopsy Punch, double packed 10 units



MediWales

**Innovation Award Winners**

# Tischler Biopsy Punches



TIME



COST

- High quality stainless steel instruments for a consistent and precise biopsy
- Sharp first time action to minimise trauma

## TISCHLER BIOPSY PUNCHES

CBT1001	Tischler Cervical Biopsy Punch 230mm, double packed	10 units
CBT1001B	Baby Tischler Cervical Biopsy Punch 230mm, double packed	10 units




5.5 x 19mm  
inner

Tischler



Tischler



4.3 x 9mm  
inner

Baby  
Tischler



Baby  
Tischler

# Bonney Uterine Polyp Forceps



- Slender strong jaw suitable for many procedures including polyp and coil removal
- Effective ratchet and robust grip



## BONNEY UTERINE POLYP FORCEPS

BPF2001 Bonney Uterine Polyp Forceps 240mm, double packed 20 units



# Forceps

- Strong jaws with effective locking ratchet and robust grip

## RAMPLEY SPONGE HOLDER

RSH2181	Rampley Sponge Holders 180mm, double packed	20 units
RSH2241	Rampley Sponge Holders 240mm, double packed	20 units



# Vulsellum Forceps

- High quality stainless steel for good tissue grasping
- Jaw with 3:4 teeth and 230mm length
- Curved to aid positioning
- Complimentary to Spackman Insufflation Cannula

## TEALE VULSELLUM FORCEPS

VFT1001	Teale Vulsellum Forceps, double packed	10 units
---------	--	----------




**i** Looking for Spackman Insufflation Cannula?  
see p5

- High quality stainless steel for good tissue grasping
- Jaw with 1:1 teeth and precise action
- Straight 265mm length to suit most procedures
- Complementary to Spackman Insufflation Cannula

## ROTUNDA VULSELLUM FORCEPS

VFR1001	Rotunda Vulsellum Forceps, double packed	10 units
---------	--	----------



# Desjardins Forceps

NEW

- Allows for better grip and hold of polyp/gallstone

## DESJARDINS

DJF2201	Desjardins Forceps 220mm, double packed	20 units
---------	---	----------




# Scissors

NEW


## EPISIOTOMY SCISSORS

WPS2018	Waldmann Perineum Episiotomy Scissors, double packed	20 units
BSS2001	Braun Stadler Episiotomy Scissors Angular 145mm, double packed	20 units



## METZENBAUM SCISSORS

ZSC2121-D	Baby Metzenbaum Scissors Curved 120mm, double packed	20 units
ZSC2141-D	Metzenbaum Scissors Curved 145mm, double packed	20 units
ZSC2161-D	Metzenbaum Scissors Curved 160mm, double packed	20 units
ZSS2181-D	Metzenbaum Scissors Straight 180mm, double packed	20 units
ZSC2181-D	Metzenbaum Scissors Curved 180mm, double packed	20 units
ZSS2201-D	Metzenbaum Scissors Straight 200mm, double packed	20 units
ZSC2201-D	Metzenbaum Scissors Curved 200mm, double packed	20 units



# Scissors

- Wide range for different procedures
- Straight and curved variants


## SCISSORS

USC2001	American Pattern Umbilical Cord Scissors 100mm, double packed	20 units
BUS2002	Busch Umbilical Scissors 160mm, double packed	20 units
MSS2141	Mayo Scissors Straight Bevelled Blade 140mm, double packed	20 units
MSC2141	Mayo Scissors Curved Bevelled Blade 140mm, double packed	20 units
MSS2171	Mayo Scissors Straight Bevelled Blade 170mm, double packed	20 units
MSC2171	Mayo Scissors Curved Bevelled Blade 170mm, double packed	20 units

*Busch Umbilical Scissors*



*American Pattern Umbilical Cord Scissors*



*Straight*



*Curved*

*Mayo Scissors*



# Delivery Packs



TIME



COST

EXTENDED RANGE

OBSTETRICS AND MATERNITY PACKS			
DPWK0505		Delivery Pack double packed 5 units	<ul style="list-style-type: none"> <li>■ Busch Umbilical Scissors 160mm</li> <li>■ Waldmann Perineum Episiotomy Scissors 180mm</li> <li>■ Kocher-Ochsner Forceps 185mm</li> <li>■ Rampley Sponge Holder 240mm x 2</li> </ul>
DPMS1004-D		Delivery Pack double packed 10 units	<ul style="list-style-type: none"> <li>■ Mayo Scissors Curved 170mm</li> <li>■ Spencer Wells Forceps 180mm x 2</li> <li>■ American Pattern Umbilical Cord Scissors 100mm</li> </ul>
DPMS2005-D		Obstetrics Delivery Pack double packed 20 units	<ul style="list-style-type: none"> <li>■ Mayo Scissors Straight 170mm</li> <li>■ Rampley Sponge Holder 240mm</li> <li>■ Spencer Wells Forceps 150mm</li> <li>■ Spencer Wells Forceps 200mm</li> <li>■ American Pattern Umbilical Cord Scissors 100mm</li> </ul>
DPMM1004		Perineal Suture Pack double packed 10 units	<ul style="list-style-type: none"> <li>■ Mayo Hegar Needle Holder 150mm</li> <li>■ Mayo Scissors Straight 140mm</li> <li>■ Mosquito Artery Forceps 122mm</li> <li>■ Adson Toothed Forceps 120mm</li> </ul>
DPIM1006-D		Delivery Pack double packed 10 units	<ul style="list-style-type: none"> <li>■ Mayo Hegar Needle Holder 130mm</li> <li>■ Fine Suture Forceps Serrated 140mm</li> <li>■ Micro Mosquito Forceps Curved 120mm</li> <li>■ Mosquito Forceps Straight 125mm</li> <li>■ Fine Iris Forceps Curved 110mm</li> <li>■ Mayo Scissors Straight 140mm</li> </ul>
MPMC1005		Maternity Pack double packed 10 units	<ul style="list-style-type: none"> <li>■ Mayo Hegar Needle Holder 150mm</li> <li>■ Mayo Scissors Straight 140mm</li> <li>■ Crile Artery Forceps Straight 140mm</li> <li>■ Adson Dissecting Forceps Toothed 120mm</li> <li>■ Adson Dissecting Forceps Non-Toothed 120mm</li> </ul>

# Dental Syringe




TIME

## SPECIAL SIZES

- Suitable for use where needle stick injury is a concern (compliant with EU directive 2010/32/EU – prevention from sharps injuries)
- Popular top loading style to accommodate 2.2ml vials
- Barrel diameter is wide enough to take most vials
- Available in extra long length for use with larger, high BMI patients

### DENTAL SYRINGE

DGS2001	Dental Syringe 85mm for 2.2ml vial, single packed	20 units
DGSL2001	Dental Syringe Long 115mm for 2.2ml vial, single packed	20 units



# Endocervical Specula

- Specula with 265mm reach and 20mm jaws
- Removes the difficulty of cleaning hinge and ratchet
- Scaled ratchet for controlled visibility

### KOGAN ENDOCERVICAL SPECULUM

ECF1001	Kogan Endocervical Specula with ratchet, double packed	10 units
---------	--	----------



COST

### KURIHARA ENDOCERVICAL SPECULUM

KSB1001	Kurihara Endocervical Specula black, double packed	10 units
---------	--	----------



TIME



COST

- Malleable for ease of use
- 300mm shaft with 10mm graduations

### SIMS UTERINE SOUND

USH2001	Sims Uterine Sound, double packed	20 units
---------	-----------------------------------	----------






# Metal Cusco Specula



- Extended range of sizes for all cases
- Strong metal construction with improved insulation performance for patient and user safety
- Insulated and non-insulated options available
- Smooth rounded design for improved patient comfort
- Robust self-retaining screw

## CUSCO INSULATED SPECULA

CSI200S	Cusco Specula Insulated with Smoke Tube Small, double packed	20 units
CSI200M	Cusco Specula Insulated with Smoke Tube Medium, double packed	20 units
CSI200ML	Cusco Specula Insulated with Smoke Tube Medium Long, double packed	20 units
CSI200L	Cusco Specula Insulated with Smoke Tube Large, double packed	20 units
CSI100LL	Cusco Specula Insulated with Smoke Tube Large Long, double packed	10 units
CSI100LE	Cusco Specula Insulated with Smoke Tube Large Extra Long, double packed	10 units
CSI100B	Cusco Specula Insulated with Smoke Tube Broad, double packed	10 units





## CUSCO NON-INSULATED SPECULA

CSM200S	Cusco Specula Metal with Smoke Tube Small, double packed	20 units
CSM200M	Cusco Specula Metal with Smoke Tube Medium, double packed	20 units
CSM200ML	Cusco Specula Metal with Smoke Tube Medium Long, double packed	20 units
CSM200L	Cusco Specula Metal with Smoke Tube Large, double packed	20 units
CSM100LL	Cusco Specula Metal with Smoke Tube Large Long, double packed	10 units
CSM100LE	Cusco Specula Metal with Smoke Tube Large Extra Long, double packed	10 units
CSM100B	Cusco Specula Metal with Smoke Tube Broad, double packed	10 units



Width x Length (mm)	Insulated Smoke Tube	Metal Smoke Tube	Metal w/out Smoke Tube
<b>SMALL</b> 25 X 85mm	CSI200S	CSM200S	CMM200S
<b>MEDIUM</b> 30 X 85mm	CSI200M	CSM200M	CMM200M
<b>MEDIUM LONG</b> 30 X 95mm	CSI200ML	CSM200ML	CMM200ML
<b>LARGE</b> 35 X 95mm	CSI200L	CSM200L	CMM200L
<b>LARGE LONG</b> 35 X 110mm	CSI100LL	CSM100LL	CMM100LL
<b>LARGE XL</b> 35 X 115mm	CSI100LE	CSM100LE	CMM100LE
<b>BROAD</b> 40 X 95mm	CSI100B	CSM100B	CMM100B

### MEASUREMENT INDICATOR



### CUSCO NON-INSULATED WITHOUT SMOKE TUBE

CMM200S	Cusco Specula Metal Small, double packed	20 units
CMM200M	Cusco Specula Metal Medium, double packed	20 units
CMM200ML	Cusco Specula Metal Medium Long, double packed	20 units
CMM200L	Cusco Specula Metal Large, double packed	20 units
CMM100LL	Cusco Specula Metal Large Long, double packed	10 units
CMM100LE	Cusco Specula Metal Large Extra Long, double packed	10 units
CMM100B	Cusco Specula Metal Broad, double packed	10 units



# Endometrial Sampling



TIME



COST

**SIGNIFICANT COST SAVINGS**

- Two and Four Hole Pipette provides improved sampling
- Firm Pipette construction improves insertion
- Endosampler for cases where larger sample size required
- Reliable supply, supported by strong service


## MEDGYN PRODUCTS

022734	Medgyn Pipette II hole	25 units
022730	Medgyn Pipette IV hole	25 units
022720	Medgyn Endosampler	25 units



For cost benefit information please contact us  
+44 (0)1792 797910

Pipette II



Pipette IV



Pipette IV

Pipette II

Endosampler



## NEW RANGE

# LLETZ Loops and Balls



COST

- Standard shaft designed to fit most electro-surgical pencils
- High cutting performance with reduced thermal damage to the tissue sample
- Tungsten loop provides ideal tension for greater control of penetration
- Overcomes cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes



Ball electrode

Loop electrode

### BALL

BE0503	Ball Electrode 150 x 3mm	5 units
BE0505	Ball Electrode 150 x 5mm	5 units

### STANDARD LOOP

028067	Loop Electrode 5 x 5mm	5 units
028069	Loop Electrode 10 x 5mm	5 units
028071	Loop Electrode 10 x 8mm	5 units
028055	Loop Electrode 10 x 10mm	5 units
028057	Loop Electrode 15 x 5mm	5 units
028058	Loop Electrode 15 x 8mm	5 units
028068	Loop Electrode 15 x 10mm	5 units
028078	Loop Electrode 15 x 12mm	5 units
028059	Loop Electrode 15 x 15mm	5 units
028060	Loop Electrode 20 x 8mm	5 units
028063	Loop Electrode 20 x 10mm	5 units
028061	Loop Electrode 20 x 12mm	5 units
028064	Loop Electrode 20 x 15mm	5 units
028065	Loop Electrode 20 x 20mm	5 units
028066	Loop Electrode 25 x 10mm	5 units

## “Supporting visibility, quality and safety in your supply chain”\*\*

Our instruments' packaging features the latest machine readable GS1-standard barcoded labels for product traceability.

The newly incorporated labelling technology is a timely development with an increasing number of health providers worldwide adopting an 'eProcurement' strategy and will provide machine readable product information on all labels to meet expected EU and other regulators demands.

### Associated benefits include:

- Improve patient care and reduce costs
- Improve stock control
- One global language
- Accurately track equipment

*“Unique device identification (UDI) compliance represents a significant development for our customers. Automation of information flow will have influences across the whole purchasing and supply process, providing the benefits of superior instrument traceability and improved patient safety.”*

*“We are working closely with Kodit UDI software to provide our labelling solutions to meet the latest international regulations and we serve as a reference site for other manufacturers to learn from.”*

**Andrew Davidson**  
(Managing Director of DTR Medical)

For more information on our new product labelling, please contact Customer Services on +44 (0) 1792 797910 or [info@dtrmedical.com](mailto:info@dtrmedical.com)

\*\* - GS1 – Global Traceability Standards





## Matching reusable standards

We benchmark against reusable British and international standards:

### **BS 5194-1:1991**

(Inc ISO 7153 – 1:2001) Surgical Instruments - (Metallic Materials and Stainless Steel)

### **BS 5194-2:1989**

Specification for Instruments with Pivot Joints

### **BS 5194-3:1985**

Specification for Dissecting Forceps

### **BS 5194-4:1985**

Specification for Scissors, Shears and other Jointed Instruments



We operate our own UK production facilities equipped with ISO class 7 cleanrooms. Our company is fully accredited to:

**EN ISO 9001:2008**

**EN ISO 13485:2012**

**MDD 93/42/EEC**



**A:** 38 Hector Street OSBORNE PARK WA 6017

**T:** 08 9361 0843 | **F:** 08 9470 3647

**E:** [medsales@medsales.com.au](mailto:medsales@medsales.com.au) | **W:** [www.medsales.com.au](http://www.medsales.com.au)