




medical
sales &
service

ENT/MAXFAX



DTR Medical[®]



YOUR STERILE SINGLE-USE SOLUTION
WITHOUT COMPROMISE

INNOVATION • QUALITY • SERVICE



Matching reusable standards

Save Time - Reprocessing

- First time sharpness
- Instrument damage/blockage is limited



Suitable for high risk procedures

Save Lives - Contamination

- Cross-infection is eliminated
- Risk is avoided



Cost effective

Save Money – Increase Productivity

- Theatre and clinic time IS REDUCED
- Diminished re-processing, repairs and replacement

YOUR STERILE SINGLE-USE SOLUTION

Contents

Ear Specula	4-5
Clearway Suction Handle (for foreign body removal)	6
Tongue Depressors & Nasal Specula	7
Mirrors	7-8
Crocodile Forceps & Scissors <i>Including Fine Jaw Micro Crocodile Forceps</i>	8-10
Forceps <i>Including Adson Toothed Forceps, Fine Suture & Tilley Forceps</i>	10-15
Clamps	15
Bone Trepine	16
Scissors <i>Including Iris & Wire Cutting Scissors</i>	17-19
Tracheal Dilators	19
Suture Packs <i>Including Fine Suture, Stevens Maxillo, & Skin Biopsy Packs</i>	20-23
Probes <i>Including Jobson & Formby</i>	23
Hooks & Points <i>Including Cawthorne</i>	24
Procedure Trays & Packs <i>Including Myringotomy, Ear & Nasal</i>	25-26
Curettes	27
Retractors & Directors <i>Including Maxillofacial Retractor</i>	28
Rhinology & Otology	29-30
Reusable Tympanoplasty Mould	30
Nasal Septal Button & Other	31
Special Suction Handles & Dissectors <i>Including Negus Aspirating Dissector</i>	31-32
Suction Handles <i>Including Frazier, Sinus, Micro, Zoellner (with Fine Ends) & Yankauer Suction</i>	33-38
Micro Needle Electrodes	39
Endoscopic Ear Range	40
Kerrison Rongeur	41
Scope Warmer	41

Ear Specula



TIME

- Sterile and non-sterile options available in a range of sizes
- Matt black finish improves light transmission and visibility
- Smooth edges for patient comfort and safety
- Large sizes (wide and short) suitable for giving injections
- Robust Single-use Specula

EAR SPECULA - OVAL 100 (NON-STERILE)

OSP1001	Oval Ear Specula size 1 (2.5mm), single packed	100 units
OSP1002	Oval Ear Specula size 2 (3.5mm), single packed	100 units
OSP1003	Oval Ear Specula size 3 (4mm), single packed	100 units
OSP1004	Oval Ear Specula size 4 (6mm), single packed	100 units

EAR SPECULA - OVAL 50

OSP0501	Oval Ear Specula size 1 (2.5mm), single packed	50 units
OSP0502	Oval Ear Specula size 2 (3.5mm), single packed	50 units
OSP0503	Oval Ear Specula size 3 (4mm), single packed	50 units
OSP0504	Oval Ear Specula size 4 (6mm), single packed	50 units

EAR SPECULA - OVAL 20

OSP0201	Oval Ear Specula size 1 (2.5mm), single packed	20 units
OSP0202	Oval Ear Specula size 2 (3.5mm), single packed	20 units
OSP0203	Oval Ear Specula size 3 (4mm), single packed	20 units
OSP0204	Oval Ear Specula size 4 (6mm), single packed	20 units



For ear, nasal treatment and Myringotomy trays see p26

EAR SPECULA - ROSEN 100 (NON-STERILE)

ESP1003	Rosen Round Ear Specula 3mm, single packed	100 units
ESP1004	Rosen Round Ear Specula 4mm, single packed	100 units
ESP1005	Rosen Round Ear Specula 5mm, single packed	100 units
ESP1006	Rosen Round Ear Specula 6mm, single packed	100 units
ESP1007	Rosen Round Ear Specula 7mm, single packed	100 units

EAR SPECULA - ROSEN 50

ESP0503	Rosen Round Ear Specula 3mm, single packed	50 units
ESP0504	Rosen Round Ear Specula 4mm, single packed	50 units
ESP0505	Rosen Round Ear Specula 5mm, single packed	50 units
ESP0506	Rosen Round Ear Specula 6mm, single packed	50 units
ESP0507	Rosen Round Ear Specula 7mm, single packed	50 units

EAR SPECULA - ROSEN 20

ESP0203	Rosen Round Ear Specula 3mm, single packed	20 units
ESP0204	Rosen Round Ear Specula 4mm, single packed	20 units
ESP0205	Rosen Round Ear Specula 5mm, single packed	20 units
ESP0206	Rosen Round Ear Specula 6mm, single packed	20 units
ESP0207	Rosen Round Ear Specula 7mm, single packed	20 units



EAR SPECULA - SHEA 100 (NON-STERILE)

SSP1001	Shea Ear Specula size 1 (3.5 x 4.1mm), single packed	100 units
SSP1002	Shea Ear Specula size 2 (4.5 x 5.5mm), single packed	100 units
SSP1003	Shea Ear Specula size 3 (5.5 x 7.5mm), single packed	100 units
SSP1004	Shea Ear Specula size 4 (6.5 x 8.5mm), single packed	100 units
SSP1005	Shea Ear Specula size 5 (7.5 x 9.5mm), single packed	100 units

EAR SPECULA - SHEA 50

SSP0501	Shea Ear Specula size 1 (3.5 x 4.1mm), single packed	50 units
SSP0502	Shea Ear Specula size 2 (4.5 x 5.5mm), single packed	50 units
SSP0503	Shea Ear Specula size 3 (5.5 x 7.5mm), single packed	50 units
SSP0504	Shea Ear Specula size 4 (6.5 x 8.5mm), single packed	50 units
SSP0505	Shea Ear Specula size 5 (7.5 x 9.5mm), single packed	50 units

EAR SPECULA - SHEA 20

SSP0201	Shea Ear Specula size 1 (3.5 x 4.1mm), single packed	20 units
SSP0202	Shea Ear Specula size 2 (4.5 x 5.5mm), single packed	20 units
SSP0203	Shea Ear Specula size 3 (5.5 x 7.5mm), single packed	20 units
SSP0204	Shea Ear Specula size 4 (6.5 x 8.5mm), single packed	20 units
SSP0205	Shea Ear Specula size 5 (7.5 x 9.5mm), single packed	20 units



EAR SPECULA - TUMARKIN 100 (NON-STERILE)

TSP10000	Tumarkin Ear Specula size 0 (3.5mm), single packed	100 units
TSP10005	Tumarkin Ear Specula size 0.5 (4.5mm), single packed	100 units
TSP10010	Tumarkin Ear Specula size 1 (5.5mm), single packed	100 units
TSP10020	Tumarkin Ear Specula size 2 (6mm), single packed	100 units
TSP10030	Tumarkin Ear Specula size 3 (6.5mm), single packed	100 units
TSP10040	Tumarkin Ear Specula size 4 (7mm), single packed	100 units

EAR SPECULA - TUMARKIN 50

TSP05000	Tumarkin Ear Specula size 0 (3.5mm), single packed	50 units
TSP05005	Tumarkin Ear Specula size 0.5 (4.5mm), single packed	50 units
TSP05010	Tumarkin Ear Specula size 1 (5.5mm), single packed	50 units
TSP05020	Tumarkin Ear Specula size 2 (6mm), single packed	50 units
TSP05030	Tumarkin Ear Specula size 3 (6.5mm), single packed	50 units
TSP05040	Tumarkin Ear Specula size 4 (7mm), single packed	50 units

EAR SPECULA - TUMARKIN 20

TSP02000	Tumarkin Ear Specula size 0 (3.5mm), single packed	20 units
TSP02005	Tumarkin Ear Specula size 0.5 (4.5mm), single packed	20 units
TSP02010	Tumarkin Ear Specula size 1 (5.5mm), single packed	20 units
TSP02020	Tumarkin Ear Specula size 2 (6mm), single packed	20 units
TSP02030	Tumarkin Ear Specula size 3 (6.5mm), single packed	20 units
TSP02040	Tumarkin Ear Specula size 4 (7mm), single packed	20 units



Clearway Suction Handle

Foreign body removal




- Removes difficult foreign objects from ear and nose
- Ideal for A&E/ER and MIU
- Less traumatic for patient and family
- Can remove the need for ENT surgeon intervention or theatre time

“After 60 minutes of trying to retrieve a small ball bearing from a child’s ear, we decided to use the Clearway Suction Handle which removed it in seconds! We now stock these as a matter of course and have used on many occasions since, with excellent feedback from our ENT surgeons”.

East Kent Hospitals University NHS Foundation Trust, 2013



special suction cup



rotates for improved clinician visibility



Training poster available on request



See our video demonstration of the Clearway at: www.dtrmedical.com

FOREIGN BODY REMOVAL

SVT1001-S Clearway Suction Handle, single packed 10 units

Tongue Depressors

INCLUDES BRUENINGS

- Robust alternative to wood

LACK

LAC1001S	Lack Tongue Depressor, Small, single packed	20 units
LAC1001M	Lack Tongue Depressor, Medium, single packed	20 units
LAC1001L	Lack Tongue Depressor, Large, single packed	20 units

BRUENINGS

BRU2001	Bruenings Tongue Depressor 190mm, single packed	20 units
---------	---	----------



Nasal Specula



TIME

- Strong flexible stainless steel construction with smooth edges
- Balanced sprung blades for ease of use
- Multipack containing three individual sizes
- Particularly suitable for clinic use

THUDICHUM NASAL SPECULUM

TNS2003	Thudichum Nasal Specula, multipack (sizes 1, 2 & 3), single packed	10 x 3 units
TNS2000	Thudichum Nasal Specula size 0 (Extra Small), single packed	20 units
TNS2011	Thudichum Nasal Specula size 1 (Small), single packed	20 units
TNS2022	Thudichum Nasal Specula size 2 (Medium), single packed	20 units
TNS2033	Thudichum Nasal Specula size 3 (Large), single packed	20 units



RANGE OF SIZES

- Strong stainless steel construction with smooth edges

HARTMANN NASAL SPECULUM

HNS2011-D	Hartmann Nasal Specula, Small, double packed	20 units
HNS2021-D	Hartmann Nasal Specula, Medium, double packed	20 units
HNS2031-D	Hartmann Nasal Specula, Large, double packed	20 units



Mirrors



TIME

DENTAL MIRRORS

DMM2010	Dental Mirror 10mm diameter, single packed	20 units
DMM2015	Dental Mirror 15mm diameter, single packed	20 units
DMM2020	Dental Mirror 20mm diameter, single packed	20 units
DMM2025	Dental Mirror 25mm diameter, single packed	20 units
DMM2030	Dental Mirror 30mm diameter, single packed	20 units



Mirrors



TIME

FULL RANGE

- Robust metal construction with smooth edges and secure mirrors

LARYNGEAL MIRRORS

LYM2010	Laryngeal Mirror 10mm diameter, single packed	20 units
LYM2015	Laryngeal Mirror 15mm diameter, single packed	20 units
LYM2020	Laryngeal Mirror 20mm diameter, single packed	20 units
LYM2025	Laryngeal Mirror 25mm diameter, single packed	20 units
LYM2030	Laryngeal Mirror 30mm diameter, single packed	20 units
LYM5010	Laryngeal Mirror 10mm diameter, single packed	50 units
LYM5015	Laryngeal Mirror 15mm diameter, single packed	50 units
LYM5020	Laryngeal Mirror 20mm diameter, single packed	50 units
LYM5025	Laryngeal Mirror 25mm diameter, single packed	50 units
LYM5030	Laryngeal Mirror 30mm diameter, single packed	50 units



TIME

POST NASAL MIRRORS

PNM2010	Post Nasal Mirror 10mm diameter, single packed	20 units
PNM2012	Post Nasal Mirror 12mm diameter, single packed	20 units
PNM2014	Post Nasal Mirror 14mm diameter, single packed	20 units



You may also be interested in the ClearView Scope Warmer see p41 for more details

Crocodile Forceps



TIME




COST

5MM JAW

- Single-use avoids costly repairs to damaged or worn jaws
- Available in steel and matt black finish
- Eliminates risk of contamination from occluded surfaces
- High quality stainless steel with precise jaw action
- Single or double packed to suit procedure

CROCODILE MICRO FORCEPS

CF2001-S	Crocodile Micro Forceps steel finish, single packed	20 units
CFB2001-S	Crocodile Micro Forceps black finish, single packed	20 units
CF2001-D	Crocodile Micro Forceps steel finish, double packed	20 units
CFB2001-D	Crocodile Micro Forceps black finish, double packed	20 units



Crocodile Fine Jaw Micro Forceps

Fine specification



- Single-use avoids costly repairs to damaged or worn jaws
- Eliminates risk of contamination from occluded surfaces
- Precise jaw action
- Single or double packed to suit procedure



CROCODILE FINE JAW MICRO FORCEPS

CFBF2001-S	Crocodile Fine Jaw Micro Forceps black finish, single packed	20 units
CFBF2001-D	Crocodile Fine Jaw Micro Forceps black finish, double packed	20 units

Crocodile Scissors




TIME



COST


STRAIGHT AND CURVED

- Precise cutting action




CROCODILE SCISSORS

CSS1001-S	Crocodile Scissors Straight stainless steel, single packed	10 units
CSC1001-S	Crocodile Scissors Curved stainless steel, single packed	10 units
CSS1001-D	Crocodile Scissors Straight stainless steel, double packed	10 units
CSC1001-D	Crocodile Scissors Curved stainless steel, double packed	10 units



Curved



4.5mm
blade
length


Straight

Forceps

- Fenestrated cup jaws to aid wax removal
- Eliminates risk of contamination from difficult to clean jaw

HENCKEL AURAL FORCEPS

HAF1001-D	Henckel Aural Forceps, double packed	10 units
-----------	--------------------------------------	----------




8.5mm
top jaw length

7mm
opening of jaw

- High quality stainless steel with precise jaw action
- Matt black finish to minimise reflection
- Single or double packed to suit procedure
- Eliminates risk of contamination from difficult to clean jaw

ORMEROD FORCEPS

OFB1001-S	Ormerod Forceps black finish, single packed	10 units
OFB1001-D	Ormerod Forceps black finish, double packed	10 units



6mm
top jaw length

5mm opening
of jaw

Forceps



TIME




COST

- Elongated shaft to aid procedure

ORAL BIOPSY FORCEPS

OBF1001-S	Oral Biopsy Forceps 90mm stainless steel, single packed	10 units
OBFL1001-S	Oral Biopsy Forceps Large 120mm stainless steel, single packed	10 units




MAGILL CATHETER FORCEPS

MCF11601	Magill Catheter Forceps 160mm, double packed	10 units
MCF12001	Magill Catheter Forceps 200mm, double packed	10 units
MCF12501	Magill Catheter Forceps 250mm, double packed	10 units



NEW SMALLER SIZE

MCF11601



Forceps

- Effective locking ratchet and robust grip


CRILE AND MOSQUITO FORCEPS

CFS2001-D	Crile Artery Forceps 140mm Straight, double packed	20 units
CFC2001-D	Crile Artery Forceps 140mm Curved on Flat, double packed	20 units
MQT5001-S	Mosquito Forceps 125mm Straight, single packed	50 units
MQC2012	Micro Mosquito Forceps 120mm Curved, double packed	20 units
MQC2025-D	Mosquito Forceps Curved 125mm, double packed	20 units
MQC2125-D	Mosquito Forceps Curved Toothed 1x2 125mm, double packed	20 units



KOCHER ARTERY FORCEPS

KAF2014-D	Kocher Artery Forceps 140mm, double packed	20 units
-----------	--	----------



NEW

- Effective at clamping off bleeding from blood vessels


KOCHER-OCHSNER FORCEPS

KOF2181	Kocher-Ochsner 1x2 Toothed Forceps Straight 185mm, double packed	20 units
---------	--	----------



SPENCER WELLS FORCEPS


SWF2010	Spencer Wells Artery Forceps 120mm Straight, double packed	20 units
SWF2015	Spencer Wells Artery Forceps 150mm Straight, double packed	20 units
SWF2018	Spencer Wells Artery Forceps 180mm Straight, double packed	20 units
SWF2020	Spencer Wells Artery Forceps 200mm Straight, double packed	20 units



- Precision instrument with firm grip
- Good spring action

DENTAL FORCEPS

CDF2001 College Dental Forceps 160mm, single packed 20 units




NEW RANGE

ADSON AND McINDOE FORCEPS

ADS5001-S	Adson Non-Toothed Forceps 120mm, single packed	50 units
ATS5001-S	Adson Toothed Forceps 120mm, single packed	50 units
ATP2012-D	Adson Forceps 1x2 with plateau 120mm, double packed	20 units
ATP2015-D	Adson Forceps 1x2 with plateau 150mm, double packed	20 units
ATS2012F	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm, double packed	50 units
MDF5001-S	McIndoe Dissecting Forceps 150mm, Non-Toothed, single packed	50 units



Adson Toothed



Adson Non-Toothed

1.9mm tip width



We also do Adson Bipolar Forceps see general catalogue p15

Adson Toothed Forceps 120mm




Adson Forceps 150mm with plateau




Adson Forceps 120mm with plateau



Fine Tip Adson Toothed Forceps



Adson Forceps



McIndoe Forceps 150mm




Forceps




TIME

FINE SUTURE FORCEPS

SPF2001-S	Fine Suture Forceps Plain 105mm, single packed	20 units
SSF2001-S	Fine Suture Forceps Serrated 105mm, single packed	20 units
SPF2001-D	Fine Suture Forceps Plain 105mm, double packed	20 units
SSF2001-D	Fine Suture Forceps Serrated 105mm, double packed	20 units



Plain




Serrated




Looking for a Fine Suture Pack?
see p20

DEBAKEY FORCEPS

DAF2016-D	DeBakey Atraumatic Dissecting Forceps 160mm, double packed	20 units
-----------	--	----------




2 x 17mm jaws



GILLIES FORCEPS

GDF2001-S	Gillies Toothed Dissecting Forceps 160mm, single packed	20 units
-----------	---	----------



20mm
2.4mm tip width




ALLIS TISSUE FORCEPS

ATF2341	Allis Tissue Forceps 3:4 teeth, 160mm, double packed	20 units
---------	--	----------




4.8mm width of head



NEW

STRASSBURG FORCEPS


SFT2015-D Strassburg Forceps Toothed 1x2 150mm, double packed 20 units



- Strong slender jaws with good grip and tips that always meet
- Single-use avoids damage or difficult to clean box joints
- Particularly suitable for clinic use

PACKING FORCEPS


TAF2001-S	Tilley Aural Forceps 55mm, single packed	20 units
TNF2001-S	Tilley Nasal Dressing Forceps 75mm, single packed	20 units
WAF2001-S	Wilde Aural Packing Forceps 150mm, single packed	20 units



Clamps

CLAMPS

BTC2013-D	Backhaus Towel Clamp 135mm, double packed	20 units
TCF2051	Tubing Clamps 150mm, double packed	20 units



i You may also be interested in Silicone Tubing see general catalogue p28-31

NEW

Bone Trephine

- Used for bone grafting, including taking circular bone graft biopsies for paediatric cleft surgery
- Less invasive
- First time sharpness
- Less patient trauma

T-Handle

- Reusable
- Used to hold Bone Trephine Burr to allow for coring

“I am delighted to be using the trephine with T bar handle for alveolar bone grafting now improved and produced by DTR Medical. It allows a relatively minimal access to harvesting bone from the iliac crest yet with a high bone yield and minimises the external scar. My patients are up and walking the same day after surgery”.

*Alistair RM Cobb. Consultant Surgeon.
Southwest UK Cleft service*

robust handle

easy insertion for Bone Trephine



BONE TREPHINE

BTB1002 Sterile Single-use Bone Trephine Burr Procedure Pack, inc Bone Trephine Burr & No.11 Disposable Scalpel, double packed 10 units

T-HANDLE

LX203R T-Handle with Jacobs Chuck Diameter 6.5mm 1 unit

Scissors




COST

- Stainless steel sharp/sharp scissors
- First time sharpness which cut even at the tips



KILNER SCISSORS

KSS2001-S Kilner Scissors 116mm, single packed 20 units



IRIS SCISSORS

ISS2111-S	Iris Scissors Straight 90mm, single packed	20 units
ISC2111-S	Iris Scissors Curved 90mm, single packed	20 units
ISS2110-D	Iris Scissors Straight 110mm, double packed	20 units
ISC2110-D	Iris Scissors Curved 110mm, double packed	20 units



METZENBAUM SCISSORS


ZSC2121-D Baby Metzenbaum Scissors Curved 120mm, double packed 20 units




Scissors

STEVENS TENOTOMY SCISSORS


TSS2001-S	Stevens Tenotomy Scissors 115mm Straight, single packed	20 units
TSC2001-S	Stevens Tenotomy Scissors 115mm Curved, single packed	20 units
TSS2001-D	Stevens Tenotomy Scissors 115mm Straight, double packed	20 units
TSC2001-D	Stevens Tenotomy Scissors 115mm Curved, double packed	20 units



Straight



Curved



We also do Stevens Maxillo Suture packs see table p22



COST

WIRE CUTTING SCISSORS

WCS2001-S	Wire Cutting Scissors 120mm, single packed	20 units
-----------	--	----------



NEW

REYNOLDS SCISSORS

RSC2014-D	Reynolds Scissors Curved 145mm, double packed	20 units
-----------	---	----------




Scissors




- Both straight options have 0.4mm tip width
- Curved option offers greater flexibility for delicate procedures

VANNAS SCISSORS


VSS2081	Vannas Scissors Straight 85mm, double packed	20 units
VSS2121	Vannas Scissors Straight 125mm, double packed	20 units
VSC2081	Vannas Scissors Curved 85mm, double packed	20 units



VSS2081




VSS2121




VSC2081

MAYO SCISSORS

MSS2141	Mayo Scissors Straight Bevelled Blade 140mm, double packed	20 units
MSC2141	Mayo Scissors Curved Bevelled Blade 140mm, double packed	20 units
MSS2171	Mayo Scissors Straight Bevelled Blade 170mm, double packed	20 units
MSC2171	Mayo Scissors Curved Bevelled Blade 170mm, double packed	20 units



Straight



Curved




Tracheal Dilators



- Smooth mechanism to ensure high quality performance
- Adult and paediatric sizes

TRACHEAL DILATORS

TTD2101	Trouseau Tracheal Dilator 105mm, double packed	20 units
TTD2141	Trouseau Tracheal Dilator 140mm, double packed	20 units









Suture Packs









■ Customised and double packed options available







NEW LINES FOR 2016

SUTURE PACKS			
FSP2002-S		Fine Suture Pack single packed 20 units	<ul style="list-style-type: none"> ■ Iris Scissors Straight 90mm ■ Fine Suture Forceps Serrated 105mm
SPA1004		Suture Pack single packed 10 units	<ul style="list-style-type: none"> ■ Iris Scissors Curved 90mm ■ Adson Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 150mm ■ No3 Blade Handle
SPIA2002-S		Suture Pack single packed 20 units	<ul style="list-style-type: none"> ■ Adson Non-Toothed Forceps 120mm ■ Iris Scissors Curved 90mm
SPIA2003-S		Suture Pack single packed 20 units	<ul style="list-style-type: none"> ■ Iris Scissors Straight 90mm ■ Adson Non-Toothed Forceps 120mm ■ Mosquito Forceps Straight 125mm
SPMA1005		Suture Pack double packed 10 units	<ul style="list-style-type: none"> ■ Mayo Scissors Straight 140mm ■ Adson Toothed Forceps 120mm ■ Adson Non-Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 130mm ■ Micro Mosquito Forceps 120mm Curved
SPMM1005		Suture Pack double packed 10 units	<ul style="list-style-type: none"> ■ Mayo Scissors Straight 140mm ■ Adson Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 130mm ■ McIndoe Dissecting Forceps Non-Toothed 150mm ■ Micro Mosquito Forceps Curved 120mm



SUTURE PACKS

<p>SPIA1004</p>		<p>Suture Pack single packed 10 units</p>	<ul style="list-style-type: none"> ■ Iris Scissors Straight 90mm ■ Adson Toothed Forceps 120mm ■ McIndoe Dissecting Forceps Non-Toothed 150mm ■ Mayo Hegar Needle Holder 150mm
<p>SPIA1005</p>		<p>Suture Pack single packed 10 units</p>	<ul style="list-style-type: none"> ■ Iris Scissors Straight 90mm ■ Adson Toothed Forceps 120mm ■ Mosquito Forceps Straight 125mm ■ Mayo Hegar Needle Holder 130mm ■ No3 Blade Handle
<p>SPKA2002-S</p>		<p>Suture Pack single packed 20 units</p>	<ul style="list-style-type: none"> ■ Kilner Scissors 120mm ■ Adson Non-Toothed Forceps 120mm
<p>SPKA1004</p>		<p>Suture Pack single packed 10 units</p>	<ul style="list-style-type: none"> ■ Kilner Scissors 120mm ■ Adson Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 150mm ■ No3 Blade Handle
<p>STKA1004</p>		<p>Suture Tray double packed 10 units</p>	<ul style="list-style-type: none"> ■ Kilner Scissors 120mm ■ Adson Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 150mm ■ No3 Blade Handle
<p>SPKA1136</p>		<p>Suture Pack double packed 10 units</p>	<ul style="list-style-type: none"> ■ Kilner Scissors 120mm ■ Adson Toothed Forceps 120mm ■ Adson Non-Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 130mm ■ Mosquito Forceps Straight 125mm ■ Micro Mosquito Forceps Curved 120mm

Suture Packs

SUTURE PACKS			
SPKA1156		Suture Pack double packed 10 units	<ul style="list-style-type: none"> ■ Kilner Scissors 120mm ■ Adson Toothed Forceps 120mm ■ Adson Non-Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 150mm ■ Mosquito Forceps Straight 125mm ■ Micro Mosquito Forceps 120mm Curved
SPKM2003-S		Suture Pack single packed 20 units	<ul style="list-style-type: none"> ■ Kilner Scissors 120mm ■ McIndoe Dissecting Forceps Non-Toothed 150mm ■ Mosquito Forceps Straight 125mm
SPSA1006		Stevens Maxillo Suture Pack double packed 10 units	<ul style="list-style-type: none"> ■ Stevens Tenotomy Scissors Straight 115mm ■ Stevens Tenotomy Scissors Curved 115mm ■ Kilner Scissors 120mm ■ Adson Toothed Forceps 120mm ■ Adson Non-Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 130mm
SPSA1007		Stevens Maxillo Suture Pack double packed 10 units	<ul style="list-style-type: none"> ■ Stevens Tenotomy Scissors Straight 115mm ■ Stevens Tenotomy Scissors Curved 115mm ■ Kilner Scissors 120mm ■ Adson Toothed Forceps 120mm ■ Adson Non-Toothed Forceps 120mm ■ Mayo Hegar Needle Holder 130mm ■ No3 Blade Handle
SBIG1004		Skin Biopsy Pack double packed 10 units	<ul style="list-style-type: none"> ■ Iris Scissors Straight 90mm ■ Gillies Toothed Forceps 160mm ■ Kilner Needle Holder 135mm ■ No3 Blade Handle
SBIG1005		Skin Biopsy Pack double packed 10 units	<ul style="list-style-type: none"> ■ Iris Scissors Straight 90mm ■ Gillies Toothed Forceps 160mm ■ Kilner Needle Holder 135mm ■ No3 Blade Handle ■ Gillies Skin Hook Sharp

Suture Packs

SUTURE PACKS			
MFP1003-S		Maxillofacial Pack single packed 10 units	<ul style="list-style-type: none"> ■ Iris Scissors Straight 110mm ■ Crilewood Needle Holder 150mm ■ Adson Forceps 1x2 Toothed with Plateau
TTP2002		Tongue Tie Pack double packed 20 units	<ul style="list-style-type: none"> ■ Baby Metzenbaum Scissors Curved 120mm ■ Grooved Director 130mm



Probes



- Easy to bend to shape to match anatomy and improve vision
- Fine working ends for use in even the smallest ear specula
- Jobson double ended loop and cotton wool grip
- Single-use avoids need to re-straighten or risk breakage during use





JOBSON		
JHPS5001-S	Jobson Horne Probe 127mm (5 inch), single packed	50 units
JHPS10001-S	Jobson Horne Probe 127mm (5 inch), single packed	100 units
JHPL5001-S	Jobson Horne Probe 178mm (7inch), single packed	50 units
JHPL10001-S	Jobson Horne Probe 178mm (7inch), single packed	100 units
JHPL5001-D	Jobson Horne Probe 178mm (7inch), double packed	50 units



i Looking for Ear Specula?
see p4-5

FORMBY		
FCH2001-S	Formby Cerumen Hook and Probe 178mm, single packed	20 units



Hooks and Points



TIME



COST

- Choice of 90 degree hook or point

CAWTHORNE


CTH2901-D	Cawthorne Hook 1.5mm 90 degree black, double packed	20 units
CTP2001-D	Cawthorne Straight Point black, double packed	20 units



Hook



Point



GILLIES

GHS2001	Gillies Skin Hook Sharp 175mm, double packed	20 units
GHB2001	Gillies Skin Hook Semi Sharp 175mm, double packed	20 units



Sharp




Semi-Sharp



SCHUKNECHT PICK






SPS2001	Schuknecht Pick Straight 160mm, double packed	20 units
---------	---	----------



Procedure Trays and Packs




■ Customised options available

NASAL TREATMENT			
NTP0301		Nasal Treatment Pack double packed, 10 units	<ul style="list-style-type: none"> ■ Tilley Nasal Forceps 75mm ■ Lack Tongue Depressors Large ■ Thudichum Nasal Specula size 2
NPP1006		Nasal Plastics Procedure Pack double packed 10 units	<ul style="list-style-type: none"> ■ Black Crocodile Forceps ■ Thudichum Nasal Specula size 2 ■ Jobson Horne Probe 5" ■ Kilner Scissors 120mm ■ Adson Non-Toothed Forceps 120mm ■ Zoellner Suction Handle
EAR TREATMENT			
EUMT4001		EUM Procedure Pack single packed 10 units	<ul style="list-style-type: none"> ■ Crocodile Micro Forceps Black ■ Tumarkin size 0.5 & 1 ■ Jobson Horne Probe 7"
ETP0601		Ear Treatment Pack inc Black Crocodile Forceps & Oval Ear Specula double packed 10 units	<ul style="list-style-type: none"> ■ Crocodile Micro Forceps Black ■ Jobson Horne Probe 7" ■ Formby Cerumen Hook ■ Oval Ear Specula sizes 2, 3, & 4
ETP2006		Ear Treatment Pack inc Crocodile Fine Jaw Micro Forceps & Shea Ear Specula double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Suction Regulator ■ Crocodile Fine Jaw Micro Forceps ■ Shea Ear Specula sizes 1, 3, & 5

Procedure Trays and Packs



■ Customised options available

EAR TREATMENT - MYRINGOTOMY			
MMS1017-F	 <p style="text-align: center;">MMT1017</p>	Myringotomy and Minor Ear Tray No 1 inc Shea Ear Specula sizes 2 & 3 double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Crocodile Fine Jaw Micro Forceps ■ Shea Ear Specula sizes 2 & 3 ■ Zoellner Suction Handle ■ Zoellner Fine End 18swg ■ Narrow Arrow Myringotome Blade
MMT1017		Myringotomy and Minor Ear Tray No 1 double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Crocodile Micro Forceps Black ■ Tumarkin Ear Specula sizes 0 & 0.5 ■ Zoellner Suction Handle ■ Zoellner Fine End 18swg ■ Narrow Arrow Myringotome Blade
MMT1017-F		Myringotomy and Minor Ear Tray No 1 inc Extra Fine Crocodile Forceps double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Crocodile Fine Jaw Micro Forceps ■ Tumarkin Ear Specula sizes 0 & 0.5 ■ Zoellner Suction Handle ■ Zoellner Fine End 18swg ■ Narrow Arrow Myringotome Blade
MMT1028-F		Myringotomy and Minor Ear Tray No 2 inc Extra Fine Crocodile Forceps double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Crocodile Fine Jaw Micro Forceps ■ Oval Ear Specula size 2 ■ Zoellner Suction Handle ■ Zoellner Fine End 20swg ■ Narrow Arrow Myringotome Blade ■ Formby Cerumen Hook ■ Cawthorne Point
MMT1037C		Myringotomy Tray No 3 Children's double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Crocodile Micro Forceps Black ■ Oval Ear Specula size 1 ■ Zoellner Suction Handle ■ Zoellner Fine End 18swg ■ Narrow Arrow Myringotome Blade ■ Cawthorne Point
MMT1017A		Adult Myringotomy Kit double packed 10 units	<ul style="list-style-type: none"> ■ Jobson Horne Probe 7" ■ Crocodile Micro Forceps Black ■ Tumarkin Specula sizes 1 & 3 ■ Zoellner Suction Handle ■ Zoellner Fine End 18swg ■ Narrow Arrow Myringotome Blade

Curettes

NEW LINES FOR 2016

- Precise and sharp scoop for effective curettage

BECKMANN ADENOID CURETTE

BAC2010	Beckmann Adenoid Curette, 10mm Curette Edge, 220mm, double packed	20 units
BAC2013	Beckmann Adenoid Curette, 13mm Curette Edge, 220mm, double packed	20 units
BAC2015	Beckmann Adenoid Curette, 15mm Curette Edge, 220mm, double packed	20 units
BAC2017	Beckmann Adenoid Curette, 17mm Curette Edge, 220mm, double packed	20 units



COMING SOON

HOUSE CURETTE

HCD20515	House Curette 180mm, double ended 0.5/1.5mm Scoops, double packed	20 units
----------	---	----------



- First time sharpness

ROSEN NEEDLE

RSN2016	Rosen Needle 160mm, double packed	20 units
---------	-----------------------------------	----------



Contact us for more information on
+44 (0) 1792 797910


Retractors



- Stainless steel dual hook with blunt points
- Latex-free 2mm silicone 180mm ties
- Pack of 4 for wound size flexibility

MAXILLOFACIAL RETRACTOR



MFR0410 Maxillofacial Retractors x 4 with silicone loop, double packed 10 units



- Rigid stainless steel construction
- Four hook proximal screw action

ALMS RETRACTOR

ALR1071 Alms Retractor Small 78mm length, double packed 10 units
 ALR1101 Alms Retractor Large 103mm length, double packed 10 units



- Double ended
- First time sharpness for precision

KILNER "CATSPAW" RETRACTOR

KCP2001 Kilner "Catspaw" Retractor, double packed 20 units



Directors

DIRECTORS

GDR2001 Grooved Director 130mm, double packed 20 units



i Included in Tongue Tie Pack see p23

Rhinology and Otolaryngology




- Range of Single-use blade formats
- First time sharpness for precision
- Ability to bend to required angle
- Individually packed to avoid cross contamination

MYRINGOTOME BLADES AND HANDLES

BL-1150	MicroEdge Blade with adapter, Arrow	12 units
BL-1151	MicroEdge Blade with adapter, Bevel	12 units
BL-1152	MicroEdge Blade with adapter, Narrow Arrow	6 units
BL-1425	MicroEdge Blade Handle, reusable, non-sterile	1 unit
DBL-1425	Myringotome Handle used with MicroEdge Blades, double packed	12 units
MNA2001	Myringotome Narrow Arrow, double packed	20 units



Myringotome Handle and MicroEdge Blade



Myringotome Narrow Arrow




For range of Ear Treatment Packs see p25-26

- Complies with current UK decontamination policies
- First time sharpness for precision

TILLEY LICHTWITZ SINUS PUNCH

TLP1001	Tilley Lichtwitz Sinus Punch, double packed	10 units
---------	---	----------



1mm x 2.5mm blade




Looking for suction handles? see p33-38

- Comprehensive range of designs, sizes and materials
- Available with or without wires

VENTILATION TUBES AND GROMMETS

VT-0201-01	Shepard vent tubes with wire 1.14mm lumen fluoroelastic	6 units
VT-0202-01	Shepard vent tubes 1.14mm lumen fluoroelastic	6 units
VT-0203-01	Shepard vent tubes with wire 1.0mm lumen fluoroelastic	6 units
VT-0204-01	Shepard vent tubes 1.0mm lumen fluoroelastic	6 units
VT-0802-01	Shah with tab 1.14mm lumen fluoroelastic	6 units
VT-0802-50	Shah with tab 1.14mm lumen fluoroelastic	50 units
VT-0803-01	Shah with wire 1.14mm lumen fluoroelastic	6 units
VT-0803-50	Shah with wire 1.14mm lumen fluoroelastic	50 units
VT-0804-01	Shah 1.14mm lumen fluoroelastic	6 units
VT-0804-50	Shah 1.14mm lumen fluoroelastic	50 units
VT-0400-01	T-tube 12 x 10 x 1.14mm lumen silicone	6 units
VT-0401-01	T-tube 5 x 8 x 1.14mm lumen silicone	6 units
VT-0403-01	T-tube Grommet 1.27mm lumen silicone	6 units
VT-1207-01	Reuter Bobbin Short Flange 1.0mm lumen fluoroelastic	6 units
VT-0300-01	Paparella 1.14mm lumen silicone	6 units
VT-0501-01	Armstrong bevelled straight shank 7.5 x 1.14mm lumen fluoroelastic	6 units
VT-0501-50	Armstrong bevelled straight shank 7.5 x 1.14mm lumen fluoroelastic	50 units



Reuter bobbin

Shah (with wire)



Armstrong bevelled straight shank

T-tube



Rhinology and Otology



- Effective method to deliver drug treatment to inner ear
- Standard length wick for Ménière's disease
- Longer wick for eustachian tube dysfunction
- Minimises expensive hospital time for repeat doses



SILVERSTEIN MICROWICK

OT-1700	Silverstein MicroWick (9 x 1 mm)	1 unit
OT-1720	Silverstein MicroWick ET (20 x 1 mm)	1 unit

Tympanoplasty

REUSABLE

- Reusable moulds for faciaform tympanoplasty
- Anatomically shaped for left and right tympanic membranes
- Labelled for ease of use

TYMPANOPLASTY MOULDS

TYM0102R	Reusable Tympanoplasty Mould, double packed 1 Set of Right & Left	1 Set
----------	---	-------



Splints

- Non-stick splint reduces adhesion
- Slit improves insertion and removal
- Pre-punched holes provided to ease suturing
- Available in two sizes and thicknesses

NASAL SPLINTS

SP-27000	Bivalve Nasal Splint, Standard Thin	2 x 6 units
SP-27005	Bivalve Nasal Splint, Standard Thick	2 x 6 units
SP-27010	Bivalve Nasal Splint, Oversize Thin	2 x 6 units
SP-27015	Bivalve Nasal Splint, Oversize Thick	2 x 6 units
SP-1500-L	Denver Nasal Splint Kit Large	5 units
SP-1500-S/M	Denver Nasal Splint Kit Small/Medium	5 units



EASY INSERTION AND REMOVAL

- Funnel-shaped airway enhances patient comfort
- Unique airway geometry creates pressure gradient increasing airflow by 50%
- Easier breathing and minimises blockage
- Suture indicators on anterior tip aids suture placement

BREEZE NASAL AIRWAY SPLINTS

SP-79300	Breeze Nasal Airway Splint	2 x 1 units
----------	----------------------------	-------------



Nasal Septal Button and Other



COST EFFECTIVE

- Closure of septal perforations to reduce or eliminate symptoms
- Cost effective alternative to existing products
- Soft medical grade silicone for ease of placement and patient comfort
- Easily trimmed to ensure proper fit


NASAL SEPTAL BUTTONS

SP-78100	Nasal Septal Button Standard (3cm)	1 unit
SP-78105	Nasal Septal Button Oversize (5cm)	1 unit

- Doyle patient friendly packing slit to aid removal and reduce discomfort

OTHER RHINOLOGY PRODUCTS

SP-80020	Shikani Stent Medium	6 units
RH-7411-10	Anatomical Nasal Packing 8.0 x 1.0 x 3cm	10 units
RH-7412-10	Doyle Patient Friendly Nasal Packing 8 x 1.5 x 3cm with string	10 units
RH-7413-10	Doyle Patient Friendly Nasal Packing with airway tube & string 8 x 1.5 x 3cm	10 units
RH-7414-10	Doyle Patient Friendly Nasal Packing, Thin 8 x 1.0 x 3cm with string	10 units
RH-7530-10	Post-Op Sinus Packing with String 3.5 x 0.6 x 1.2cm	10 units
RH-2315-20	Small Epistaxis Pack with String 5.5 x 1.5 x 2.5cm	20 units
OT-5200-25	Ear Wick 9 x 15mm	25 units
NGC1001	Nasal Gate Clamp, single packed	10 units



Special Suction Handles and Dissectors



- Negus, Freer "Hockey Stick" for nasal procedures
- Tonsil Dissector based on the Gwynne Evans blade pattern
- Optimum blade sharpness for precision work
- Removes risk from very difficult to clean suction blades

ASPIRATING DISSECTORS

FDS0501	Freer Aspirating "Hockey Stick" Dissector, double packed	5 units
TDS0501	Tonsil Aspirating Dissector, double packed	5 units



Freer



Tonsil



IMPROVED NEGUS

Negus Aspirating Dissector



- Optimum blade sharpness for precision work
- Removes risk from very difficult to clean suction blades
- Enlarged hole size increases the rear aperture size
- Front edge of hole is 0.75mm closer to the blade tip

Improved suction flow at the blade tip giving better visibility

Less prone to blockage and subsequent theatre delays



BEFORE
1.5mm hole



AFTER
2.25mm hole



ASPIRATING DISSECTORS

NDS0501 Negus Aspirating Dissector, double packed 5 units

Suction Handles



- Complete range of sizes from 3 to 7fg with suction control
- Working length of 77mm
- Mini Barron 7fg size suitable for use with fine end

BARRON SUCTION HANDLES

BSH2003	Barron Suction Handle 1mm (3fg), double packed	20 units
BSH2004	Barron Suction Handle 1.35mm (4fg), double packed	20 units
BSH2005	Barron Suction Handle 1.67mm (5fg), double packed	20 units
BSH2006	Barron Suction Handle 2mm (6fg), double packed	20 units
BSH2007	Barron Suction Handle 2.3mm (7fg), double packed	20 units
BLS2101	Mini Barron Suction Handle 98 x 2.33mm OD (7fg), double packed	20 units

Mini Barron Suction Handle



Barron Suction Handle



Looking for sterile silicone tubing for use with the following? general catalogue see p29



- Removes need to clean and reprocess the long narrow lumen

BELLUCCI SUCTION HANDLES

BLS2001	Bellucci Suction Handle 175 x 2mm OD (6fg), double packed	20 units
---------	---	----------



- Effective alternative to reprocessing very long narrow lumen
- Range of sizes from 240mm (6fg) to 550mm (18fg)

LONG SUCTION HANDLES

240mm		
BLSF2240	Micro Laryngeal Suction Handle 240 x 2mm OD (6fg), double packed	20 units
280mm		
BLS2280	Micro Laryngeal Suction Handle 280 x 2mm OD (6fg), double packed	20 units
330mm		
BLS2330	Micro Laryngeal Suction Handle 330 x 2mm OD (6fg), double packed	20 units
BLSF2339	Micro Laryngeal Suction Handle 330 x 3mm OD (9fg), double packed	20 units
360mm		
BLSF2369	Micro Laryngeal Suction Handle, 360 x 3mm OD (9fg), double packed	20 units
450mm		
BLS2450	Bronchoscopy Suction Handle 450 x 4mm OD (12fg), double packed	20 units
550mm		
BLS2550	Oesophageal Suction Handle 550 x 4mm OD (12fg), double packed	20 units



Frazier Suction Handles

Available without suction control



- Strong construction suitable for bending to shape during procedure
- Single-use avoids the need to restraighten or risk breakage
- Tip options include non-fenestrated, fenestrated with two holes and cross fenestrated with four holes

Supplied straight*
Bend to shape for your procedure

*Excluding 12fg range



with suction control



cross fenestrated



without suction control



with suction control



42° bend



without suction control

FRAZIER SUCTION HANDLES

2mm (6fg)

FSH2006	Frazier Suction Handle 184 x 2mm OD (6fg) straight, double packed	20 units
FSF2006	Frazier Suction Handle 184 x 2mm OD (6fg) straight, fenestrated, double packed	20 units
FSF2046	Frazier Suction Handle 184 x 2mm OD (6fg) straight, cross fenestrated, double packed	20 units
FSN2426	Frazier Suction Handle 184 x 2mm OD (6fg), 42° bend, without suction control, double packed	20 units

2.66mm (8fg)

FSH2008	Frazier Suction Handle 184 x 2.66mm OD (8fg) straight, double packed	20 units
FSF2008	Frazier Suction Handle 184 x 2.66mm OD (8fg) straight, fenestrated, double packed	20 units
FSF2048	Frazier Suction Handle 184 x 2.66mm OD (8fg) straight, cross fenestrated, double packed	20 units

3mm (9fg)

FSH2009	Frazier Suction Handle 184 x 3mm OD (9fg) straight, double packed	20 units
FSF2009	Frazier Suction Handle 184 x 3mm OD (9fg) straight, fenestrated, double packed	20 units
FSF2049	Frazier Suction Handle 184 x 3mm OD (9fg) straight, cross fenestrated, double packed	20 units
FSN2429	Frazier Suction Handle 184 x 3mm OD (9fg), 42° bend, without suction control, double packed	20 units

4mm (12fg)

FSH2012	Frazier Suction Handle 184 x 4mm OD (12fg), double packed	20 units
FSF2012	Frazier Suction Handle 184 x 4mm OD (12fg), fenestrated, 150° bend, double packed	20 units
FSF2412	Frazier Suction Handle 184 x 4mm OD (12fg), cross fenestrated, 150° bend, double packed	20 units
FSN24212	Frazier Suction Handle 184 x 4mm OD (12fg), 42° bend, without suction control, double packed	20 units



Suction Handles



- Atraumatic tip

MAGILL SUCTION HANDLES

MGH2010	Magill Suction Handle 1mm Lumen 210 x 3mm OD, double packed	20 units
MGH2015	Magill Suction Handle 1.5mm Lumen 210 x 3.5mm OD, double packed	20 units
MGH2020	Magill Suction Handle 2mm Lumen 210 x 4mm OD, double packed	20 units



- Metal version of traditional glass suction handle
- Suitable for where wide bore is required



METAL "GLASS-TYPE" SUCTION HANDLE

GSH2001	Metal "Glass-Type" Suction handle 200 x 5mm OD (15fg), double packed	20 units
---------	--	----------



- Dual purpose suction and irrigation with suction control
- Complete range of sizes from 1.8 to 5mm
- Single-use overcomes repairs to tips damaged through drilling
- Particularly suitable for collection of bone dust

HOUSE SUCTION HANDLES


HST0118/18	House Suction and Irrigation Handle (1.8 x 1.8mm), double packed	20 units
HST0125/20	House Suction and Irrigation Handle (2.5 x 2mm), double packed	20 units
HST0130/20	House Suction and Irrigation Handle (3 x 2mm), double packed	20 units
HST0145/20	House Suction and Irrigation Handle (4.5 x 2mm), double packed	20 units
HST0145/30	House Suction and Irrigation Handle (4.5 x 3mm), double packed	20 units
HST0145/45	House Suction and Irrigation Handle (4.5 x 4.5mm), double packed	20 units
HST0150/45	House Suction and Irrigation Handle (5 x 4.5mm), double packed	20 units



Suction Handles

FENESTRATED HOUSE SUCTION HANDLE

HSFL118/18 House Suction Handle Fenestrated (1.8 x 1.8mm), double packed 20 units



EXTENDED RANGE

- Stainless steel tube with colour coded luer connector for microsuction use
- Also suitable for use with syringe as irrigation device
- Available in 86, 110 and 120mm lengths

MICRO SUCTION TUBES

MSD2001	Micro Suction Device 86 x 2mm OD (6fg), double packed	20 units
86mm		
MSV508113	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g), double packed	50 units
MSV508113R	Micro Suction (Verhoeven) Red 86 x 2.3mm OD (13g) with Suction Regulator, double packed	50 units
MSV508114	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g), double packed	50 units
MSV508114R	Micro Suction (Verhoeven) Green 86 x 2.1mm OD (14g) with Suction Regulator, double packed	50 units
MSV508116	Micro Suction (Verhoeven) Cream 86 x 1.65mm OD (16g), double packed	50 units
MSV508116R	Micro Suction (Verhoeven) Cream 86 x 1.65mm OD (16g) with Suction Regulator, double packed	50 units
MSV508117	Micro Suction (Verhoeven) White 86 x 1.47mm OD (17g), double packed	50 units
MSV508119	Micro Suction (Verhoeven) Brown 86 x 1.05mm OD (19g), double packed	50 units
110mm		
MSV50114	Micro Suction (Verhoeven) Green 110 x 2.1mm OD (14g), double packed	50 units
MSV50117	Micro Suction (Verhoeven) White 110 x 1.47mm OD (17g), double packed	50 units
MSV50119	Micro Suction (Verhoeven) Brown 110 x 1.05mm OD (19g), double packed	50 units
MSV50122	Micro Suction (Verhoeven) Black 110 x 0.7mm OD (22g), double packed	50 units
MSV50127	Micro Suction (Verhoeven) Clear 110 x 0.4mm OD (27g), double packed	50 units
120mm		
MSV50214	Micro Suction (Verhoeven) Green 120 x 2.1mm OD (14g), double packed	50 units
MSV50217	Micro Suction (Verhoeven) White 120 x 1.47mm OD (17g), double packed	50 units
MSV50219	Micro Suction (Verhoeven) Brown 120 x 1.05mm OD (19g), double packed	50 units
MSV50222	Micro Suction (Verhoeven) Black 120 x 0.7mm OD (22g), double packed	50 units
MSV50227	Micro Suction (Verhoeven) Clear 120 x 0.4mm OD (27g), double packed	50 units




UNIVERSAL ADAPTER

- House style regulator with control
- Luer connector suitable for Verhoeven and other micro suction tubes

SUCTION REGULATOR

DSRP5001 Suction Regulator Yellow, double packed 50 units



Tubing and Verhoeven not included



- Single-use removes need to clean a very narrow lumen
- Working length of 77mm



SCHUKNECHT SUCTION HANDLES

SHK2018	Schuknecht Suction Handle 18swg, double packed	20 units
SHK2020	Schuknecht Suction Handle 20swg, double packed	20 units
SHK2022	Schuknecht Suction Handle 22swg, double packed	20 units
SHK2024	Schuknecht Suction Handle 24swg, double packed	20 units



FULL RANGE

SINUS SUCTION HANDLES

SSH2606	Sinus Suction Handle (6fg) 65 degrees, double packed	20 units
SSH2609	Sinus Suction Handle (9fg) 65 degrees, double packed	20 units
SSH2909	Sinus Suction Handle (9fg) 90 degrees, double packed	20 units
SSH2612	Sinus Suction Handle (12fg) 65 degrees, double packed	20 units
SSH2912	Sinus Suction Handle (12fg) 90 degrees, double packed	20 units




90 degrees



65 degrees



For the Tilley Lichtwitz Sinus Punch see p29



- Tear drop control for adjustable suction


TAPERED SUCTION HANDLES

190mm

TDH2195	Tapered Suction Handle 190mm size 5 with tear drop control, double packed	20 units
TDH2196	Tapered Suction Handle 190mm size 6 with tear drop control, double packed	20 units
TDH2197	Tapered Suction Handle 190mm size 7 with tear drop control, double packed	20 units

220mm

TDH2225	Tapered Suction Handle 220mm size 5 with tear drop control, double packed	20 units
TDH2226	Tapered Suction Handle 220mm size 6 with tear drop control, double packed	20 units
TDH2227	Tapered Suction Handle 220mm size 7 with tear drop control, double packed	20 units



Suction Handles



CHOICE OF ROSES

- Ideal when plastic handles are inappropriate
- Detachable rose with three (YST1000) or six holed modified rose (YST1000R)
- Ideal for tonsillectomy



YANKAUER SUCTION HANDLES

YST1000	Yankauer Suction Handle 305mm with rose, double packed	10 units
YST1000R	Yankauer Suction Handle 305mm with modified rose, double packed	10 units



- Available with or without fine ends
- Ergonomic balanced design for ease of use
- Tube rotation suits left or right handed users
- Range for theatre and clinic use



Zoellner Suction Handle without Fine End

ZOELLNER SUCTION HANDLES

SMS4400-S	Zoellner Suction Handle, single packed	50 units
SMS4418-S	Zoellner Suction Handle with 18swg Fine End, single packed	50 units
SMS4400-D	Zoellner Suction Handle, double packed	20 units
SMS4400-DS	Zoellner Suction Handle with stylet, double packed	20 units
SMS4418-D	Zoellner Suction Handle with 18swg Fine End, double packed	20 units

- Robust metal taper suitable for refitting during procedure
- Extra long taper gives a good fit with any suction handle
- Range of sizes from 18 to 26 standard wire gauge



ZOELLNER FINE ENDS

SM4018	Zoellner Fine Ends size 18 swg, double packed	50 units
SM4020	Zoellner Fine Ends size 20 swg, double packed	50 units
SM4022	Zoellner Fine Ends size 22 swg, double packed	50 units
SM4024	Zoellner Fine Ends size 24 swg, double packed	50 units
SM4026	Zoellner Fine Ends size 26 swg, double packed	50 units



- Cream or ointment applicator for use with luer lock syringe

CREAM APPLICATORS

ECA5001-S	Cream Applicator 86mm (6fg), single packed	50 units
-----------	--	----------



Micro Needle Electrodes





COST

- Straight and angled options available
- Compatible with most finger switches and electro-surgical pencils
- High quality stainless steel and Tungsten needle enables fine and precise tissue dissection
- Improved cutting and coagulation performance, whilst minimising thermal spread


STAINLESS STEEL ELECTRODES

MNIS2514	Micro Needle Electrode, Insulated Stainless Steel Tip, 63 x 14mm	25 units
MNIS2514-A	Micro Needle Electrode, Insulated Angled Stainless Steel Tip, 70 x 14mm	25 units
MNIS2516	Micro Needle Electrode, Insulated Stainless Steel Tip, 150 x 16mm	25 units
MNIS2518	Micro Needle Electrode, Insulated Stainless Steel Tip, 70 x 18mm	25 units
MNIS2524	Micro Needle Electrode, Insulated Stainless Steel Tip, 70 x 22mm	25 units
MNES2528	Micro Needle Electrode, Stainless Steel Tip, 70 x 28mm	25 units



TUNGSTEN ELECTRODES

MNET1014	Micro Needle Electrode, Tungsten Tip, 70 x 14mm	10 units
MNET1024	Micro Needle Electrode, Tungsten Tip, 70 x 22mm	10 units
MNIT1014	Micro Needle Electrode, Insulated Tungsten Tip, 63 x 14mm	10 units
MNIT1016	Micro Needle Electrode, Insulated Tungsten Tip, 150 x 16mm	10 units
MNIT1024	Micro Needle Electrode, Insulated Tungsten Tip, 70 x 22mm	10 units



NEW

Endoscopic Ear Range

- New range for growth in ear endoscopy which provide clear visibility at the procedure site
- Overcomes cleaning difficulty and cross contamination with narrow lumen curved tubes
- Available with 3mm & 6mm bends
- Choice of suction control – Regulator or Zoellner
- More options coming soon



Curved Zoellner
Suction Fine End



Zoellner



CURVED ZOELLNER SUCTION FINE ENDS

SME4183	Curved Zoellner Suction Fine End, 3mm Bend Radius 18 swg, double packed	20 units
SME4186	Curved Zoellner Suction Fine End, 6mm Bend Radius 18 swg, double packed	20 units

CURVED ENDOSCOPIC SUCTION TUBES


CES2183	Curved Endoscopic Suction Tube, 3mm Bend Radius 18 swg, double packed	20 units
CES2186	Curved Endoscopic Suction Tube, 6mm Bend Radius 18swg, double packed	20 units

Kerrison Rongeur

- 200mm shaft length with 40 degree up bite
- Well balanced stainless steel construction

KERRISON RONGEUR

KRU1101	Kerrison Rongeur Punch 1mm jaw, length 200mm with 40 degree up bite, double packed	10 units
KRU1201	Kerrison Rongeur Punch 2mm jaw, length 200mm with 40 degree up bite, double packed	10 units
KRU1301	Kerrison Rongeur Punch 3mm jaw, length 200mm with 40 degree up bite, double packed	10 units
KRU1401	Kerrison Rongeur Punch 4mm jaw, length 200mm with 40 degree up bite, double packed	10 units
KRU1501	Kerrison Rongeur Punch 5mm jaw, length 200mm with 40 degree up bite, double packed	10 units



Scope Warmer



TIME



COST

THEATRE EFFICIENCY

- Simple preheating to eliminate condensation
- Secures the surgeon's vision
- Warms full scope length to 37°C within 90 seconds, staying above body temperature for 120 minutes
- Long-acting to avoid disruption or repeated defogging

SCOPE WARMER

31SW10	ClearView Anti-Fog Scope Warmer, sterile packed for single-use	10 units
31SW30	ClearView Anti-Fog Scope Warmer, sterile packed for single-use	30 units



Instructional poster available
Contact us on +44 (0)1792 797910

“Supporting visibility, quality and safety in your supply chain”**

Our instruments' packaging features the latest machine readable GS1-standard barcoded labels for product traceability.

The newly incorporated labelling technology is a timely development with an increasing number of health providers worldwide adopting an 'eProcurement' strategy and will provide machine readable product information on all labels to meet expected EU and other regulators demands.

Associated benefits include:

- Improve patient care and reduce costs
- Improve stock control
- One global language
- Accurately track equipment


“Unique device identification (UDI) compliance represents a significant development for our customers. Automation of information flow will have influences across the whole purchasing and supply process, providing the benefits of superior instrument traceability and improved patient safety.”

“We are working closely with Kodit UDI software to provide our labelling solutions to meet the latest international regulations and we serve as a reference site for other manufacturers to learn from.”

Andrew Davidson
(Managing Director of DTR Medical)

For more information on our new product labelling, please contact Customer Services on +44 (0) 1792 797910 or info@dtrmedical.com

** - GS1 – Global Traceability Standards





medical
sales &
service

A: 38 Hector Street OSBORNE PARK WA 6017

T: 08 9361 0843 | **F:** 08 9470 3647

E: medsales@medsales.com.au | **W:** www.medsales.com.au