



medical  
sales &  
service



## BritePro Solo Combination Sets

### Complete

Sterile packed with fitted batteries and changeable blade, BritePro Solo is ready to use straight out of the box. Comes complete with handle, blade and batteries.

### Reliable

With the advantage of being new every time, BritePro Solo can be tested in pack without compromising sterility.

### Performance

Unique textured grip offers excellent control and feel, while the high intensity light provides excellent visualisation.

### Ultra-Low Profile Blade

Macintosh blade design reduces the risk of dental trauma and improves visualization of the vocal cords, even during difficult intubation.

ORDER CODE	DESCRIPTION	UOM
PAFLC-040-331	Laryngoscope BritePro Solo Disposable Combo Set F/O MAC SIZE 1	EACH
PAFLC-040-332	Laryngoscope BritePro Solo Disposable Combo Set F/O MAC SIZE 2	EACH
PAFLC-040-333	Laryngoscope BritePro Solo Disposable Combo Set F/O MAC SIZE 3	EACH
PAFLC-040-334	Laryngoscope BritePro Solo Disposable Combo Set F/O MAC SIZE 4	EACH
PAFLC-040-336	Laryngoscope BritePro Solo Disposable Combo Set F/O MAC SIZE 3 Curve Blade	EACH
PAFLC-040-340N	Laryngoscope BritePro Solo Disposable Combo Set F/O MILLER SIZE 00	EACH
PAFLC-040-340	Laryngoscope BritePro Solo Disposable Combo Set F/O MILLER SIZE 0	EACH
PAFLC-040-341	Laryngoscope BritePro Solo Disposable Combo Set F/O MILLER SIZE 1	EACH
PAFLC-040-342	Laryngoscope BritePro Solo Disposable Combo Set F/O MILLER SIZE 2	EACH
PAFLC-040-343	Laryngoscope BritePro Solo Disposable Combo Set F/O MILLER SIZE 3	EACH
PAFLC-040-344	Laryngoscope BritePro Solo Disposable Combo Set F/O MILLER SIZE 4	EACH
PAFLC-040-350N	Laryngoscope BritePro Solo Disposable MINI Combo Set F/O MILLER SIZE 00	EACH
PAFLC-040-350	Laryngoscope BritePro Solo Disposable MINI Combo Set F/O MILLER SIZE 0	EACH
PAFLC-040-351	Laryngoscope BritePro Solo Disposable MINI Combo Set F/O MILLER SIZE 1	EACH

# New BriteBlade Pro Blades

BriteBlade Pro is a metal disposable fibre optic blade designed for single patient use which eliminates risk of cross contamination from both blade and handle. It has a unique hinge design that prevents the folded

### Enhanced Vision

Bright, focused high intensity illumination

### Eliminate Risk of Cross-Infection

Single-use individually packed pre-sterilized blade.

### Eliminate Risk of Cross-Contamination

Blade stop design in folded position prevents the blade from touching and transferring contaminants to the handle.<sup>1</sup>

Fits all Green Specification Handles



### Confidence & Control

100% metal construction matches the proven re-usable blade characteristics.

### Reduced Risk of Damaging Patient Teeth

Low profile design on both Mac & Miller laryngoscope blades.

### Wide Choice at a Low Cost

Full range of Mac & Miller Blades with a lower overall cost than re-usable purchase, maintenance, cleaning and sterilization costs.

ORDER CODE	DESCRIPTION	UOM
PAFLC-040-711	Laryngoscope BriteBlade Pro Disposable Blade F/O MAC SIZE 1	BOX 10
PAFLC-040-712	Laryngoscope BriteBlade Pro Disposable Blade F/O MAC SIZE 2	BOX 10
PAFLC-040-713	Laryngoscope BriteBlade Pro Disposable Blade F/O MAC SIZE 3	BOX 10
PAFLC-040-714	Laryngoscope BriteBlade Pro Disposable Blade F/O MAC SIZE 4	BOX 10
PAFLC-040-720N	Laryngoscope BriteBlade Pro Disposable Blade F/O MILLER SIZE 00	BOX 10
PAFLC-040-720	Laryngoscope BriteBlade Pro Disposable Blade F/O MILLER SIZE 0	BOX 10
PAFLC-040-721	Laryngoscope BriteBlade Pro Disposable Blade F/O MILLER SIZE 1	BOX 10
PAFLC-040-722	Laryngoscope BriteBlade Pro Disposable Blade F/O MILLER SIZE 2	BOX 10
PAFLC-040-723	Laryngoscope BriteBlade Pro Disposable Blade F/O MILLER SIZE 3	BOX 10
PAFLC-040-724	Laryngoscope BriteBlade Pro Disposable Blade F/O MILLER SIZE 4	BOX 10

# New BritePro Max LED Handles

### Energy Efficient

Delivers over 25,000 hours of light without the need to replace the bulb providing cost savings.

### Compatible with all Fibre Optic Blades

BritePro Max is compatible with all fibre optic blades conforming to ISO EN 7376 which includes, the BriteBlade Pro, EVO and reusable fibre optic blades within the Flexicare range.





### Green System Identification

BritePro Max has the universally accepted colour identification green band to signify that the handle will fit fibre optic blade systems.

### Lower Battery Cost

The energy efficiency of LED technology uses less power from the batteries than conventional halogen or krypton bulbs which reduces the need for frequent battery replacement.

True White Light



ORDER CODE	DESCRIPTION	UOM
PAFLC-040-813	BritePro Max LED Fibre Optic Laryngoscope Handle Standard Reusable	Each
PAFLC-040-814	BritePro MAX LED Fibre Optic Laryngoscope Handle Paediatric Mini Reusable	Each
PAFLC-040-815	BritePro MAX LED Fibre Optic Laryngoscope Handle Stubby Reusable	Each