

FITLEGS™ Anti-Embolism Stockings

FITLEGS™ AES (anti-embolism stockings) help to prevent the formation of Deep Vein Thrombosis (DVT) in the legs by applying clinically proven graduated compression.

FITLEGS™ AES are designed to be simple to fit and apply, which reduces sizing errors and nursing time.

IN STOCK NOW!




KEY FEATURES

- DUOFIT™ Technology provides a two-way stretch for easier application.
- Single ankle measurement reduces nursing time when measuring patients.
- Improved inspection hole design and placement for superior comfort.
- Available in below knee and thigh length.
- Four sizes within each length reduces inventory.
- Colour coded packaging reduces time identifying correct size.
- Darker colour (Teal) improves patient compliance.
- Anti-bacterial fabric for odour management.



BENEFITS OF FITLEGS™

The compression within the garment is highest at the ankle and reduces further up the leg. The pressure from the Fitlegs™ stockings forces blood back towards the heart.



ORDERING INFORMATION



FITLEGS AES BELOW KNEE

Size	Ankle Measurement	Order Code
Small	19-23cm	PAGN-401000
Medium	23-26cm	PAGN-401001
Large	26-29cm	PAGN-401002
Extra Large	29-33cm	PAGN-401003



FITLEGS AES THIGH LENGTH

Size	Ankle Measurement	Thigh Measurement	Order Code
Small	19-23cm	65cm	PAGN-401004
Medium	23-26cm	70cm	PAGN-401005
Large	26-29cm	78cm	PAGN-401006
Extra Large	29-33cm	82cm	PAGN-401007