Reynard IV access and clinical skin cleansing. Our extensive range means you're prepped for anything.

The perfect solution for every clinical skin cleansing situation.

Our broad range includes innovative packaging with 'easy-peel' and 'tear open' options. This, combined with unsurpassed value, ensure our prep pads and swabsticks are used daily throughout Australia. Providing every healthcare system provider in the country with the highest quality, TGA registered products.





Reynard[™] IV access & clinical skin cleansing.

The organisms may be micro, but the impact can be big.

Surgical site infections (SSIs) are a major challenge for patients and health sector resources. Effective cleansing and disinfection using appropriate Reynard $^{\text{TM}}$ prep pads and swabsticks can help minimise the risk of surgical site infections and improve overall patient outcomes.

- SSIs negatively impact patient physical and mental health.
- Increased patient morbidity, mortality, and loss of earning during recovery are some of the indirect costs associated with infection.
- Intangible costs may also be incurred by the patient, such as pain and anxiety. In addition, patients may experience delayed wound healing and more susceptible to secondary complications, such a bacteraemia.
- Distress may also be caused to the patient and family member if the patient is absent from home and work for a prolonged period.
- Accordingly, prolonged hospitalisation and increased morbidity as a result of developing an SSI have been shown to negatively impact patient health-related quality of life.
- Increase in healthcare costs per patient and additional stress on healthcare systems and staff.

The perfect solution for every problem.

Reynard[™] provides the broadest range of prep pads and swabsticks in a variety of formulations and configurations to suit your specific clinical need and appropriate application.

- 0.5% chlorhexidine digluconate (CHG) and 70% isopropyl alcohol (IPA)
- 2% chlorhexidine digluconate (CHG) and 70% isopropyl alcohol (IPA)
- 70% isopropyl alcohol (IPA) sterile
- 0.5% chlorhexidine digluconate (CHG)
- 2% chlorhexidine digluconate (CHG)



Reynard[™] Health Supplies.

Reynard[™] Health Supplies was established in New Zealand in 1999 and today is the number one supplier of surface cleaning & disinfection wipes, specialist patient care wipes and clinical skin cleansing products to hospitals in New Zealand and Australia.

Reynard[™] is owned by a dedicated group of Australian and New Zealand shareholders who are all actively and directly involved in the day-to-day operation and development of the business globally.

Used extensively in the private sector and every public healthcare system in Australia.

Reynard[™] is reassurance.

At Reynard[™], we strive to provide products and support services that provide you with the reassurance you need to be able to focus on optimum patient outcomes.

This reassurance stems from:

- Quality & Range a wide range of quality products to suit every need
- Cost Savings cost effective and pricing stability provide real value
- Reliability a reliable, proven supply chain that has demonstrated reliability and consistency throughout the global pandemic
- Support Dedicated technical, educational and compliance support services
- Research & Development focus drives continuous product development and a pipeline of innovative solutions
- Flexibility the capability and desire to design new products and supply chain options





Our extensive range means you're prepped for anything.

| Size and type of product | | | | | | |
|--------------------------|-----------|---------|----------------|-----------|-----------|--|
| | Prep pads | | | | | |
| | | | Substrate only | | | |
| Formula | 6 x 3cm | 6 x 6cm | 8.5 x 4.5cm | 20 x 16cm | Easy Peel | |
| 70% IPA only | RHS451G | RHS417 | RHS438 | | | |
| 2% CHG + 70% IPA | RHS461 | RHS401 | RHS437 | RHS422 | RHS422EP | |
| 0.5% CHG + 70% IPA | | | | | | |
| 2% CHG only | | | RHS421 | | | |
| 0.5% CHG only | | RHS479 | | | | |

| Size and type of product | | | | | | | | |
|--------------------------|------------|----------------|-----------|----------|-----------|--------|----------|-----------|
| | Swabsticks | | | | | | | |
| | Cotton he | ead | | | Foam Head | | | |
| Formula | Regular | Triple Pack | Easy Peel | Tint | Regular | Maxi | Reg Tint | Maxi Tint |
| 70% IPA only | | | | | | | | |
| 2% CHG + 70% IPA | RHS408 | RHS408T | RHS408EP | RHS408TR | RHS473 | RHS474 | | RHS472 |
| 0.5% CHG + 70% IPA | | | | | RHS470 | | RHS470TR | RHS471 |
| 2% CHG only | | | | | RHS475 | | | |
| 0.5% CHG only | | | | | RHS477 | RHS478 | | |

Skin cleansing considerations.

| Product formulation | Recommended product application |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 0.5% CHG + 70% IPA | General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potentially allergies, skin reactivity and risk. |
| 2% CHG + 70% IPA | General use where sensitivity etc. are not a concern. This solution provides broad capacity for most patients. |
| 70% IPA Sterile | General use, also suitable for paediatrics and areas where use of CHG is not preferred due to potentially allergies, skin reactivity and risk |
| 0.5% CHG only | General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potentially allergies, skin reactivity and risk. Alcohol (IPA) free, appropriate where use and application of IPA is not suitable. |
| 2% CHG only | General use where CHG concentration sensitivity is not a concern. Alcohol (IPA) free, appropriate where use and application of IPA is not suitable. |

All Reynard[™] clinical skin cleansing products are suitable for a wide range of skin cleansing requirements including but not limited to:

- Blood culture collections
- Central venous catheters insertion
- General skin preparation, pre-procedure, injection, or incision
- Insertion of percutaneous devices
- Percutaneous device site maintenance
- Peripheral intravascular catheters insertion
- Site maintenance of central lines

Note: Only use on intact skin before any invasive procedures

| Directions | Product should be applied with gentle friction using a repeated back-and-forth motion for manufacturer-recommended application time. |
|------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Dry time | Allow product to aur dry completely; dry time can be affected by volume of product, body sit, presence or absence of hair, humidity etc. |

Reynard[™] prep pad and Swabstick product application

Instructions for use:

- 1. Check allergy adverse reaction status with your patient.
- 2. Thoroughly clean and treat the chosen area with the prep pad using a repeated backand-forth scrubbing motion for 20 to 30 seconds.
- Use the necessary number of prep pads dependent on the clinical procedure.
- 4. Allow the skin area to fully air dry to be most effective against bacteria growth. Allow approximately 30 seconds for a dry site and up to 2 minutes for a moist site.
- Discard of prep pad after single use, do not flush.
- 6. Dispose of prep pad in adherence with your facility infection control guidelines.

Is this product flammable?

- Refer to manufacturer's instructions for use, and any associated hazard warnings
- Follow AORN guidelines for perioperative practice



Prep pads

CHLORHEXIDINE & ALCOHOL

| | Alcohol and chorhesidine Prep pads Typed T | Alcohol and chlorheridine Prep pads 15-12-2 Types Types | Model and chlorhostifice Prop page 819192 | ************************************** |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------|----------------------------------------|
| Product code | RHS401 | RHS422 | RHS422EP | RHS437 |
| Swabs per box | 200 | 100 | 50 | 100 |
| Boxes per carton | 25 | 10 | 20 | 10 |
| Wipe size | 6 x 6cm | 20 x 16cm | 20 x 16cm | 8.5 x 4.5cm |
| Solution | 70% IPA & | 70% IPA & | 70% IPA & | 70% IPA & |
| | 2% CHG | 2% CHG | 2% CHG | 2% CHG |
| Volume | 0.7ml | 4.5ml | 5ml | 1.6ml |
| GSM | 45gsm | 40gsm | 40gsm | 100gsm |

| CH | - | | | 11 - |
|----|---|-------|---|-------------|
| | | 54 51 | - | VI - |
| | | | | |


| | Charleston Prep pads For part of the par | Chlorheddine Prep pads For presented and abuptals as Colorheddine Prep pads For presented and buptals as Colorheddine Prep pads For presented and buptals as Reynard Reynard Reynard Rosstry | Acond a Chlorhesidite Prep pads Story Figure Figure |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product code | RHS421 | RHS479 | RHS461 |
| Swabs per box | 100 | 200 | 200 |
| Boxes per carton | 10 | 25 | 25 |
| Wipe size | 8.5 x 4.5cm | 6 x 6cm | 6 x 3cm |
| Solution | 2% CHG | 0.5% CHG | 70% IPA & 2% CHG |
| Volume | 0.8ml | 0.7ml | 0.5ml |
| GSM | 45gsm | 45gsm | 45gsm |

ALCOHOL

| | Alcohol STERK Prep pads Pr | Premium Alcohol Prep pads Prep pads Prepads Pr | Alcohol Prep pads Framework and and an |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------|
| Product code | RHS417 | RHS438 | RHS451G |
| Swabs per box | 200 | 100 | 200 |
| Boxes per carton | 25 | 10 | 25 |
| Wipe size | 6 x 6cm | 8.5 x 4.5cm | 6 x 3cm |
| Solution | 70% IPA | 70% IPA | 70% IPA |
| Volume | 0.7ml | 1.6ml | 0.5ml |
| GSM | 45gsm | 100gsm | 45gsm |







22cm

| Size | Product code | Coverage area |
|---------|------------------------|--------------------|
| Small | RHS451G, RHS461 | 42cm² |
| Medium | RHS401, RHS417, RHS479 | 100cm ² |
| Large | RHS421, RHS437, RHS438 | 225cm ² |
| X-Large | RHS422, RHS422EP | 484cm² |

*Note: coverage area may vary in accordance with variations in atmospheric conditions and condition of the skin. Tests we conducted on hairless male skin at 21 degrees Celsius. The skin area was wetted sufficiently to ensure it was wet for a minimum 30 seconds before allowing to dry. In each test enough liquid was left in the product to allow up to a further 20% area to be covered. Diagram is not to scale.

Swabsticks

COTTON HEAD









| | Reynard assets | Reynard scent | Reynard mesers | Reynard wester |
|-------------------|------------------|------------------|------------------|------------------|
| Product code | RHS408 | RHS408T | RHS408TR | RHS408EP |
| Sticks per sachet | 1 | 3 | 1 | 1 |
| Sachets per box | 50 | 25 | 50 | 50 |
| Boxes per carton | 10 | 10 | 10 | 10 |
| Stick type | Cotton | Cotton | Cotton | Cotton |
| Solution | 70% IPA & 2% CHG |
| Tint | No | No | Yes | No |
| Volume | 1.6ml | 4.5ml | 1.6ml | 1.6ml |

| F | O, | ¥М | н | EΑ | D |
|---|----|----|---|----|---|
| | | | | | |











| | Royand work | V Later Control | ₩ Reynard work | | |
|------------------|-----------------------|-----------------------|---------------------|--------|----------|
| Product code | RHS470 | RHS470TR | RHS473 | RHS475 | RHS477 |
| Sticks per box | 25 | 25 | 25 | 25 | 25 |
| Boxes per carton | 10 | 10 | 10 | 10 | 10 |
| Stick type | Foam | Foam | Foam | Foam | Foam |
| Solution | 70% IPA & 0.5% CHG | 70% IPA & 0.5% CHG | 70% IPA & 2% CHG | 2% CHG | 0.5% CHG |
| Tint | No | Yes | No | No | No |
| volume | 2.5ml | 2.5ml | 2.5ml | 2.5ml | 2.5ml |

LARGE FOAM HEAD










| Product code | RHS471 | RHS472 | RHS474 | RHS478 |
|------------------|--------------------|------------------|------------------|------------|
| Stick per box | 25 | 25 | 25 | 25 |
| Boxes per carton | 10 | 10 | 10 | 10 |
| Stick type | Large Foam | Large Foam | Large Foam | Large Foam |
| Solution | 70% IPA & 0.5% CHG | 70% IPA & 2% CHG | 70% IPA & 2% CHG | 0.5% CHG |
| Tint | Yes | Yes | No | No |
| Volume | 5ml | 5ml | 5ml | 5ml |



22cm

| Size | Product code | Coverage area |
|---------------------------|------------------------------------------|-------------------|
| Regular cotton swabsticks | RHS408, RHS408EP, RHS408T, RHS408TR | 72cm ² |
| Regular foam swabsticks | RHS470, RHS470TR, RHS473, RHS475, RHS477 | 144cm² |
| Large foam swabsticks | RHS471, RHS472, RHS474, RHS478 | 484cm² |

^{*}Note: coverage area may vary in accordance with variations in atmospheric conditions and condition of the skin. Tests we conducted on hairless male skin at 21 degrees Celsius. The skin area was wetted sufficiently to ensure it was wet for a minimum 30 seconds before allowing to dry. In each test enough liquid was left in the product to allow up to a further 20% area to be covered. Diagram is not to scale.

Reynard focus on killing germs so you can focus on saving lives.

Germs.

They can proliferate in the cleanest of places. A constant bugbear, literally, to be eradicated, first time, every time. Reynard™ constantly adopt new technology, new research to develop new products in all our ranges. Working to protect patients, and everyone and everything they come into contact with. An ever-growing Reynard™ product range means you can choose a product that meets every operational need, whether cleaning, disinfecting, skin cleansing or patient care. Many applications developed in consultation with health professionals, in response to urgent clinical needs. So you can concentrate on doing what you do best. Caring for patients. Visit www.reynardhealth.com.au

Reynard is reassurance.



More from Reynard[™]...

Personal Care



Our diverse range of personal care products provide gentle and effective cleansing while minimising cross-infection risks.

Cleaning & disinfection



Developed with convenience in mind, this range provides fast and effective solutions for every cleaning need.

P.P.E.



Promotes infection control and works towards minimising the spread of bacteria, dust and pathogens.

Dry Wipes



Disposable, practical and absorbent, Reynard[™] dry wipes are some of the softest and strongest wipes on the market.

Hand Hygiene



Killing 99.9% of germs while leaving hands soft and refreshed. Available with or without chlorhexidine.

Baby Care



Chemical and fragrance-free, our baby care range is designed to be ultra-soft and gentle on little ones.

