




medical
sales &
service




Leadercath Arterial

Setting High Standards to Reduce Line Infection




Safe
Simple
Successful



Product Information



Leadercath Arterial



Insertion of arterial lines into small vessels can be technically challenging, particularly if the patient is hypotensive or peripherally vasoconstricted. Leadercath Arterial offers the 'classic' Seldinger method. This technique has demonstrated a significantly higher rate of success than the direct threading method. S. C. Beards, et al recommend the use of a 'classical' Seldinger technique (method B) for arterial line insertion.⁽¹⁾



"My technique was the result of an idea... Vygon have developed that idea with skill and innovation, resulting in an excellent product which I highly recommend."

S. Seldinger
Sven-Ivar Seldinger

Features and Benefits

Classic Seldinger insertion

for higher success rate.
93% success rate compared to 76% with cannula.

Anti-kink collar

maintains patency.

Winged hub


for secure fixation.

Range of sizes

for radial and femoral placement to suit all clinical requirements.

Customisation

to match your clinical requirements when placing arterial lines and maintain maximum sterile barrier precautions.



Ordering Information

Product Codes		Description	Catheter		Needle		Guidewire		Unit of Sale
Vygon	NPC		Size (G)	Length (cm)	Size (G)	Length (cm)	Size (mm)	Length (cm)	
0115.09	FSQ049	Leadercath Arterial PE (red hub)	20	8	20	3.2	0.50	20	20
0115.11	FSQ004	Leadercath Arterial PE (red hub)	18	10	19	5.4	0.70	30	20
1211.04	FSQ322	Leaderflex	20	4	20	3.2	0.50	20	20
1211.06	FSQ323	Leaderflex	20	6	20	3.2	0.50	20	20
1211.08	FSQ324	Leaderflex	20	8	20	3.2	0.50	20	20
1211.12	FSQ325	Leaderflex	20	12	20	3.2	0.50	30	20
1212.04	FSQ326	Leaderflex with Integrated Extension	22	4	22	4.2	0.46	23	20
1212.06	FSQ327	Leaderflex with Integrated Extension	22	6	22	4.2	0.46	23	20
1212.08	FSQ328	Leaderflex with Integrated Extension	22	8	22	4.2	0.46	26	20
1213.06	FSQ331	Leaderflex (red hub)	18	6	19	5.4	0.70	20	20
1213.10	FSQ332	Leaderflex (red hub)	18	10	19	5.4	0.70	30	20
80.1170.25A		Spare Guidewire for 115.09					0.50	20	20

Product Information



Radial Arterial Pack 20G 8cm



The Vygon radial arterial pack is setting the standard for arterial line placement. The components of the pack have been specifically designed and selected by clinicians who insert arterial lines on a daily basis. Having all of the required components in one pack reinforces the use of maximum barrier precautions, thereby increasing the likelihood that the highest standard of asepsis is maintained by all clinicians.

Features and Benefits

Maximum barrier protection

helping reduce the risk of infection to the patient.

Everything you need in one pack

reduced set-up time, and reduced drop contamination.⁽²⁾

Red hub


for easy arterial line identification.

Classic Seldinger insertion

for higher success rate.⁽¹⁾

Kit Contents

1 x Leadercath Arterial catheter	1 x Hypodermic needle 25G 5/8"
1 x 20G Introducer needle	5 x Swabs 10 x 10cm
1 x Guidewire	1 x Dressing 10 x 12cm
1 x 2.5ml Luer-slip syringe	1 x Chloraprep® Frepp 1.5ml
1 x 5ml Luer-slip syringe	1 x Arterial peelable drape 45 x 70cm



Leadercath Arterial catheter

Ordering Information

Vygon Code	NPC Code	Product Description	Unit of Sale
OP0115802PH	FSQ935	Radial arterial pack 20G 8cm with Tegaderm™ dressing	40
OP0115803PH	FSU312	Radial arterial pack 20G 8cm with IV3000 dressing	40