

Precision that lasts a lifetime



ri-san® aneroid sphygmomanometer

The ri-san® high-quality, palm-style sphygmomanometer is the embodiment of ergonomic design. The air release valve is conveniently regulated by finger-tip control. Also optionally available in a latex-free version (when ordering enter LF for latex-free).

- **Specially tempered copper-beryllium diaphragm, almost non-aging**
- **Precision air release valve with fine adjustment and wear-free**
 - Latex-free bulb for pumping up the cuff
 - Tube connection at the top of the pressure gage for unhindered blood pressure measurement
 - Suitable for right or left handed users
- **Microfilter protects air release valve and measuring system**
- **Diaphragm can be pressurized up to 600 mmHg**
- **Made in Germany**
 - Easily legible scale Ø 64 from 0 to 300 mmHg
 - Maximum error tolerance of +/- 3 mmHg
 - In vinyl zippered case
 - Comprehensive selection of different sizes and versions of cuffs

- ri-san® with cotton bandage cuff for adults
- ri-san® with cotton hook cuff for adults
- ri-san® with calibrated nylon velcro cuff for adults

slate gray	blue	saffron	green
No.: 1510	No.: 1515	No.: 1520	No.: 1525
No.: 1511	No.: 1516	No.: 1521	No.: 1526
No.: 1512	No.: 1517	No.: 1522	No.: 1527