



# FarrowWrap® 4000

Patient Name \_\_\_\_\_ DOB \_\_\_\_\_ Date \_\_\_\_\_

Address \_\_\_\_\_ Gender M  F

City/State/Zip \_\_\_\_\_ Diagnosis \_\_\_\_\_


Doctor / Address \_\_\_\_\_ PO# \_\_\_\_\_

City/State/Zip \_\_\_\_\_

Original Order <input type="checkbox"/> Reorder w Changes <input type="checkbox"/>	
Exact Reorder <input type="checkbox"/>	

### 3 Easy Measuring Steps:

1. Length (B-D) - Ankle to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) - Calf (widest)
3. Circumference (B) - Ankle (smallest part)



Measurements		
	Left	Right
cC		
cB		
/B-D		

### JOBST® FarrowWrap® 4000

### JOBST® FarrowWrap® 4000 Size Chart

	XS	S	M	L
<b>C</b> Mid-Calf	25-38 cm	32-48 cm	37-56 cm	43-67 cm
<b>B</b> Ankle	17-23 cm	19-28 cm	25-37 cm	29-43 cm
<b>B-D</b> Regular	30-33 cm	30-33 cm	30-33 cm	30-33 cm
<b>B-D</b> Tall	34-37 cm	34-37 cm	34-37 cm	34-37 cm

Ref No.	Size	Length	Color	Quantity
7666100	XS	Regular	Tan	
7666101	S	Regular	Tan	
7666102	M	Regular	Tan	
7666103	L	Regular	Tan	
7666108	XS	Regular	Black	
7666109	S	Regular	Black	
7666110	M	Regular	Black	
7666111	L	Regular	Black	
7666104	XS	Tall	Tan	
7666105	S	Tall	Tan	
7666106	M	Tall	Tan	
7666107	L	Tall	Tan	
7666112	XS	Tall	Black	
7666113	S	Tall	Black	
7666114	M	Tall	Black	
7666115	L	Tall	Black	