



# FarrowWrap® Basic

Patient Name / BSN File # \_\_\_\_\_ DOB \_\_\_\_\_ Date \_\_\_\_\_

Address \_\_\_\_\_ Gender M  F

City/State/Zip \_\_\_\_\_

Diagnosis \_\_\_\_\_


Doctor / Address \_\_\_\_\_

City/State/Zip \_\_\_\_\_

<b>PO#</b>
Original Order <input type="checkbox"/> Reorder w Changes <input type="checkbox"/>
Exact Reorder <input type="checkbox"/>

### 3 Easy Measuring Steps:

1. Length (A-D) - Floor to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) - Calf (widest)
3. Circumference (B) - Ankle (smallest part)



Measurements		
	Left	Right
cC		
cB		
cA-D		

### JOBST® FarrowWrap® Basic (30-40 mmHg\*)

Ref No.	Size	Length	Quantity
7666000	XS	Regular	
7666001	S	Regular	
7666002	M	Regular	
7666003	L	Regular	
7666004	XL	Regular	
7666005	XS	Tall	
7666006	S	Tall	
7666007	M	Tall	
7666008	L	Tall	
7666009	XL	Tall	

### JOBST® FarrowWrap® Basic Size Chart

	XS	S	M	L	XL
<b>C</b> Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
<b>B</b> Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
<b>A-D</b> Regular	33-37 cm	35-39 cm	37-41 cm	39-43 cm	39-43 cm
<b>A-D</b> Tall	38-41 cm	40-43 cm	42-45 cm	44-47 cm	44-47 cm

\* The mean compression for an average size ankle  
Please refer to the product label and / or package insert for full instructions on the safe use of these products.



JOBST®,  
an Essity brand



/JOBSTUSA



@JOBST\_USA



@JOBSTforUSA



jobst-usa.com

BSN medical Inc., an Essity company • 5825 Carnegie Blvd., Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 • To Order Online: <https://order.jobst.com/us>  
To order toll-free: Tel. (+1) 877 417 5187 • [www.bsneshop.com](http://www.bsneshop.com)