

15301

## **FluidGard® 160 Anti-Fog Procedure Mask**


- Complies with ASTM F2100 Level 3 Requirements
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



### Specifications

Product	FluidGard® 160 Anti-Fog Procedure Mask
Color	Blue Diamond
Fluid Resistance	160mm Hg
Filtration - BFE	≥99.8
Filtration - PFE	≥99
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10/Boxes, 500/Cs

This product is NOT made with natural rubber latex



### Complies with ASTM F2100

**Fluid Resistance (mmHg): 160**  
**BFE (%) @ 3.0 µm: ≥ 98%**  
**PFE (%) @ 0.1 µm: ≥ 98%**  
**ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5**  
**Flame Spread: Class 1**