


ACE - FL

Single Lumen Ovum Pickup Needle with Flushing Line




allwin®
Medical Devices
www.allwinmedical.com

ACE - FL

Single Lumen Ovum Pickup Needle with Flushing Line

The ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.




Tri-facet geometry
with highly
echogenic tip



Multiple operation
functionality
by using 3 way
stopcock



Luer lock
connector



Silicone bung
for tight fit

Order No.	Description	Needle Length (cm)	Needle Gauge	Qty. / Box
A-OVM-SL16/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	16	10
A-OVM-SL17/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	17	10
A-OVM-SL18/35-FL	ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line	35	18	10

*All products may not be available in all countries. Please contact our sales department for details.



US Sales

✉ sales.usa@allwinmedical.com
☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com
☎ +91 99870-84563

Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com
☎ +91 95949-81581