


# FLEXI IUI




## Intra Uterine Insemination Catheter

**allwin**<sup>®</sup>  
Medical Devices  
[www.allwinmedical.com](http://www.allwinmedical.com)

# FLEXI IUI

## Intra Uterine Insemination Catheter

The FLEXI Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.




Stylet




Centimetre marking provided



Soft uterine catheter




Firm Support



Shapable outer sheath

Order No.	Description	Qty. / BOX
F-IUI	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	25
F-IUI-O	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip)	25
STY18	Stylet for FLEXI IUI Intra Uterine Insemination Catheter	25

\*All products may not be available in all countries. Please contact our sales department for details.



Close Tip



Open Tip

### US Sales

✉ sales.usa@allwinmedical.com  
☎ +1 714-572-1709

### India Sales

✉ sales.india@allwinmedical.com  
☎ +91 99870-84563

### Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com  
☎ +91 95949-81581