

The MegaVID Wi-Fi 12MP camera is an incredible addition to our camera line. Operating at 5GHz frequency for maximum speed (3840x2160@30fps) the MegaVID can be mounted to nearly any microscope that supports trinocular C-mount and connects to any iOS or Android device through the LW Scientific app. Utilize the MegaVID Wi-Fi to its fullest capabilities using the included software disc to connect to any desktop. The MegaVID Wi-Fi has many applications but is best suited for teaching, training, and educating clients and students.

SPECIFICATIONS

- Sony IMX412 CMOS Image Sensor Compatible with Nikon & Olympus
- 12MP Camera 3840 x 2160 30fps
- 5G Wi-Fi
- Real time auto and manual white balance adjustment
- DC 5V 2A Power Supply

FEATURES

- · Live stream video clips on any device
- High-res images in BMP, JPG, and PNG formats
- LW Scientific app available on iOS and Android or desktop through the included software disc

REQUIREMENTS

System Requirements:

- iOS11 or higher (Mac OS supported through the LW Scientific app)
- Android 3.0 or higher
- Microsoft Windows 7 or higher (32 & 64 bit)
- Compatible C-Mount Trinocular pipe



SKU MVC-12MP-WIFI



PACKAGE INCLUDES

- [1] MegaVID Wi-FI camera
- [1] HDMI Cable
- [1] Ethernet Cable
- [1] Software Disk
- [1] Power Adapter

