

PRO EMBRYO TRANS



Embryo Transfer Catheter

PRO EMBRYO TRANS

Embryo Transfer Catheter

The PRO EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

PRO EMBRYO TRANS Embryo Transfer Catheter is also available with compatible stylet or soft obturator which assist the guide catheter in gaining uterine access in the embryo transfer procedure.




Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter




Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter



Soft inner catheter and smooth tip



Silicone depth positioner



Centimetre marking



Orientation arrow on hub

Order No.	Description	Qty. / Box
PEB-ETC	PRO EMBRYO TRANS Embryo Transfer Catheter	10
PEB-ETC-SB	PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
PEB-ETC-STY	PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10
STY-PE	Stylet for Pro Embryo Trans / Pro Echo Trans / Twinkle Pro Embryo Trans Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.



US Sales

✉ sales.usa@allwinmedical.com
☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com
☎ +91 99870-84563

Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com
☎ +91 95949-81581