Revelation III

The Revelation III professional-grade microscope is a best-seller into physician clinics, as well as universities and medical schools. The Revelation III comes with LED illumination for daylight color, low heat, low power consumption, and super-long bulb life (30,000 hours). Equipped for performance, its features include titanium-finished DIN achromatic or DIN Plan optics and a 30-year anti-fungal coating. The Revelation III is the choice for superior performance at a great value!

HEAD

Binocular (Seidentopf)
Inclined 30°, rotates 360°
Diopter adjustment
10X/18 wide field eyepieces (10X/20HP available)
Monocular, trinocular, dual binocular and
45° sliding heads available
Interpupillary distance range:
Binoc head: 48-75mm
Trinoc head: 55-75mm



NOSEPIECE

Quadruple forward-facing nosepiece Multiple ball bearing mounted

OBJECTIVES

DIN Achromatic or DIN Plan 30 year anti-fungal coating 4X, 10X, 40X, 100XR (oil) Parfocal, parcentric, color-coded

STAGE

Mechanical stage (140mm x 130mm) Coaxial drive controls Range of traverse is 78mm x 50mm Stage locking lever Slow-close hydraulic slide finger

FOCUS

Coarse adjustment: range of 30mm Fine adjustment: graduation of 2µm Tension control knob

ILLUMINATION

Moveable Abbe condenser, NA 1.25, rack & pinion
Iris diaphragm
100-240V AC / 50-60Hz Auto-Switching Input
Light Options:
Variable halogen light source (12V/20W bulb) - Optional
Variable LED light source (3W bulb) - Standard

Revelation III

Model # Description

R3M-BN4A-DAL3 Binoc head, 4-10-40-100 DIN Achromatic objectives, LED illumination

R3M-TN4A-DAL3 Trinoc head, 4-10-40-100 DIN Achromatic objectives, LED illumination

R3M-BN4A-DPL3 Binoc head, 4-10-40-100 DIN Plan objectives, LED illumination

R3M-TN4A-DPL3 Trinoc head, 4-10-40-100 DIN Plan objectives, LED illumination

3M-MN4A-DAL3 Monoc head, 4-10-40-100 DIN Achromatic objectives, LED illumination

- •Medical-Grade Microscope
- Anti-Fungal Coating on All Lenses
- •Field-Ready Portable Options

WEIGHT AND DIMS

Height: 15" (380mm) Length: 9" (230mm) Width: 7" (178mm) Weight: 14 lbs. (6.4kg)



P 770.270.1394 F 770.270.2389

865 Marathon Parkway Lawrenceville GA 30046

LWScientific.com

SLS-7.5.3-L-012 | Rev 6