

MXU Centrifuge

With the choice of either an 8-place fixed-angle rotor or a 6-place swing-out rotor, the MXU centrifuge from LW Scientific is the perfect solution for both laboratories and doctors' offices. This USA-made centrifuge features a ZERO-RPM locking lid for safety and digital programming of time and speed ranging from 800-3400rpm. Achieve up to 1450g with the fixed-angle rotor and 1809g with the swing-out rotor resulting in clean, easy-to-read lines of separation. The LW Scientific MXU centrifuge has a 110-220v auto-switching power adapter.

ADVANTAGES

Accurate speeds for specific g-force
Ideal for blood, urine, and other fluids
Removable bowl for easy cleaning

DESIGN FEATURES

Rubber feet to prevent slipping
Easy removal of rotors for accessibility
Zero-rpm locking lid
Service Indicator Light
Stores last four programmed speeds

SPECIFICATIONS

800-3400 rpm
12-volt maintenance-free brushless motor
Programmable Digital Time and Speed
Fixed-angle rotor: 8-place 3-15ml tubes
Swing-out rotor: 6-place 3-15ml tubes
110/220v auto-switching power adapter
53 dBA +/- 3 dBA @ 18"-24"

INCLUDES

Power adapter (CE, UL, cUL approved)
8-place fixed-angle rotor
-8 small tube shields
(accommodates min. tube size of 10.25x65mm)
-8 large tube shields
(accommodates max. tube size of 17x130mm)
or
6-place metal rotor
-6 small tube shields
(accommodates min. tube size of 10.25x65mm)
-6 large tube shields
(accommodates max. tube size of 17x130mm)

CONTROLS

Touch Pad time and speed controls
800-3400 rpm programmable



MXU Centrifuge

Model #	Description
MUC-08AD-15T3	8-place, 3-15ml angled rotor, digital
MUC-06SD-15T3	6-place, 3-15ml swing-out rotor, digital

- Removable bowl for easy cleaning
- Programmable Digital Time and Speed
- Made in the USA

WEIGHT AND DIMS

Height: 9.12" (231 mm)
Length: 16" (406 mm)
Width: 14" (355 mm)
Weight without rotor: 17.1 lb (7.76 kg)
Weight with rotor: 18.35 lb (8.32 kg)

6-place swing-out rotor



8-place angled rotor



Designed, Engineered &
Assembled Proudly in America

LW
Scientific

P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com