# MX5 Centrifuge

With many configurations, the high-performance USA MX5 is the most versatile, large-volume centrifuge manufactured by LW Scientific. It spins many sizes of test tubes and microtubes at speeds up to 5,000rpm (4,220g), with microprocessor speed control. The digital menu is simple-to-use, and the digital interface ensures accurate processing. The zero-rpm interlocking lid will stay locked until the rotor has come to a stop, to ensure safety. The all-metal rotors are built for strength and durability, and the wide-mouth metal bowl allows for easy cleanup. The whisper-quiet, heavy-duty, brushless DC motor can handle high volume use and high-speed full loads for years to come.



### **ADVANTAGES**

9 programmable settings
Digital rotor selection
Ability to display rpm or g-force
High-strength stainless steel body
Adjustable acceleration and braking
Straight-line separation for accurate processing
when using swing-out rotors
Reduces high potassium results
Tighter gel line reduces re-mixing and re-spins
Very high g-forces for fast separation
Fast and accurate veterinary fecal flotations

# DESIGN FEATURES

LCD display
Rubber feet to prevent slipping
Impact-resistant lid
Zero-rpm interlocking lid
Large opening for easy cleanup
Metal rotor and tube shields
Cooler chamber for cooler samples

# **SPECIFICATIONS**

100-5,000 rpm (4,220g)
Brushless DC motor
9 Programmable settings
Dual-fuse input module (5a/250v)
110 volt (220 volt available)

## CONTROLS

Digital 1-99 minute programmable timer 100-5,000 rpm programmable speed

<b>#</b>	Description
15T1	MX5 Swing-Out, 8-Place Rotor, 1-15ml, 110 volt
15T1	MX5, 24-Place Angled Rotor, 3-15ml, 110 volt

MX5 Swing-Out Bucket Rotor, 24-Place, 3-10ml, 110 volt

MX5 Swing-Out Bucket Rotor, 4-Place, 50ml, 110 volt

• For 220 volt units, change the last digit of the part number from 1 to 2.

MX5 Centrifuge

### WEIGHT AND DIMS

#### Unit (No rotor):

Model #

M5C-08SU-

M5C-24AU-

M5C-BKSU-15T1

M5C-BKSU-50T1

Height: 11.5" (292 mm) Length: 16.5" (419 mm) Width: 16" (419 mm) Weight: 38.25 lb (17.34 kg) Packaged (No rotor): Height: 16" (406 mm) Length: 20" (508 mm) Width: 20" (508 mm) Weight: 43.2 lb (19.59 kg)



P 770.270.1394 F 770.270.2389

865 Marathon Parkway Lawrenceville GA 30046

SLS-7.5.3.L-117 | Rev 3 LWScientific.com

# **MX5** Configurations

#### MX5 centrifuge with 8-Place 15ml Swing-Out Rotor

#### M5C-08SU-15T1

- Spins 8, 1ml-15ml tubes
- Maximum 5,000rpm / 4,220g



#### MX5 centrifuge with 24-Place 15ml Angled Rotor

#### M5C-24AU-15T1

- Spins 24, 3ml-15ml tubes
- Maximum 5,000rpm / 3,829g



#### MX5 centrifuge with Swing-Out Bucket Rotor & 24-place test tube inserts

#### M5C-BKSU-15T1

- Includes 4, 6-place 3ml-10ml inserts
- Spins 24, 3ml-10ml tubes
- Spins 8, 15ml conical tubes
- Maximum 5,000rpm / 4,220g



#### MX5 centrifuge with Swing-Out Bucket Rotor & 4-place 50ml test tube inserts

#### M5C-BKSU-50T1

- Spins 4, 2-place 50ml tubes
- Maximum 5,000rpm / 4,220g



# **Optional Configurations**

#### **Crit Carrier Inserts**

#### CNP-04CC-PR77

- Spins up to 8, 75mm standard microhematocrit tubes
- Includes EZ Reader Card



#### **IRAP Angled Rotor for Equine applications**

#### C5P-RT04-IRFT & (4) C5P-TBS7-60IR

• Spins 4, 60ml IRAP tubes

