

# GAMMEX<sup>®</sup> Non-Latex

S U R G I C A L   G L O V E S

## Neoprene, chemical accelerator-free surgical glove for advanced allergy protection and peace of mind

### KEY FEATURES AND BENEFITS:

- Non-latex and 100% chemical accelerator-free to eliminate the risk of latex Type I and minimize chemical Type IV allergies and skin sensitivities
- Extended cuff length for extra protection
- Anti-slip, smooth external finish for secure instrument handling and grip
- Tested for use with chemotherapy drugs and a wide range of chemicals used in hospitals\*

*\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations.*

### Recommended For

- Damp donning
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Tested for use with chemotherapy drugs

### FEATURED TECHNOLOGIES

 **A.R.T.**  
Pinhole Barrier Technology

 **DERMASHIELD™**  
Donning Aid Technology



## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Non-Latex

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	Neoprene (Polychloroprene)
Powder Content	Powder-Free
Color	Green
Thickness	Standard
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I), Chemical (Type IV)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F1670, ASTM F739, EN 16523-1, EN 374:2003, EN 420:2003 + A1:2009, EN 421:2010, EN 455 1-4, EN ISO 374-5:2016, EN ISO 374-1:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	None

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	8511	8512	8513	8514	8515	8516	8517	8518

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com