



Endure Lab Coat

INSTRUCTION OF USE

Selection: Choose the appropriate size for a comfortable fit. The coats are available in small, medium, large, and extra-large sizes.

Features and Application: The lab coats are knee-length with long sleeves, providing ample coverage for protection against light fluid contact and particles. They have a single pocket on each side at hip level and a left-side breast pocket, useful for carrying small items.

Collar Style: Select between knit and traditional collar styles based on preference and specific requirements of your work environment.

Donning the Coat: Wear the coat by slipping your arms into the sleeves and securing the front with the centered snap closures.

Usage: The coats are designed for single use, ensuring hygiene and reducing the risk of cross-contamination in sensitive environments.

Disposal: After use, dispose of the lab coat responsibly in accordance with your facility's guidelines for medical waste.

Material and Comfort: Made from Spunbond-Meltblown-Spunbond (SMS) material, these lab coats provide good air permeability, balancing protection with wearer comfort.

Storage: Store unused lab coats in a clean, dry place to maintain their integrity.

Endure Industries, Inc.
45 N Fulton St
Homer, NY 13077
www.endureind.com
sales@endureind.com