



Endure Lab Jacket

INSTRUCTION OF USE

Selecting the Right Size:

Choose from small, medium, large, and extra-large sizes for a proper fit.

Features:

The jackets have a knit collar and cuffs, are hip-length, and made from Spunbond-Meltblown-Spunbond (SMS) material. They're designed for single use, offering protection against light fluid contact and particles while maintaining breathability.

Wearing the Jacket:

Put on the jacket and secure it using the five front-centered snap closures. The long sleeves provide additional protection.

Utility:

Use the side pockets at hip level and the left-side breast pocket for carrying small items.

Disposal:

After use, dispose of the jacket according to your facility's guidelines for medical waste.

Endure Industries, Inc.
45 N Fulton St
Homer, NY 13077
www.endureind.com
sales@endureind.com