



Etac Swift toilet seat raiser

Designed to assist users and carers with every day routines.

 **etac**[®]

Creating Possibilities

Etac Swift free standing toilet seat raiser

Swift free standing is a high quality toilet seat raiser which also can be used as a shower stool.

High stability

Swift free standing toilet seat raiser provides a very stable support. Its large base area and soft, large ferrules provide excellent floor grip.

Flexibility

The seat is height adjustable and can be set with a slight forward tilt. This facilitates standing up and relieves the hip joints. Its light weight makes the Swift Free Standing easy to move, while the legs are removable for easy transport.

3 in 1 product

Besides the toilet function, Swift free standing can be used as a shower stool and with complementing accessories it may also be used as a commode.



Swift free standing has smooth surfaces which enables easy cleaning.

Accessories



Description Back support, grey

Item no. 84005069



Description Soap dish, grey

Item no. 80209266



Description Soft Seat

Item no. 84005064 ✓



Description Soft Back

Item no. 84005065



Description Pan with Lid

Item no. 84005066 ✓

Description	Weight	Item no.
Swift Free Standing Toilet Seat Raiser	3.9 kg 8.6 lbs	81702020 ✓

Dimensions

Total width: 56 cm 22"
 Seat width: 54 cm 21¼"
 Seat height min-max. 42-57 cm 16½"-22½"
 Width between the arm supports: 45 cm 17¾"

Material

Polypropylene and aluminium




Clean the product with a solvent-free cleaning agent (pH level 6-9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



 **etac**[®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item