

Bed Rail Protector



Suits majority of bed rails

Reduces risk of injury in relation to bed rails

Soft to touch

These super soft, foam Bed Rail Protectors are designed to protect patients from injury through limbs striking or getting caught in the metal bed rails.

They are made from soft medical grade foam cushions covered with a soft two way stretch PU cover including an easy secure Velcro closure for secure attachment.

The cover can be wiped clean with mild detergent or Hypochlorite solution (up to 95°).

Designed for use on both static and collapsible bed rails, the Bed Rail Protectors allow easy access to the bed rail mechanism.

Sold as a pair, the Bed Rail Protectors are available in two sizes – standard and large – to fit the majority of bed rails.



SIZE	STANDARD	LARGE
CODE	BL5247	BL5254
HEIGHT (FOLDED)	335 mm	490 mm
WIDTH	670 mm	980 mm
LENGTH	1260 mm	1580 mm
THICKNESS	25 mm	
CONFORMS TO	ARTG155475	
WARRANTY	2 year	

