


## Pressure Relief Cushion



This two layer foam cushion provides improved pressure redistribution when seated.

This multi-purpose product is suitable for users with mild postural needs and comprises a high density foam support base with a super soft memory foam top layer.

Enclosed by a zippered two-way stretch waterproof polyurethane fabric cover for enhanced anti-shear and anti-friction benefits, with 360° zip and waterfall skirt to protect from moisture ingress.

Seams are stitched for strength, then welded for or maximum infection control.

o Formerly called the Comfort Cushion

<b>CODE</b>	BL7258
<b>HEIGHT</b>	75 mm
<b>WIDTH</b>	410 mm
<b>DEPTH</b>	410 mm
<b>CAPACITY</b>	100 kg
<b>WARRANTY</b>	1 year

