

Pressure Surface Range



BetterLiving offers a broad range of static pressure relieving surfaces to suit a range of pressure care needs.

Made from premium multi-density and/or memory foam material combined with high quality vapour permeable, water resistant polyurethane cover, all mattress products are vacuum packed for economical shipping and compact storage.

All models feature antibacterial, double strength welded covers. Seams are first stitched for strength, then welded or tape sealed for maximum infection control.

Mattresses conform to ARTG 176837 and BS7177:2008 (Medium Hazard) fire retardency.

Pressure Surface Range

TripleLayer^{King}



With all the benefits of the Triple Layer, this surface has increased depth, width and weight capacity.

TripleLayer^{Concave}



With all the benefits of the TripleLayer, combined with the added safety of raised, tapered sidewalls.

TripleLayer^{Plus}



With all the benefits of the TripleLayer, convoluted memory foam increases the weight capacity of this pressure surface.

TripleLayer



The three layered, zoned design of this high risk pressure surface is ideal for patients with poor mobility or uncertain skin sensitivity.

DIMENSIONS (W X H)	1050 x 1960 mm	865 x 1960 mm	865 x 1960 mm	865 x 1960 mm
DIMENSIONS DEPTH	180 mm	150 mm	150 mm	150 mm
LAYERS	3	3	3	3
SIDE SUPPORT	✓	✓	✓	✓
TURNABILITY	1 way	2 way	1 way	1 way
CAPACITY (MAX WEIGHT)	250 kg	220 kg	160 kg	160 kg
INFECTION CONTROL SEAMS	✓	✓	✓	✓
HINGED	✓	✗	✓	✓
CONVOLUTED (EGGSHELL)	✓	✗	✓	✗
RISK CATEGORY	very high	very high	very high	high
WARRANTY	6 year body (foam) 2 year cover	6 year body (foam) 2 year cover	6 year body (foam) 2 year cover	6 year body (foam) 2 year cover
CODE	BSPS0469	BSPS0473	BSPS0454	BSPS0453
REPLACEMENT TOP COVER	BSPS0470	BSPS3210	BSPS0466	BSPS0465

DoubleLayer



This dual layer surface for medium risk users is an excellent multi-use solution for hospitals, residential or home care environments.

SingleLayer



An economical, low risk pressure surface for general use in residential or home care.

Overlay



Made of memory foam, an affordable way to add pressure relief to a firm mattress or add extra comfort for frail and/or lightweight patients.

Underlay



Designed to provide “bottoming out” protection for alternating pressure overlays, while still minimising the height of the overall support system on the bed.

Extension



200 mm foam extension for use with the TripleLayer / DoubleLayer high pressure surface. Available in Standard (**S**) and King Single (**K**) widths

865 x 1960 mm	865 x 1960 mm	865 x 1960 mm	865 x 1960 mm	S 865 x 200 mm K 1050 x 200 mm
150 mm	125 mm	50 mm	50 mm	S 150 mm K 180 mm
2	1	1	1	2
X	X	X	X	X
2 way	2 way	4 way	4 way	NA
130 kg	130 kg	130 kg	130 kg	NA
✓	✓	✓	✓	✓
X	X	X	X	X
X	X	X	X	X
medium	low	Underlay or Overlay foam surfaces are designed to complement existing static or alternating mattresses to provide enhanced comfort or safety.		
6 year body (foam) 2 year cover	6 year body (foam) 2 year cover	6 year body (foam) 2 year cover	6 year body (foam) 2 year cover	6 year body (foam) 2 year cover
BSPS0452	BSPS0450	BSPS0455	BSPS0456	S BSPS0472 K BSPS0475
BSPS0464	BSPS0462	BSPS0467	BSPS0468	NA

