

# BetterLiving

## PressureSurface Overlay



**This memory foam overlay provides an affordable way to add pressure relief to a firm mattress or add extra comfort for frail and/or lightweight patients.**

This memory foam overlay is constructed from long wearing visco-elastic foam designed to conform to the body's shape. The overlay can be used on top of the existing mattress to provide additional user comfort and pressure redistribution. Foam is covered by an inner cotton cover for added protection.

The polyurethane top cover features multi-stretch material for reduced shear and friction. With waterfall skirt to provide protection from moisture ingress, seams are stitched for strength, then welded or tape sealed for maximum infection control. A PVC slip-resistant base helps maintain product position on bed. The cover is also pre treated for protection against dust mites, mould, mildew and bacteria and includes an asset management panel for recording facility details.

**Not to be used as standalone product**

- ✓ memory foam overlay for increased patient comfort
- ✓ immediate pressure relief
- ✓ fits over any standard single sized mattress
- ✓ max user weight 130 kg

**Formerly known as  
Memory Overlay**

*Specifications detailed over*

# BetterLiving

## PressureSurface Overlay



Vapour permeable, bacteriostatically treated, fire retardant and water resistant cover



Seams are stitched, then welded or tape sealed for maximum infection control



Quality fire retardant foam protected by an inner cotton cover

ITEM	<b>Overlay</b>
CODE	<b>BSPS0455</b>
CAPACITY	130 kg
HEIGHT	50 mm
WIDTH	865 mm
LENGTH	1960 mm
TOP COVER	Latex free polyurethane
	6 year body, 2 year cover
WARRANTY	Product must be removed from vacuum packaging prior to the 'Reconstitute By' date printed on the carton. Failure to do so will void the product warranty.
TURN	4 way
COMPLIANCE	<b>ARTG 176837</b>
FIRE RETARDANCY	Conforms to BS7177:2008 (medium hazard)
CARE	To clean, wipe down the top cover and PVC base with warm water and a solution of sodium hypochlorite or similar (up to 10,000 ppm available chlorine). Top cover can be machine washed up to 80°C.
ACCESSORIES	Replacement top cover <b>BSPS0467</b>
CARTON SIZE	950 x 260 x 260 mm