

# BetterLiving

## PressureSurface TripleLayer<sup>KING</sup>



**This foam pressure surface has increased depth and width for added comfort and capacity for the larger user. The surface is sized to fit a standard size king single bed, and supports user up to a maximum capacity of 250 kilograms.**

The three layered, zoned design is ideal for patients with poor mobility or uncertain skin sensitivity. Made from premium material, each foam layer varies in density to maximise the redistribution of pressure.

A top layer of high density, temperature sensitive visco-elastic memory foam conforms to the body's shape, working in combination with the higher resistance support layer to disperse pressure evenly. The core layer and edges provide form, strength and stability. Foam is covered by an inner cotton cover for added protection.

The polyurethane top cover features multi-stretch material for reduced and friction. With waterfall skirt to provide protection from moisture ingress, seams are stitched for strength, then welded or tape sealed for maximum infection control. A PVC slip-resistant base helps maintain product position on bed. The cover is also pre treated for protection against dust mites, mould, mildew and bacteria and includes an asset management panel for recording facility details.

- ✓ sized to fit a standard size king single bed/bed base
  - ✓ three layer, zoned design with hinge for profiling beds
  - ✓ multi-stretch top cover with welded / taped seams for infection control
  - ✓ max user weight 250 kg
- Formerly known as King Single**

*Specifications detailed over*

# BetterLiving

## PressureSurface TripleLayer<sup>KING</sup>



Vapour permeable, bacteriostatically treated, fire retardant and water resistant cover



Seams are stitched, welded and taped for maximum infection control



Multi-density foam layers for maximum strength and comfort; protected by an inner cotton cover



Six hinge system to conform to profiling beds

|                 |   |
|-----------------|---|
| ITEM            | <b>TripleLayer King</b>   |
| CODE            | <b>BSPS0469</b>   |
| CAPACITY        | 250 kg  |
| HEIGHT          | 180 mm  |
| WIDTH           | 1050 mm   |
| LENGTH          | 1960 mm   |
| TOP COVER       | Latex free polyurethane   |
| WARRANTY        | 6 year body, 2 year cover<br>Product must be removed from vacuum packaging prior to the 'Reconstitute By' date printed on the carton. Failure to do so will void the product warranty.                                  |
| TURN            | 1 way   |
| COMPLIANCE      | <b>ARTG 176837</b>  |
| FIRE RETARDANCY | Conforms to BS7177:2008 (medium hazard)   |
| CARE            | To clean, wipe down the top cover and PVC base with warm water and a solution of sodium hypochlorite or similar (up to 10,000 ppm available chlorine).<br>Dry thoroughly before use. Do not machine wash the top cover. |
| ACCESSORIES     | Replacement top cover <b>BSPS0470</b>   |
| CARTON SIZE     | 1150 x 380 x 380 mm   |