



# Etac HandyGrip

Tools to help and support in mobilization and transfer.

 **etac**<sup>®</sup>

Creating Possibilities

# Adjustable grip to aid turning in bed









Etac HandyGrip provides help with independent turning in bed. It can also be used when self transferring from wheelchair to bed. Soft, compact and easy to hide away. The adjustable buckle allows the length and position of HandyGrip to be varied, providing the user with a good grip at the required position.

*Soft, compact and easy to hide away.*



## Intended Use


- Support independent turning in bed
- Support independent transfer to bed

Description	Width	Length	Order Code
HandyGrip	50 mm	1600 mm	IM461 
<b>Specifications</b>			
Polyester			 max 150 kg  max 330 lbs
Belt Buckle: Acetal and nylon			
<b>Washing Instructions</b>			
Do not use fabric softener or bleach.			
			 



 **etac**<sup>®</sup>  
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item