

Etac Optima toilet arm supports

Designed to assist users and carers with every day routines.



Etac OptimaL toilet arm supports

Modern and elegantly designed arm supports.

Function and simplicity in an elegant combination

OptimaL is an elegant wall mounted support rail. The leg secures the stability and the small wall attachment enables mounting close to the toilet. OptimaL stays safely folded until it is flipped down with the grey handle.

Accessories

See our wide range of accessories at the end of this brochure



Etac Optima toilet arm supports

Optima toilet arm supports provides a firm and comfortable grip and folds smoothly to the wall.

Comfortable and well-fitting

Optima arm support cover has a softly rounded end which provides a comfortable grip. The supporting leg has an adjustable foot which ensures a perfect floor contact. The large wall plate gives a firm fixation to the wall.



Etac Optima Lockable toilet arm supports

Optima lockable offers extra security and support.



Optima Lockable is available in two different lengths.

Description Optima Lockable	Weight 2.4 kg 5.3lbs	Item no. 80303004
Dimensions Length arm support Wall bracket	71.5 cm 281/4" 16.5 x 30.4 cm 61/2" x	12"
Material Wall bracket Powder-painted aluminium Arm supports tubing Powder-painted steel tubing Locking device Polyamide Arm supports washer Polypropylene		
	solvent-free cleaning agent	
or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.		



Etac Optima accessories



Description Item no. Attachment for signal device 80303027

Suitable for Optima Optima Lockable



DescriptionItem no.Floor pillar80303028

For floor-mounting, positions arm supports height at 80 cm $31 \frac{1}{2} ^{\prime\prime}$

Suitable for Optima Optima Lockable



DescriptionItem no.Supporting leg80303040

For relieving strain on walls.

Suitable for Optima Lockable only



 Description
 Item no.

 Extension, 10 cm 4"
 80303030

 Extension, 20 cm 734"
 80303031

 Extension, 30 cm 1134"
 80303032

Suitable for Optima Optima Lockable



DescriptionItem no.Toilet paper holder83030001Complete. fitted to the arm supports.

Suitable forOptimaL
Optima
Optima Lockable







"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

