



# Etac SlingOn

Systems and tools for horizontal positioning, turning, pressure care, moving higher up in bed.

 **etac**<sup>®</sup>

Creating Possibilities

# SlingOn

## Multi-tool for application of lifting sling and positioning

Etac SlingOn is primarily designed for easy insertion of a hoist sling in bed without the need for rolling or moving the user. SlingOn consists of two separate MultiGlide ultra low friction nylon slide sheets held together by snap buttons in the corners and handles along the sides. This holds the slide sheets well in place and facilitates application. Etac SlingOn is an optimal tool for users that needs frequent hoisting, heavy users, or users that may not be rolled for application of a sling. And since the material is the same as used for the MultiGlide slide sheets, it can also be used for all positioning in bed.

SlingOn is put in place under the user using a simple folding technique, whether for application of sling or for positioning. The hoist sling is then inserted between the two layers of MultiGlide. The ultra low friction facilitates sliding the sling in place underneath the user, and adjusted if necessary. When the lifting sling is in place, remove the top part of SlingOn.

**Note: The user is always lifted in the sling, never in SlingOn!  
Always remove the top part of SlingOn before hoisting!**



### Intended Use

- Application of a hoist sling without rolling
- Minimise friction under pressure points
- Turning
- Positioning
- Move higher up in bed

*Etac SlingOn is the optimal solution for insertion of a hoist sling without rolling.*



### Description

| Description   | Width        | Length       | Item no.  |
|---------------|--------------|--------------|-----------|
| SlingOn       | 80 cm 31.5"  | 170 cm 66.9" | IM437/S ✓ |
| SlingOn Large | 120 cm 47.2" | 150 cm 59.1" | IM437     |

### Specifications

Nylon  
Polyester straps  
Push buttons

### Washing Instructions


Can be wiped with disinfectant.  
Do not use fabric softener or bleach.





 **etac**<sup>®</sup>  
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item