



Etac Sling

Tools to help and support in mobilization and transfer.

 **etac**[®]

Creating Possibilities

Sling with handles

Etac Sling with handles acts as extension of the arms, allowing the caregiver to offer support and assistance whilst maintaining a good ergonomic position and grip, without encroaching on the service user's personal space.

Etac Sling is a flexible tool that allows the caregiver to maintain a more ergonomic grip. It can be used to support many different transfers such as, moving higher up in bed, transfer between sitting positions, etc.

Sling has 4 handles on each side. The inside has a non-slip material to help prevent sliding out of position during the transfer, whilst the outside has a low friction surface to aid placement and removal



Etac Sling is a flexible tool offering a number of different functions.

Intended Use

- Support moving higher up in bed
- Support transfer between bed, wheelchair, chair, etc
- Support stand up / sit down
- Moving back / forward in chair / wheelchair



Description

Sling, 2 x 4 Handles

Width

160 mm

Length

1900 mm

Order Code

IM428

Specifications

Backside: Nylon
 Padding: Polyester
 Inside: Nylon with polyurethan
 Handles: Polyester

Washing Instructions


Do not use fabric softener or bleach.





 **etac**[®]
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item