

Etac Swift commode

Designed to assist users and carers with every day routines.



Etac Swift Commode

Swift Commode is a three in one product with many smart functions.

Comfortable and functional

Swift Commode is discreet and comfortable. It can be set with a slight forward tilt to facilitate sitting down and standing up. The integrated pan has a tightly securable lid that locks automatically when lifted. Swift commode is supplied complete with soft seat pads.

3 in one

Swift commode can also be used as a shower chair or as a freestanding toilet seat raiser.

Easy cleaning

With smooth surfaces and removable padding, the Swift commode has been designed to enable easy cleaning. Its light weight and the integrated handle make it easy to transport and clean.

Easy handling

Swift Commode comes flat packed and is easy to assemble and height adjust, without using tools. The legs can be removed for transport.



Integrated pan lid

The pan with lid is integrated in the seat surface.



Order Code Description Capacity Swift Commode 130 kg 81702030 Pan with Lid 84005066 Soap dish, grey 80209266

Weight

5.4 kg 12 lbs

Dimensions

Total width: 56 cm 22" Seat width: 54 cm 211/4" Seat height min-max. 42-57 cm 161/2-221/2" Width between the arm supports: 45 cm 173/4"

Material

Polypropylene and aluminium







Clean the product with a solvent-free cleaning agent (pH level 6-9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.

With removable padding, the Swift commode has a discrete design and is also comfortable for dressing and undressing.

Accessories

Soap dish



Pan with Lid







"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances."

