

# 925 980 9108 WWW.SAPLACOR.COM



# **AEGIS** Neonate Medical Wrap

AEGIS (The Seatbelt of Maternity) crash tested, award winning neonate transport system, is the solution for safe couplet-care transports to minimize risk. AEGIS provides upright, hands-free, skin to skin proximity for optimal assessment & monitoring throughout transport. Transporting together enhances optimal physiological outcomes, alleviates the need of an additional food source, heat source and ambulance.

A baby delivery is one of the most memorable calls you will get, Choose AEGIS.



#### **Features & Benefits**

- Decrease mother & newborn stress.
- Newborns face is visible to assess & monitor.
- Nano-guard supports newborn's head & neck.
- Properly positioned Skin to Skin per AAP guidelines.
- 4-way stretch material harnesses entire newborn baby & absorbs energy from impact.
- Patient use on or off the gurney, in or out of the ambulance.
- Quick & easy assembly, wearability & storage.
- Allows tubes/lines to pass through top or bottom easily
- Open bottom to allow umbilical cord & placenta to remain intact during transport
- X-Static knitted throughout for antimicrobial benefits.
- X-Static & Skin to Skin reduce the risk of hypothermia.
- Skin to Skin / Oxytocin reduce the risk of hemorrhage.
- Option to wash & re-use or single patient use.

### Size / Weight Chart / SKU#

Size 1	29-32 inches	4023267618
Size 2	32-37 inches	4023267619
Size 3	37-42 inches	4023267620
Size 4	42-47 inches	4023267621
Size 5	47-52 inches	4023267622

#### Recommended for newborns 4-14 lbs ( $1.8\text{-}6.3\ kg$ ).

## **Product Testing**

FDA cleared General Wellness

Tested @ Calspan to adult SAE crash testing - Passed

Weight & chemical testing @ Intertek- Passed

X-Static Silver Technology @ Noble labs - Certified



# ABE() Seatbelt attachment

## ONLY pregnancy seatbelt for ambulance use.

- One size fits all
- Made of Biothane, Reusable, wipe clean, easy to store.
- The H-harness design reduces the risk of fetal distress & risk of abruption.
- Passed frontal, rear & side crash testing at Calspan.
- ABEO must be used in combination with AEGIS to transport mom/baby in the ambulance.



• Pediatric

SKU # 4023202697