




Project Type IB - EN	Date: November 20 2020	PO: 24802-00	 Contact: Graphic's Department graphics@thermor-ins.com
	Customer: Thermor Ltd.	Revision: A	
UPC (For Reference Only)	Job Number: TH6747		
	Item Number: 60058		
	Designer: Natalia Siguenza		
	Reason For Project: New Item: New Branding: Design Update: X N/A:	Artwork is not to be amended or scaled. If any changes are required please send through Thermor office with instructions, and we will amend here and resend artwork through. <b>PLEASE RETURN WITH SIGNATURE</b>	
		<b>Colour Breakdown</b> 	

# BIOS living

## 60058 Bamboo Shower Bench Instructions



- When screwing the bolts in, do not tighten all the way until all 8 bolts are approximately 80% tight .
- Then proceed to fully tighten all 8 bolts.

### MAINTENANCE AND CLEANING:

- After use remove the bench from the shower as to not get moldy.
- Use a towel to wipe down the bench if it gets wet.
- To clean use a mild soap and luke warm water. Make sure to dry completely.

### ONE YEAR WARRANTY:

The BIOS Living Bamboo Shower Bench Backrest has a one year warranty to be free of manufacturing defects in materials and workmanship under normal applications. If this product becomes inoperable due to defect and requires repair, return the product with all component pieces and proof of purchase to the address listed below. This warranty does not cover any shipping/transport costs. This warranty does not apply if the product is subject to misuse, neglect, rough handling or damage.

Ship the unit prepaid and insured (at owners option) to:

**Thermor Ltd.**  
 Attn: Repair Department  
 16975 Leslie Street  
 Newmarket, ON L3Y 9A1

www.biosmedical.com  
 Email: support@biosmedical.com

BIOS | Medical  
 MANUFACTURED BY / FABRIQUÉ PAR :  
**THERMOR LTD.**  
 16975 LESLIE STREET  
 NEWMARKET, ON L3Y 9A1  
 MADE IN CHINA / FABRIQUÉ EN CHINE  
 WWW.BIOSMEDICAL.COM