


## IMPLANT ONE BALL ABUTMENT

Implant One Ball Abutments are used to secure dentures using Rhein 83 Standard 2.5mm Attachments.

- Made from 6AL4V E.L.I. Titanium
- Available for all Implant One platforms in 4 different cuff heights.
- Rhein 83 components can be ordered directly from Rhein 83 USA ([www.rhein83usa.com](http://www.rhein83usa.com))

### Sizes Available and Part Numbers

CUFF HEIGHT	300 series	400 series	500 series
2 mm	ABL-0320-00	ABL-0420-00	ABL-0520-00
3 mm	ABL-0330-00	ABL-0430-00	ABL-0530-00
4 mm	ABL-0340-00	ABL-0440-00	ABL-0540-00
5 mm	ABL-0350-00	ABL-0450-00	ABL-0550-00
PROFILE DIAMETER	4 mm	5 mm	5.5 mm



Use the DEH-0000-00 Ball Abutment Driver to install into implant (2.3mm hex)

**Maximum Torque**  
**20Ncm (300 Series)**  
**30Ncm (400 & 500 Series)**

Ball compatible with Rhein 83 standard 2.5mm components

Specifications
Material: 6AL4V E.L.I. Titanium
Profile Diameter: 4 mm – 5.5 mm
Cuff height: 2 mm – 5 mm

