




IMPLANT ONE **TEMPORARY ABUTMENT**

Designed to be temporarily seated with a temporary attachment or crown while the final restoration is being made.

- Made from 6AL4V E.L.I. Titanium.
- Available in 6 different cuff heights ranging from .5 mm to 5 mm.
- Titanium retaining screw included.
- Temporary crown must be configured to prevent loading.



| Specifications |
|---------------------------------|
| Material: 6AL4V E.L.I. Titanium |
| Profile Diameter: 4 mm – 5.5 mm |
| Cuff height: .5 mm – 5 mm |



Sizes Available and Part Numbers

| CUFF HEIGHT | 300 series | 400 series | 500 series |
|------------------|-------------|-------------|-------------|
| 0.5 mm | ATM-0305-00 | ATM-0405-00 | ATM-0505-00 |
| 1 mm | ATM-0310-00 | ATM-0410-00 | ATM-0510-00 |
| 2 mm | ATM-0320-00 | ATM-0420-00 | ATM-0520-00 |
| 3 mm | ATM-0330-00 | ATM-0430-00 | ATM-0530-00 |
| 4 mm | ATM-0340-00 | ATM-0440-00 | ATM-0540-00 |
| 5 mm | ATM-0350-00 | ATM-0450-00 | ATM-0550-00 |
| PROFILE DIAMETER | 4 mm | 5 mm | 5.5 mm |