


**Important!** The [Implant One 300 Series Ti-base abutments](#) have been redesigned to improve the single-tooth restoration/abutment connection in narrow spaces.

Please [download](#) and use the latest libraries for 3Shape and exocad, available on our website (details below), in order to restore with these new, smaller diameter abutments. Failure to do so will result in restorations being made that will not fit on the abutments.

### New 300 Series Ti-base

Part IDs:  
ATB-0305-01  
ATB-0315-01  
ATB-0325-01


**Now available**



### Old 300 Series Ti-base

Part IDs:  
ATB-0305-00  
ATB-0315-00  
ATB-0325-00

**No longer available**



Download the latest Implant One 3Shape and exocad libraries here:  
<https://implantlogistics.com/pages/lab-downloads>

**Please Note:** It is necessary to remove the previous Implant One library *before* installing the new one.

Questions? Call us at 608-498-4855, or  
email [clinicalsupport@implantlogistics.com](mailto:clinicalsupport@implantlogistics.com)