



Implant-One® Temporary Abutment


Implant-One Temporary Abutments are used to temporarily restore an implant at the time of placement

Made of Grade 5 Titanium.

Available in 5 different cuff heights ranging from 1 mm to 5mm, for both 100 and 200 series.

Titanium retaining screw included.

Temporary Crown must be configured to prevent loading



100 Series
200 Series


IT121-10
IT221-10

IT121-20
IT221-20

IT121-30
IT221-30

IT121-40
IT221-40

IT121-50
IT221-50



Temporary Abutment
with Screw
and Implant

Implant-One Temporary Abutment Part Numbers