



## Implant-One Ti Base (hybrid abutment interface) 3-Shape and Sirona

**Implant-One Non-Prepable Ti Base is used for 3-Shape designed hybrid abutments and are compatible with Sirona in-office milling system**


Made of Grade 5 Titanium.

Available in .5mm cuff height for both 100 and 200 series.


Titanium retaining screw included.

Intended for use as a Titanium base for Titanium / Zirconium custom abutments either as a generic Ti Base or in conjunction with Sirona

Sirona Scan abutment also available as shown



### Implant-One Ti Base Part Numbers



Ti Base shown with Screw and Implant