

# READY TO WEAR

COMPRESSION & POST-SURGICAL



866.251.0076  
[www.wearease.com](http://www.wearease.com)



**Wear Ease** is the leader in providing “Ready To Wear” therapeutic compression garments. Our products are designed to provide mild compression to help reduce swelling and discomfort caused by medical conditions such as lymphedema, lipedema, edema, and post-surgical recovery. We have consulted with medical professionals, fitters, and patients alike to create our garments from high performance fabrics to provide superior support with our unique designs.

This year, we are proud to announce the addition of several new items to our collections, so that more people can experience the relief of Wear Ease compression garments! With Wear Ease, you can trust that you will be receiving the highest quality garment available on the market.



#### **Karena Textured Compression Bra -**



**Style 770** is a revolutionary new product for comfortable, all-day support. The textured design also provides compression with micro massage to help reduce swelling and soften radiation fibrosis and scarring.

#### **High Waist Compression Short in Two Lengths -**

**Styles 614S** (short) and **614L** (long; 3” longer than the 614S), and a wider size range (S, M, L, XL, 1X, 2X, 3X, 4X, 5X, 6X). Provides comfortable compression and great containment for the abdomen, hips, buttocks, and thighs. These shorts are the same sizing as our High Waist Compression Capri, Style 615, just different lengths.

#### **Taylor T Compression Shirt -**

**Style 914** Now available in black. Women’s long sleeve compression shirt designed with a double layer of our popular high performance nylon spandex fabric providing compression for the torso, abdomen, back and shoulders.



#### **Katy T -**

**Style 918** Now available in black. Pocketed for additional relief from axillary swelling. This shirt is now available in two colors, nude and black.

Listening to fitters, medical professionals and patients, we know what is important.

- Superior post-surgical and compression products that improve quality of life by being comfortable, functional, attractive and affordable.
- Continual product updates based upon medical professionals’ feedback and patient needs.
- Personalized attention and unforgettable customer service. Our knowledgeable and professional customer service team is always available to assist you.

# TABLE OF CONTENTS



Style 785

Style 611

<b>PRODUCT AND COMPANY INFORMATION</b> .....	2
<b>MEASUREMENT GUIDELINES</b> .....	3
<b>FITTING GUIDELINES</b> .....	3
<b>COMPRESSION SHAPEWEAR</b> .....	4
<b>STEP IN • STEP OUT Instructions</b> .....	5
Katy T - 918 .....	7
Taylor T - 914 .....	9
Ava Compression Camisole - 911 .....	9
Slimmer - 910 .....	10
Compression Camisole - 912 .....	10
Compression T - 915 .....	11
Compression Crop Top - 785 .....	11
Sydney Bra - 780 .....	13
Crisscross Shaper - 970 .....	13
High Waist Compression Short - 614S & 614L.....	15
High Waist Compression Capri - 615 .....	15
Compression Short - 612 .....	16
Compression Capri - 611 .....	16
Short Compression Capri - 613 .....	16
Compression Bodysuit - 1000.....	17
Ellen Compression Bodysuit - 1020 .....	17
<b>COMPRESSION BRAS &amp; TOPS</b> .....	19
Compression Bra - 790/791/792 .....	20
Sternum & Genital Pad - 309 .....	20
Sizing & Fitting for Style 790 .....	21
Karena Bra - 770 .....	23
Andrea Compression Shirt - 960/962.....	24
Women's Torso Compression Vests - 950 & 951 .....	25
Women's Compression Vest - 796/797 .....	25
<b>MEN'S COMPRESSION VESTS &amp; SHIRTS</b> .....	27
Andrew Shirt - 961/963 .....	27
Eric Compression Shirt - 916 .....	27
High Waist Compression Short - 614S & 614L.....	28
High Waist Compression Capri - 615 .....	28
Men's Torso Compression Vests - 953 & 954 .....	29
Men's Compression Vest - 798/799 .....	29
<b>ACCESSORIES - 400/402/403/410/510/514/516</b> .....	30
<b>USING SWELL SPOTS WITH WEAR EASE</b> .....	31
<b>POST-SURGICAL CAMISOLES</b> .....	32
Beth Camisole - 903/904 .....	32
Dawn Camisole - 900/901 .....	33
<b>POST-SURGICAL BRAS</b> .....	34
Allyson Bra - 750/751/752 .....	34
Grace Bra - 740/741/742 .....	35
Compression Bra - 790/791/792 .....	35
<b>PRESCRIPTION FORM</b> .....	36-37

# Product and Company Information

## Garment Care

Laundry garments according to the care tag attached to the garment. For best results, wash in a lingerie bag. Do not use chlorine bleach or fabric softener. Tumble dry low or hang to dry.

## Garment Replacement

Wear Ease Compression garments should be replaced every three to six months, depending on wear, to ensure consistent compression levels. Garments should be replaced when a patient gains or loses weight.

## Processing & Shipping

All parcels ship via USPS or FedEx at the discretion of Wear Ease, Inc. Upon request, Wear Ease will overnight, 2nd day or drop ship. For standard orders received before 3pm US Mountain Standard Time Monday through Friday, in-stock items will be processed and shipped that same day. Standard orders for in-stock items received after 3 pm will ship the next business day. We will notify you if an item is out of stock and will take longer than 3 days before shipping. Wear Ease is not responsible for negligent acts by delivery service, including lost, stolen or damaged parcels or weather related shipment delays.

## 4 WAYS TO ORDER

Online 24/7 - Call to set up.

By email: [info@wearease.com](mailto:info@wearease.com)

By phone or fax

Phone: 1-866-251-0076

Fax: 1-208-445-0957

(MST) Mon-Fri 8am-5pm

## Returns & Exchanges

Our goal is that you are completely satisfied with your order, but in case you need to exchange or return an item we will accept returns and exchanges within 30 days of receipt. We will issue you a credit, less shipping and restocking fees. Please contact customer service at 866-251-0076 or email [info@wearease.com](mailto:info@wearease.com) to obtain a Return Authorization (RA) number. Due to health reasons, special care needs to be taken when trying on garments. **Please note:** All returns are subject to inspection before refunds or exchanges are issued. Refunds or exchanges will not be issued for garments that are soiled, altered, washed, worn, damaged, emitting odors (including detergent, lotion, fragrances, sweat, smoke or food) or indicating that they have come into contact with other irritants/allergens such as pet hair or dander. Items returned for manufacturing defect are accepted regardless of condition. Please inspect your order within 10 days of receipt.

**Please refer to our website: [www.wearease.com](http://www.wearease.com) for additional details and current updates.**




## Measurement Guide for Wear Ease Garments

To ensure the best fit of our garments, follow these guidelines to determine correct size.

### Bust

Measure at the fullest point of the breast/chest (around nipple line).\*



### Waist

Measure at the navel, between the lowest rib and hips, identified by bending the body to one side.

### Hips

Measure the fullest part of the hips and buttocks.



## Fitting Guidelines for Wear Ease Compression Garments

With the right fit, your garments should feel snug and comfortable. **Please note that your regular bra, shirt or pant size will not be the same size as a Wear Ease garment.** \*Please measure yourself or have someone measure you and use our size charts.

### Tops, Shirts and Bodysuits

- Styles 910, 911, 912, 914, 915, 916, 918, 95X, 96X, 970, 1000 and 1020 are sized by fitting the breast/chest, waist and hip measurements. These garments should fit snug but not tight. Compare your measurements to the closest fit per the size chart. If your measurements fall between sizes, choose the smaller size. (If unsure of the fit consult your lymphedema therapist.) When fitted properly the relief is evident.

### Bras and Vest

- Styles 770, 780, 785, 790, 796 and 798 are sized by fitting the breast/chest measurement. \*Take the measurement around the breast/chest at the fullest point (nipple line) with a bra on and compare to the size chart for each item. If your measurement is between sizes, choose the smaller size. For additional information on fitting the Compression Bra, style 790, please refer to page 17.

### Capris and Shorts

- Styles 611, 612, 613, 614 and 615 are sized by fitting the waist and hips. Measure the waist and hips and compare to the closest fit per the size chart. If between sizes, choose the smaller size. (If unsure of the fit consult your lymphedema therapist.) To don, place your feet into the leg openings and pull the garment up, similar to pantyhose. Align crotch opening as you pull the shorts or capris all the way up, comfortably positioning the rise.

\*If you will be wearing a prosthesis, measure with it on to get your full breast measurement.

Please call customer service with questions: 866-251-0076

# Compression Shapewear

## BY WEAR EASE®

Compression garments designed to provide support for the body and can be used in the management of medical conditions such as lipedema, lymphedema and edema. Specifically designed to help reduce swelling caused by these conditions while providing comfort. Compression shapewear is recommended for long-term wear during daily activities including sleep, travel, active wear and even swimming. By wearing compression shapewear regularly throughout the day, it can help improve comfort and quality of life.

### Compression Levels

Wear Ease garments provide mild compression, 10-20 mmHG, depending on the garment and the fit. This provides support to the affected area without being too restrictive or uncomfortable. The compression helps promote proper circulation which can help reduce swelling in the affected area.

### Features & Benefits:

- "Step into" garment designed for ease of donning
- High cut underarms and back targets swelling
- **Built-in pockets** for prostheses if needed
- Snug enough to hold chip pads and swell spots in place
- Hidden seams and sturdy construction
- UPF 50+
- Comfortable garments for uninterrupted wear
- Lightweight, supportive & moisture resistant
- Firm fabric for control, compression and containment
- High performance 80% Nylon 20% Spandex fabric, **Latex free**



Style 780



Style 785



Style 910



Style 911



Style 912



Style 914



Style 915



Style 916



Style 918



Style 970



Style 1000



Style 1020



Style 611



Style 612



Style 613



NEW - Style 614S



NEW - Style 614L



Style 615

# STEP IN • STEP OUT

## WEAR EASE SHAPEWEAR

Donning and doffing our shapewear is simple and painless!

### STEP IN



1

Stretch neckline and gather sides of garment.



2

Step in to garment through neckline.



3

If limited balance, sit on a chair or stool - step into garment.



4

Slide garment up over hips.



5

Slip arms in to straps.



6

It's on!

### STEP OUT



1

Slip arms out of straps.



2

Slide garment down over hips. If limited balance, sit on a chair or stool.



3

Slide garment down and step out.

COMPRESSION +  
EVERY DAY STYLE



*Tops 918, Bottoms 614S & 615*



## Katy T (Axilla Compression T) 918

**Pocketed for additional relief from axillary swelling. Compression for underarm (axilla), arm, lateral chest wall, shoulder and back. Can accommodate the Wear Ease Axilla pad, style 305.**

This new T comes with pockets sewn into the lateral chest wall and the inside of the upper arm for placement of chip pads, swell spots or foam. Keep those items in place with a pocket. Great for travel, sleeping and active wear. UPF 50+.

\*See page 4 for more features.

**Colors:** Black ● Nude ○

Size	S	M	L	XL	1X	2X
Bust	36-38	38-40	40-43	43-46	46-49	49-53
Waist	31-34	34-37	37-40	40-43	43-46	46-49
Hip	38-40	40-42	42-45	45-48	48-51	51-54



Shirt is shown inside-out in photos below.  
Shown with SP-06 Small Circular Pad and Style 305 Wear Ease Axilla Pad.



*I featured the Katy T compression shirt and the 2 Axilla Pads by Wear Ease in the single arm MLD demonstration for the NLN Lymphedema Awareness Month celebration. I also wore mine to sleep for 2 nights with the pads under my breast and into the axilla. Although I do not have breast cancer Lymphedema, I do have fibrosis in my breasts and will often pool fluid. I was so pleased to see the difference in my breasts just by sleeping in the garment and using the 2 pads. I think there may be a whole new market for you to consider in treating fibrotic breasts!*  
- Louisa Boyd OTR, CLT-LANA

## Axilla Pad BY WEAR EASE® 305



**Treats swelling and fibrosis in the underarm and the lateral chest wall.**

This pad is held in place by the pockets in each side underarm area of the Katy T shirt. One pocket on the lateral chest wall and the other pocket on the underside of the sleeve.

**Color:** Tan ●



# COMPRESSION + STYLISH COVERAGE



*Tops 914, Bottoms 614S & 615*

# Taylor T 914

Compression for underarm (axilla), back, breast/chest and abdomen.

**Features:**

- Snug fit for consistent and full torso compression
- Double layer fabric for body, single layer for sleeves
- Seamless front, no seam under the bust
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Nude ○ Cappuccino ● Lt Coral ●



# Ava Compression Camisole 911

Compression for underarm (axilla), back, breast/chest and abdomen.

**Features:**

- Snug fit for consistent and full torso compression
- Wider strap than style 910
- Seamless front, no seam under the bust
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Nude ○ Lilac ●

\* 1- Model is wearing size Small.  
Bust: 34.5", Waist 28", Hips 39.5"



\* 2- Model is wearing size Medium.  
Bust: 40", Waist 31.5", Hips 42"


Size chart for styles 911 & 914

Size	XS	S	M	L	XL	1X	2X	3X
Bust	32-34	35-37	38-40	41-43	44-46	47-50	51-54	55-58
Waist	29-30	31-32	33-34	35-37	38-40	41-43	44-47	48-51
Hip	34-36	37-39	40-42	43-45	46-48	49-52	53-56	57-60

Size chart is based on actual measurements in inches and not bra size. If measurement falls in between sizes, always size down. When determining size, fit breast/chest measurement first.

*VERY impressed with the Ava Compression Camisole. I chose no reconstruction and no prosthetic after my mastectomy, so having a flatter camisole without extra fabric is perfect for me. I'll be buying more in the future. And I love the color Lilac!*

- Kathy



*You guys are terrific! And the camisole ~ it's not just a garment, it improves my quality of life. Thanks and long life to you and the company.*

BKH

## Slimmer 910

**Compression for underarm (axilla), back, breast/chest and abdomen.**

For trunk compression, the Slimmer provides immediate relief from swelling discomfort. Therapists recommend this compression garment for sleeping. Also great for travel and sportswear.

**Features:**

- Snug fit for consistent and full torso compression
- A-line style that won't ride up
- High cut front, underarms and back
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Navy ● Nude ○ White ○

**Fitter's Notes for styles 910 & 912**

- Although these garments are pocketed, they may not provide enough support for all breast forms, in which case a bra is required.



**Size chart for styles 910 & 912**

Size	XS	S	M	L	XL	1X	2X	3X
Bust	31-33	34-36	37-39	40-42	43-45	46-48	49-51	52-54
Waist	26-28	29-31	32-33	34-36	37-39	40-42	43-45	46-48
Hip	34-36	37-39	40-42	43-45	46-48	49-51	52-54	55-57

Size chart is based on actual measurements in inches and not bra size. If measurement falls in between sizes, always size down. When determining size, fit breast/chest measurement first.

FOR THE PETITE FIGURE



## Compression Camisole 912

**Compression for underarm (axilla), back, breast/chest and abdomen.**

Similar to our popular Slimmer, style 910, but with wider, shorter shoulder straps, a shorter length and higher fitting arm holes. The variety of beautiful colors make it a "wear alone" camisole for summer and layering. (If taller than 5'4" consider ordering style 910 instead.)

**Features:**

- Snug fit for consistent and full torso compression
- Good option for the petite figure under 5'4".
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Nude ○ Rose Pink ● Teal ●



# Compression T 915

**Compression for underarm (axilla), back, breast/chest and abdomen.**

This T is specifically designed to provide compression for the entire torso. Great for travel and sportswear.

**Features:**

- Double layer of fabric for comfortable compression
- Snug enough to hold chip pads in place
- Sleeves for additional coverage (not compression)
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Ivory ○ White ○

**Fitter's notes**

- Although this garment is pocketed, it may not provide enough support for all breast forms, in which case a bra is required.



Photo shows how Style 915 overlaps the edge of your compression sleeve.



**Size chart for styles 915 & 785**

Size	XS	S	M	L	XL	1X	2X	3X
Bust	32-34	35-37	38-40	41-43	44-46	47-50	51-54	55-58
Waist	29-30	31-32	33-34	35-37	38-40	41-43	44-47	48-51
Hip	34-36	37-39	40-42	43-45	46-48	49-52	53-56	57-60

Size chart is based on actual measurements in inches and not bra size.  
 If measurement falls in between sizes, always size down.  
 When determining size, fit breast/chest measurement first.

# Compression Crop Top 785

**Compression for underarm (axilla), upper back and breast/chest.**


This top is a short version of our popular Compression T with the addition of a 1 1/2" soft band under the bust for support. An alternative for the woman who may be hard to fit in style 915, but still needs upper torso compression.

**Features:**

- Sleeves for additional coverage
- Band provides support for breast forms
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Nude ○





COMPRESSION +  
SWIMWEAR UPF 50+



*Tops 780, Bottoms 614S & 615*





## Sydney Bra 780

**Compression for underarm (axilla) and upper back.**

The Sydney Bra by Wear Ease® is designed for women to wear as a compression bra. The wide band support and the surplice style of the bust provides a great bra to wear while managing tissue expanders. Can be worn as a swim top or an activity bra.

**Features:**

- Removable bust cups for additional shape
- Wide band for exceptional support
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Fuchsia ● Nude ○ White ○



**Size chart for styles 780 & 970**

Size	XS	S	M	L	XL	1X	2X	3X
Bust	31-33	34-36	37-39	40-42	43-45	46-48	49-51	52-54
Waist	26-28	29-31	32-33	34-36	37-39	40-42	43-45	46-48
Hip	34-36	37-39	40-42	43-45	46-48	49-51	52-54	55-57

Size chart is based on actual measurements in inches and not bra size.  
 If measurement falls in between sizes, always size down.  
 When determining size, fit breast/chest measurement first.

## Crisscross Shaper 970

**Compression for underarm (axilla), back and abdomen.**


Our Crisscross Shaper is a pretty alternative to a compression bra. It features a surplice neckline and a double layer of supportive fabric from the bra line to the bottom front.

**Features:**

- Removable bust cups for smoothness and balance
- Criss-cross style flatters the figure
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Colors:** Black ● Ivory ○



COMPRESSION +  
FULL CONTAINMENT



*Tops 914, Bottoms 614S & 615*

# Compression Shorts & Capris

Ideal for anyone who needs a great fit and comfortable compression. Frequently layered with bandages, wraps or stockings. Whether you prefer shorts or capris, the "High Waist" compression style provides the perfect combination of comfort and containment to help you look good and feel great. UPF 50+.

*These compression bottoms are excellent for women with fuller figures, provides good tissue containment, are easy to get on and off, and can be worn underneath jeans or pants. Our clinic has experienced excellent compliance with these garments.*

*Christine Rice, OTR/L, CLT,  
North Star Cancer Care, Yakima*



## High Waist Compression Short 614S

3" shorter than style 614L

Size	S	M	L	XL	1X	2X	3X	4X	5X	6X
Waist	31-33	34-36	37-39	40-44	45-49	50-54	55-60	61-66	66-72	73-78
Hip	35-36	37-39	40-42	48-52	53-58	59-64	65-70	71-76	77-82	83-88
Front Rise	13½	14¾	15	15½	16	16½	17	18	20	22
Back Rise	15	16	17	18	19	20	21	22	23	24
Inseam	8¼	8¾	9	9	9	9½	9½	10	13	13½

## High Waist Compression Short 614L

Size	S	M	L	XL	1X	2X	3X	4X	5X	6X
Waist	31-33	34-36	37-39	40-44	45-49	50-54	55-60	61-66	66-72	73-78
Hip	35-36	37-39	40-42	48-52	53-58	59-64	65-70	71-76	77-82	83-88
Front Rise	13½	14¾	15	15½	16	16½	17	18	20	22
Back Rise	15	16	17	18	19	20	21	22	23	24
Inseam	11¼	11¾	12	12	12	12½	12½	13	16	16½



## High Waist Compression Capri 615

Size	XL	1X	2X	3X	4X
Waist	40-44	45-49	50-54	55-60	61-66
Hip	48-52	53-58	59-64	65-70	71-76
Front Rise	15½	16	16½	17	18
Back Rise	17	19	20	21	22
Inseam	19¼	19¼	20	20½	21

Model is wearing size 1X in Styles 614 & 615.  
Height 5' 11", Waist: 45", Hip: 59", Thigh 31",  
Front Rise 15½", Back Rise 21½"

## Compression Short 612

Compression for waist, hips, buttocks and thighs.

Try our compression short for relief from swelling caused by lymphedema, lipedema, edema and post-surgical recovery. This short provides just the right amount of support and comfort for long term wear and compliance.

**Features:**

- Double layer of fabric for comfortable compression
- Hidden seams for no contact with skin
- UPF 50+

\*See page 4 for more features.

**Color:** Black ●



Size	S	M	L	XL
Waist	31-33	34-36	37-39	40-42
Hip	35-36	37-39	40-42	42-44



## Compression Capri 611

Compression for waist, abdomen, hips, buttocks, thighs and calves.

**Features:**

- Double layer of fabric for comfortable compression
- Hidden seams for no contact with skin
- UPF 50+

\*See page 4 for more features.

**Color:** Black ●

<sup>1</sup> Model is wearing size 1X in Styles 614 & 615. Height 5' 11", Waist: 45", Hip: 59", Thigh 31", Front Rise 15½", Back Rise 21½"

<sup>2</sup> Model is wearing size M. Height 5' 6", Waist: 31½", Hip: 42"

**Size chart for styles 611 & 613**

Size	S	M	L	XL	1X	2X	3X	4X
Waist	31-33	34-36	37-39	40-42	43-46	47-50	51-54	55-58
Hip	35-36	37-39	41-43	45-47	49-51	53-55	57-59	61-63



## Short Compression Capri 613

Compression for waist, abdomen, hips, buttocks and thighs.

**Features:**

- 6" shorter than style 611
- Double layer of fabric for comfortable compression
- Hidden seams for no contact with skin
- UPF 50+

\*See page 4 for more features.

**Color:** Black ●

<sup>3</sup> Model is wearing size M. Height: 5'7", Waist: 32", Hip: 40", Front Rise: 11", Back Rise: 15", Thigh: 21"

<sup>4</sup> Model is wearing size 1X. Height: 5'1", Waist: 42", Hip: 54", Front Rise: 12 1/2", Back Rise: 17", Thigh: 30"



## Compression Bodysuit 1000

Compression for underarm (axilla), arm, breast/chest, shoulder, back and abdomen.

Total torso coverage for a variety of procedures. Take our popular Compression T, style 915, and add a ¾ sleeve and a panty with a snap crotch. (If taller than 5'4" consider ordering style 1020 instead.)



**Features:**

- Fastens at crotch. Will not ride or roll up on fuller figures
- Comfortable ¾ sleeves
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Color:** Black ●


Size	XS	S	M	L	XL	1X	2X	3X
Bust	32-34	35-37	38-40	41-43	44-46	47-50	51-54	55-58
Waist	29-30	31-32	33-34	35-37	38-40	41-43	44-47	48-51
Hip	34-36	37-39	40-42	43-45	46-48	49-52	53-56	57-60

Size chart is based on actual measurements in inches and not bra size. If measurement falls in between sizes, always size down. When determining size, fit breast measurement first.

## Ellen Compression Bodysuit 1020

Compression for underarm (axilla), breast/chest, back and abdomen.

Provides comfortable compression for the total torso. This bodysuit can also be worn with compression shorts and capris for additional abdominal compression.



**Features:**

- Seamless front, no seam under the bust
- Fastens at crotch. Will not ride or roll up on fuller figures
- Double layer fabric for firm fit
- Pocketed for prostheses
- UPF 50+

\*See page 4 for more features.

**Color:** Nude ○

Size	S	M	L	XL	1X	2X
Bust	35-37	38-40	41-43	44-46	47-50	51-54
Waist	31-32	33-34	35-37	38-40	41-43	44-47
Hip	37-39	40-42	43-45	46-48	49-52	53-56

Size chart is based on actual measurements in inches and not bra size. If measurement falls in between sizes, always size down. When determining size, fit breast/chest measurement first.



# COMPRESSION+ COOL SUPPORT



Tops 790; Bottoms 614S, 614L & 612

# Compression Bras and Tops

## BY WEAR EASE®

Wear Ease Compression Bras and Tops are ideal for treating swelling from lipedema, lymphedema, edema and surgery. Constructed of a double layer of breathable and moisture management fabric, developed for medical compression garments, 83% nylon/17% lycra. A snug and comfortable fit will result in mild compression, 10-20 mmHg.

**\*\*Features:**

- Cushioned front zipper closure
- Adjustable and detachable hook and eye shoulder straps
- Soft "latex free" elastic band encircles torso, provides support and allows attachment of drainage tube pouches (styles 791, 797, 799)
- High cut underarm for maximum compression on the axilla
- Worn comfortably under clothing for 100% compliance



Style 770



Style 790



Style 796



Style 950



Style 951



Style 960/962



Style 798



Style 961/963



Style 953



Style 954



# Compression Bra

## 790/791/792

### Compression for underarm (axilla), breast and upper back

Designed to provide compression for the upper torso. Double layer bra with high sides and built-in slings that lift the bust. Great coverage for the underarm area (axilla) where fluid tends to collect.

#### Features:

- Ready to wear and available in a wide range of sizes
- Cool and soft to the skin for continuous comfort
- Inside pockets to hold and keep breast forms in place
- Soft "latex free" elastic band encircles the torso, provides excellent support and allows attachment of drainage tube pouches (style 791)
- Adjustable shoulder straps for ease of fit
- Easy and flattering fit for long term wear

\*\*See page 18 for more features.

**Style 790:** Compression Bra

**Style 791:** 790 + a pair of removable drainage tube pouches

**Style 792:** 791 + a pair of fiberfill breast forms

**Colors:** Black ● Latte ● White ○

**Shoulder Strap Extenders**  
See page 30



A, B, C & D Cup Sizes	Size	XXS	XS	S	M	L	XL	1X	2X	3X	4X
	Bust (at fullest in inches)	29-31	32-34	35-37	38-40	41-43	44-46	47-49	50-52	53-58	59-64

***DD/DDD Cup Sizes	Size	SD	MD	LD	XLD	1XD
	Bust (at fullest in inches)	38-39	40-41	42-43	44-47	48-51

See page 21 for detailed fitting.

## Sternum & Genital Pad

### BY WEAR EASE® 309

Insert under the bra to treat swelling and fibrosis of the sternum and the inner breast area. Shown above tucked inside Wear Ease Compression Bra, Style 790. May also be used as a genital pad.



*I was really happy to see that you all have a sternum pad. When I had my surgery all of the post-surgical bras (not Wear Ease bra) hurt my sternum and it was very uncomfortable to sleep and sit up throughout the day. I can definitely let customers know how useful it will be after surgery.*



Sonya Evans, Love Is Pink

“ Thank you so much for the sample. I tried it on, and it is one of the nicest bras I have worn. I am recovering from breast reconstruction and have tried all the manufacturers, and this is my favorite one. I feel like it was made for me. I am going to order more today for myself. This one will go in our fitting room for our patients.

Terri McClanahan,  
Owner Hope 4 Healing, Inc.

## Sizing and Fitting Compression Bra 790/791/792

### How To Determine Your Compression Bra Size

#### 1. BUST SIZE

Measure circumference at the fullest part of the bust.  
Measure with a bra on, and prosthesis if applicable.

#### 2. BAND SIZE

Measure the circumference directly under the bust.



#### 3. CUP SIZE

Subtract the band size from the bust.

2" or less = A/B = regular 790

2" - 4" = C/D = regular 790

4" - 6" = DD/DDD = 790D\*\*\*



### Fitting the Compression Bra

1. When you receive your new Compression Bra by Wear Ease® take it out of the package and unzip the front zipper all the way. If wearing a prosthesis, slip it into the pocket inside the bra.
2. Slip the bra on as you would any front-closure bra, attach the hook and eyes inside then zip up the zipper.
3. Lean over, adjusting your breasts comfortably in the bra. Be sure that the band fits comfortably under your breasts, supporting them. The bra should fit snug like a sports bra.
4. The bra should fit high under your arms to provide full coverage of the axilla area. If the bra fits too tight under the arms you can lengthen the bra straps by adjusting the hook and eye row on the shoulder straps. **You can also attach our Shoulder Strap Extender, Style 400, for additional length. See inset on page 20.**



COMPRESSION+  
TEXTURED MASSAGE



Tops 770; Bottoms 614S, 614L & 612



# Karena Textured Compression Bra 770

The Karena Compression Bra from Wear Ease is a revolutionary product for comfortable, all-day support. The unique textured design also provides compression with micromassage to help reduce swelling and soften radiation fibrosis and scarring.

**NEW!**



inside out

**Features:**

- Zipper front allows for easy donning and doffing
- Adjustable straps enables fitting various sizes
- Breathable and lightweight fabric makes it comfortable to wear day and night
- Double thickness of fabric (placket) provides coverage for the sternum and separation of the breasts
- Textured lining helps soften scars and radiation fibrosis at the breast, thorax and axilla. It also helps new scars heal with better results.

**Note:** *this garment is not appropriate to use immediately after surgery because the micromassage action may be too aggressive on skin. Wait at least 3-4 weeks until surgical and radiation wounds are healed before wearing.*

- Combination of compression and massaging provides relief from pain, swelling, inflammation and overall discomfort. Helps massage the inflamed area while compressing it for more comfortable healing

**Color:** Beige ●

Size	S	M	L	XL	1X
Bust	34-36	37-39	40-42	43-45	46-48

*I have many breast cancer patients who struggle with the fibrosis aftermaths of treatment: surgical scarring and radiation fibrosis. This bra is lined with a micromassage fabric that provides treatment to help soften the fibrosis of the breast and thorax as well as being a great compression garment. I originated the concept by sewing fabric from another micromassage compression garment into my patients' bras: Wear Ease has taken this bra a step further with a truly wonderful product!*

Karen Ashforth, MS, OTR, CLT-LANA



## Andrea Shirt 960/962

**Compression for underarm (axilla), breast/chest, abdomen, back and hips.**

A new way to wear compression pads. Sewn in axilla pads eliminate the hassle of keeping pads from moving around. May be purchased with 2 pads or 1 pad (right pad only or left pad only). Or no pads, style 962.

**Features:**

- Set in sleeves for additional coverage
- Sewn in pads provide compression for lateral chest wall under the axilla

\*\*See page 18 for more features.

**Axilla Pad Options**

- Style 960: Both
- Style 960R: Right
- Style 960L: Left
- Style 962: No pads

**Color:** White ○



Style 960 (inside out)

Size	XS	S	M	L	XL	1X	2X	3X	4X
Bust	31-33	33-35	35-37	37-39	39-41	41-44	44-47	47-50	50-53
Waist	33-35	35-37	37-39	39-40	40-42	42-45	45-48	48-51	51-54
Hip	36-38	38-40	40-42	42-44	44-46	46-49	49-52	52-55	55-58

Size chart is based on actual measurements in inches. For help with sizing please call customer support 866-251-0076.



## Torso Compression Vests 950 & 951

Compression for underarm (axilla), breast/chest, abdomen, back and hips.

For women who need total torso compression. Ideal solution for the bilateral woman who cannot wear other compression garments due to rolling. Features include a scoop neck, zipper opening with hook and eye fasteners and soft placket underneath. **Style 951:** Features a V-neck

**Features:**

- "A shaped" top for figures fuller in the waist and hips
- Seamless body contact, inside and out for comfort
- Full length
- No pockets

\*\*See page 18 for more features.

**Colors:** Black (951) ● White (950) ○



Shoulder Strap Extenders  
See website

**Size chart for styles 950/951**

Size	XS	S	M	L	XL	1X	2X	3X	4X
Bust	31-33	33-35	35-37	37-39	39-41	41-44	44-47	47-50	50-53
Waist	33-35	35-37	37-39	39-40	40-42	42-45	45-48	48-51	51-54
Hip	36-38	38-40	40-42	42-44	44-46	46-49	49-52	52-55	55-58

Size chart is based on actual measurements in inches and not bra size. The size chart is meant to be a guideline. When determining size, fit breast/chest measurement first.



## Women's Compression Vest 796/797

Compression for underarm (axilla), breast/chest and upper back.

Designed for compression all around the torso for women with post-surgical discomfort, breast lymphedema or edema. The vest is recommended for women who may need more compression than provided by our Compression Bra, style 790. The 4-way stretch of this power mesh provides the fit and compression for the torso being fitted. This feature makes this garment an ideal solution for bilateral surgery or uneven breasts.

**Features:**

- Soft "latex free" band encircles torso, provides support and allows attachment of drainage tube pouches (**style 797**)
- **No cups or pockets for prostheses**

**Style 796:** Compression Vest

**Style 797:** 796 + 2 removable drainage tube pouches

**Colors:** White ○ Latte ●

\*\*See page 18 for more features.



Shoulder Strap Extenders  
See page 30

**Size chart for styles 796/797**

Size	XXS	XS	S	M	L	XL	1X	2X	3X
Bust	30-31	32-33	34-35	36-38	39-41	42-44	45-47	48-50	51-53

Size chart is based on actual measurements in inches and not bra size. The size chart is meant to be a guideline. If measurement falls in between sizes, always size down. When determining size, fit breast/chest measurement first.

# COMPRESSION+ FOR ALL SHAPES & SIZES



Top 916; Bottoms 612

# Men's Compression Vests and Shirts

Styles 961/963, 953/954, 798/799

Ideal for managing swelling after chest procedures (male breast cancer surgery, male chest reconstruction, gynecomastia treatment and thoracic surgery). Featuring a front zipper with hook and eye and comfortable placket underneath to protect the skin from the zipper.

**Styles 798/799, 953/954, 961/963 feature:**

- High cut underarm for maximum compression on the axilla
- Seamless body contact, inside and out
- Latex Free



## Andrew Shirt 961/963

**Torso Compression for underarm (axilla), chest, back and abdomen.**

Full torso compression and upper arm support. Sewn in pads provide extra compression and coverage for the axilla. No need for additional fasteners to hold pads in place. Model is also wearing Style 611 Capris.

### Axilla Pad Options

**Style 961:** Both

**Style 961R:** Right

**Style 961L:** Left

**Style 963:** No pads

\*\*See page 18 for more features.

**Color:** White

Size	XS	S	M	L	XL	1X	2X	3X	4X
Chest	31-33	33-35	35-37	37-39	39-41	41-44	44-47	47-50	50-53
Waist	33-35	35-37	37-39	39-40	40-42	42-45	45-48	48-51	51-54
Hip	36-38	38-40	40-42	42-44	44-46	46-49	49-52	52-55	55-58

Size chart is based on actual measurements in inches.

## Eric Compression Shirt 916

**Compression for the underarm (axilla), back, chest and abdomen.**

A stylish masculine crew neck shirt that provides mild and comfortable compression for the chest, shoulders and overall torso. For post surgery recovery, relief from swelling and comfortable support.

### Features:

- Snug fit for consistent and full torso compression
- Double layer of fabric for strong support
- UPF 50+

\*See page 4 for more features.

**Color:** Black

Size	XS	S	M	L	XL	2XL	3XL	4XL
Chest	34-36	38-40	42-44	46-48	50-52	54-56	58-60	62-64
Waist	30-32	34-36	38-40	42-44	46-48	50-52	54-56	58-60
*Length	26.5	27	27.5	28	28.5	29	29.5	30

\*Measure neck to hemline.





# Men's Compression Shorts & Capris

The High Waist Compression Shorts are perfect for anyone needing a great fit and comfortable compression. The two lengths, 614S (short) and 614L (long), provide a variety of options depending on your preference. The 614L provides an extra 3" in length, in addition to a wider size range, for a better fit. Whether you prefer short or long compression shorts capris, the High Waist Compression Shorts provide the perfect combination of comfort and containment to help you look good and feel great.



## High Waist Compression Short 614S

3" shorter than style 614L

Size	S	M	L	XL	1X	2X	3X	4X	5X	6X
Waist	31-33	34-36	37-39	40-44	45-49	50-54	55-60	61-66	66-72	73-78
Hip	35-36	37-39	40-42	48-52	53-58	59-64	65-70	71-76	77-82	83-88
Front Rise	13½	14¾	15	15½	16	16½	17	18	20	22
Back Rise	15	16	17	18	19	20	21	22	23	24
Inseam	8¼	8¾	9	9	9	9½	9½	10	13	13½



## High Waist Compression Short 614L

Size	S	M	L	XL	1X	2X	3X	4X	5X	6X
Waist	31-33	34-36	37-39	40-44	45-49	50-54	55-60	61-66	66-72	73-78
Hip	35-36	37-39	40-42	48-52	53-58	59-64	65-70	71-76	77-82	83-88
Front Rise	13½	14¾	15	15½	16	16½	17	18	20	22
Back Rise	15	16	17	18	19	20	21	22	23	24
Inseam	11¼	11¾	12	12	12	12½	12½	13	16	16½

## Compression Short 612

Compression for waist, hips, buttocks and thighs.

Try our compression short for relief from swelling caused by lymphedema, lipedema, edema and post-surgical recovery. This short provides just the right amount of support and comfort for long term wear and compliance.

### Features:

- Double layer of fabric for comfortable compression
- Hidden seams for no contact with skin

\*See page 4 for more features.

Color: Black ●

Size	S	M	L	XL
Waist	31-33	34-36	37-39	40-42
Hip	35-36	37-39	40-42	42-44





## Men's Torso Compression Vests 953 & 954

Torso Compression for underarm (axilla), chest, back and abdomen.

**Features:**

- Detachable and adjustable shoulder straps.

**Style 954:** Features a V-neck

\*\*See page 18 for more features.

**Colors:** Black (954) ● White (953) ○

Size	XS	S	M	L	XL	1X	2X	3X	4X
Bust	31-33	33-35	35-37	37-39	39-41	41-44	44-47	47-50	50-53
Waist	33-35	35-37	37-39	39-40	40-42	42-45	45-48	48-51	51-54
Hip	36-38	38-40	40-42	42-44	44-46	46-49	49-52	52-55	55-58

Size chart is based on actual measurements in inches.



**Shoulder Strap Extenders**  
See website



## Men's Compression Vest 798/799

Compression for underarm (axilla), chest and upper back.

Short length ends right above the navel which eliminates having to fit the abdomen.

**Style 799:** Packaged with pair of drain tube pouches to accommodate surgical drains.

\*\*See page 18 for more features.

**Color:** White ○

Size	XXS	XS	S	M	L	XL	1X	2X	3X
Bust	30-31	32-33	34-35	36-38	39-41	42-44	45-47	48-50	51-53

Size chart is based on actual measurements in inches. If measurement falls in between sizes, always size down.



**Shoulder Strap Extenders**  
See page 30

# Accessories

## Bust Cups 514/516

Lightweight fabric-covered oval and triangular shaped molded foam bra cups used for modesty and to establish symmetry. Can be used to shape fiberfill leisure form, style 510.

### Benefits:

- Support for swimwear, dresses & formalwear
- Can be cut to desired shape without fraying
- Can be worn in any Wear Ease Bra or Camisole
- Machine washable & dry cleanable

**Colors:** Black ● Nude ○

### Packaging Options:

- 514: 1 pair per pack
- 514Z: 1 dozen (one size per dozen)
- 516: 1 pair per pack
- 516Z: 1 dozen (one size per dozen)



### Triangle Bust Cups, Part 514

#### Sizes:

- |            |              |
|------------|--------------|
| - 6 (XS-S) | - 12 (XL)    |
| - 8 (M)    | - 14 (1X)    |
| - 10 (L)   | - 16 (2X-3X) |



### Oval Bust Cups, Part 516

#### Sizes:

- |       |       |       |
|-------|-------|-------|
| - 30B | - 36B | - 38C |
| - 32B | - 38B | - 40C |
| - 34B | - 36C | - 40D |



## Fiberfill Breast Forms 510

Wear Ease fiberfill forms are triangular shaped and may be easily adjusted to any specific size and shape by adding or removing the fiberfill. Constructed of polyester and spandex microfiber, they are filled with 100% hypoallergenic polyester fluff. They may be worn in any Wear Ease Post-Surgical bra or camisole and are soft enough to be worn after surgery. Cool option for summer and leisure wear. Can be worn with oval and triangle bust cups. (See styles 514 and 516)

### Benefits:

- Restores shape without discomfort
- Provides additional padding over surgical site
- Triangular shape
- Easily adjusted to become any specific size required
- Sold by the pair

**Colors:** Black ● Nude ○ White ○

**Sizes:** S, M, L, XL

## Shoulder Strap Extenders 400

Shoulder Strap Extenders are designed to add length to the shoulder straps.

**Use with Styles 790/796/798**

### Features:

- One pair per package
- Adds 1" - 3" to shoulder straps

**Colors:** Black ● Nude ○ White ○



## Solaris Collection™ SwellSpots and Wear Ease Compression Bras, Camisoles & Pads

SwellSpots, part of the Solaris Collection by L&R USA INC., are professionally designed edema and fibrosis pads available in a wide range of versatile shapes and sizes that tuck into compression garments or under compression bandages.

**Designed to address swelling and/or fibrosis in the breast and chest region.**

Below are some examples of how various swellsots can be used with Wear Ease garments.



People.Health.Care. **SwellSpots™**



### Curved Pad SP-01

Placed under the arm or inside the bra to aid in reducing swelling. Can be used to direct lymph flow around scars in any areas of the body.



Shown with Wear Ease Compression Bra, Style 790

### Breast Pad SP-09

Inserted under the bra to treat swelling and fibrosis of the breast area due to a lumpectomy or reconstructive surgery.



Shown with Wear Ease Compression Bra, Style 790

### Circular Pad SP-06

Recommended to treat fibrotic tissue on the chest wall due to radiation and or surgical incision scarring.



Shown with Wear Ease Slimmer, Style 910

### Lateral Bra Pad SP-16-A

Tucked inside the compression bra pockets, regular bra or camisole to reduce scarring and protect delicate areas from strap pressures.



Shown with Wear Ease Slimmer, Style 910

### Kidney Pad SP-08

Placed under the arm to aid in reducing the swelling in the underarm area (axilla).



Shown with Wear Ease Compression Crop Top, Style 785

### Full Bra Pad SP-16-B

Selected for similar use as the SP-16-A above but for front opening bras.



Shown with Wear Ease Compression Bra, Style 790



### Axilla Pad by Wear Ease 305

Shown with Wear Ease Katy T, Style 918



### Sternum & Genital Pad by Wear Ease 309

Insert under the bra to treat swelling and fibrosis of the sternum and the inner breast area. Shown above tucked inside Wear Ease Compression Bra, Style 790.

May also be used as a genital pad

**Consider several pads due to perspiration and hygiene. Wash one, wear one.**



# Post-Surgical Camisoles

BY WEAR EASE®

Wear Ease offers a unique style of post-surgical camisoles with built-in support of a shelf bra. Our camisoles have different donning and fastening methods depending on the needs of the patient and the requirements of the surgeon. These camisoles can be worn after a wide range of surgical procedures and treatments including: mastectomy, lumpectomy, biopsy, breast augmentation, breast reduction, breast reconstruction, radiation treatment and for transition during mastectomy recovery.

**Features:**

- Built-in shelf bra provides under bust support
- Pocketed for breast forms if needed
- Fine, no-pill cotton, 95% cotton, 5% spandex
- Packaged with a pair of detachable drainage tube pouches
- Optional: packaged with a pair of fiberfill breast forms
- Latex free

**Fitters notes:**

- Inside soft, latex free elastic encircles the entire torso allowing positioning of the removable drainage tube pouches from front, sides or back according to needs. The drainage tube pouches can be removed and the camisole can be worn for leisure wear after surgery, during radiation and recovery.
- If surgery requires more than one pair of pouches, additional pairs can be purchased under part 410.
- Should fit comfortably but not too snug. Need room for holding and hiding the drain pouches.



## Beth Camisole 903/904

Beth offers all the features and benefits of a post-surgical bra with the additional security of a camisole.


**Features:**

- Built-in shelf bra
- Soft zipper front closure with protective placket
- Adjustable straps for easy fitting

**Style 903:** includes a pair of removable drainage tube pouches

**Style 904:** 903 + a pair of fiberfill breast forms

**Color:** White ○



Size	XS	S	M	L	XL	1X	2X	3X	4X
Bust	34	36	38	40	42	45	48	51	54
Waist	26	28	30	32	34	37	40	43	50
Hip	36	38	40	42	44	47	50	53	58





## Dawn Camisole 900/901

Dawn is a post-surgical camisole that does not look surgical. With a little bit of lace and choice of colors (white, nude, black), women appreciate wearing this attractive and supportive top during recovery. Lingerie satin elastic and stretch lace add to the feminine feel. This no-closure camisole is easy to slide on and off (step in/step out). Often worn for years after surgery as a replacement for a bra.

**Features:**

- Built-in shelf bra
- Color-matched accessories

**Style 900:** includes a pair of removable drainage tube pouches

**Style 901:** 900 + a pair of fiberfill breast forms

**Colors:** Black ● Nude ○ White ○



## Drain Tube Pouches 410

**Benefits:**

- Eliminates the need for lanyards and pins to hold surgical drains
- Attaches to soft elastic band in Wear Ease bras and camisoles
- Helps to control pain by easing stress or pulling on incision
- Convenient and secure
- One size, one pair per package
- Some procedures may require more than one pair

**Colors:** Black ● Nude ○ White ○

Size	XS	S	M	L	XL	1X	2X	3X
Bust	34	36	38	40	42	45	48	51
Waist	26	28	30	32	34	37	40	43
Hip	36	38	40	42	44	47	50	53

# Post-Surgical Bras

## BY WEAR EASE®

Wear Ease has designed post-surgical soft cup bras to meet the needs of women through recovery and transition. Our bras are designed with attractive gathering but with different closures and adjustability. Designed to be worn after a wide range of surgical procedures and treatments, including: Mastectomy, Lumpectomy, Biopsy, Radiation Therapy, Breast Augmentation, Breast Reduction, Breast Reconstruction and Thoracic Surgery.



## Allyson Bra 750/751/752

A supportive bra with zipper front closure and adjustable straps. Comfortable choice to wear during transition, mastectomy recovery, daily wear and as a sleep bra.

**Features:**

- Adjustable shoulder straps for easy fitting
- Wide band for exceptional support
- Pocketed and can be worn as a mastectomy bra
- Comfortable and soft cotton spandex fabric, 95% cotton, 5% spandex

**Style 750:** Allyson Bra

**Style 751:** 750 + a pair of removable drainage tube pouches

**Style 752:** 751 + a pair of fiberfill breast forms

**Colors:** Black ● Nude ○ White ○

Size	XS	S	M	L	XL	1X	2X	3X
Bust	34	36	38	40	42	45	48	51
Waist	26	28	30	32	34	37	40	43





## Grace Bra 740/741/742

Grace Post-Surgical Bra is a front closure bra that is extremely comfortable and is one of our most popular specialty bras. It is also comfortable for daily wear and as a sleep bra.

**Features:**

- Cushioned front zipper closure
- Comfortable and gentle moisture wicking microfiber, 92% polyester, 8% spandex
- Wide band for exceptional support
- Pocketed and can be worn as a mastectomy bra

**Style 740:** Grace Bra

**Style 741:** 740 + a pair of removable drainage tube pouches

**Style 742:** 741 + a pair of fiberfill breast forms

**Colors:** Black ● Nude ● White ○

**Size chart for styles 740/741/742 & 700/701/702**

Size	XS	S	M	L	XL	1X	2X	3X
Bust	34	36	38	40	42	45	48	51
Waist	26	28	30	32	34	37	40	43



## Compression Bra 790/791/792

**Compression for underarm (axilla), breast and upper back**

Designed to provide compression for the upper torso. Double layer bra with high sides and built-in slings that lift the bust. Great coverage for the underarm area (axilla) where fluid tends to collect.

**Style 790:** Compression Bra

**Style 791:** 790 a pair of removable drainage tube pouches

**Style 792:** 791 a pair of fiberfill breast forms

**Colors:** Black ● Latte ● White ○

See page 19 for features and page 20 for size chart.



# Prescription Form

Date \_\_\_\_\_

Prescriptions for upper body compression and diagnosis of lymphedema in the surgical quadrant.



Wear Ease, Inc.  
5491 W Kendall St.  
Boise, ID 83706  
866.251.0076  
www.wear ease.com

## TOPS

Name	Compression T	Slimmer	Katy T	Crisscross Shaper	Taylor T	Ava Compression Camisole	Compression Camisole	Compression Bodysuit
Style #	915	910	918	970	914	911	912	1000
Shoulder	•		•		•			•
Underarm (axilla)	•	•	•	•	•	•	•	•
Arm			•		•			•
Breast/Chest	•	•	•		•	•	•	•
Upper Back	•	•	•	•	•	•	•	•
Back	•	•	•	•	•	•	•	•
Abdomen	•	•	•	•	•	•	•	•

Name	Ellen Comp. Bodysuit	Sydney Bra	Crop Top	Compression Bra	Karena Compression Bra	Compression Vest	Torso Comp. Vest	Torso Comp. Vest	Andrea Shirt
Style #	1020	780	785	790, 791, 792	770	796, 797	950	951	960, 962
Shoulder			•						•
Underarm (axilla)	•	•	•	•	•	•	•	•	•
Breast/Chest	•	•	•	•	•	•	•	•	•
Upper Back	•	•	•	•	•	•	•	•	•
Back	•								•
Abdomen	•								•
Hips									•

Patient Name \_\_\_\_\_ Date of Birth \_\_\_\_\_  
 Dr's Name \_\_\_\_\_ Dr's Sig \_\_\_\_\_ NPI# \_\_\_\_\_  
 (Dx) \_\_\_\_\_ Written Description: \_\_\_\_\_



MADE IN USA



# Prescription Form

Date \_\_\_\_\_

Prescriptions for lower body compression.

Prescriptions for upper body compression and diagnosis of lymphedema in the surgical quadrant.



Wear Ease, Inc.  
5491 W Kendall St.  
Boise, ID 83706  
866.251.0076  
www.wearase.com

## BOTTOMS

Name	Compression Capri	Compression Short	Short Comp. Capri	High Waist Comp. Short	High Waist Comp. Short	High Waist Comp. Capri
Style #	611	612	613	614S	614L	615
Abdomen	•	•	•	•	•	•
Hips	•	•	•	•	•	•
Buttocks	•	•	•	•	•	•
Thighs	•	•	•	•	•	•
Calves	•	•	•	•	•	•



## MEN'S

Name	Chest Comp. Vest	Torso Comp. Vest	Torso Comp. Vest	Andrew Shirt	Eric Shirt
Style #	798, 799	953	954	961, 963	916
Shoulder	•	•	•	•	•
Underarm (axilla)	•	•	•	•	•
Breast/Chest	•	•	•	•	•
Upper Back	•	•	•	•	•
Back	•	•	•	•	•
Abdomen	•	•	•	•	•
Hips	•	•	•	•	•



Patient Name \_\_\_\_\_ Date of Birth \_\_\_\_\_

Dr's Name \_\_\_\_\_ Dr's Sig \_\_\_\_\_ NPI# \_\_\_\_\_

(Dx) \_\_\_\_\_ Written Description: \_\_\_\_\_



MADE IN USA



5491 W. KENDALL ST.  
BOISE, ID 83706

# READY TO WEAR

COMPRESSION & POST-SURGICAL



FIND A RETAILER:

Visit [wearease.com](http://wearease.com) >>> store locator

BECOME A DEALER:

1.866.251.0076 or  
[wearease.com](http://wearease.com) >>> for wholesalers



866.251.0076  
[www.wearease.com](http://www.wearease.com)

 MADE IN THE USA

7\_23\_4k