

00:01:00

ONE MINUTE PRODUCT MEETING

LuxaFlow® Ultra

LuxaFlow Ultra is the ideal add-on resin for use with Luxatemp. Specifically designed for excellent compatibility with Luxatemp Ultra and other bisacryl provisional materials, it is formulated with fluorescence to yield a completely natural looking result in all lighting conditions.

Delivered from a high precision syringe tip dispensing system, LuxaFlow Ultra provides excellent flowability in a non-slumping formula. It easily blocks out undercuts, repairs enamel defects, voids, margins, and seals gingival margin interfaces.

It is also perfectly suited for small class III, IV and V restorations and minimally invasive restorations.



FEATURES & BENEFITS

- › Natural-looking results
- › Add-on resin for bisacryl provisionals
- › Optimum flowability
- › Non-slumping formula
- › Light-cured
- › High-precision, syringe-tip delivery
- › Precise provisional color matching
- › Radiopaque

PROCEDURES

- › Class III, IV and V restorations

COMPANION PRODUCTS

- › Luxatemp product family

ITEM # & DESCRIPTION

- 224000 | Intro Kit
- 224001 | A1-Refill Package
- 224002 | A2-Refill Package
- 224003 | A3-Refill Package
- 224004 | A3.5-Refill Package
- 224005 | B1-Refill Package
- 211759 | Syringe Tips (Bag of 20)

LOOKING FOR MORE INFO?

Visit www.dmg-america.com/products.