



# SALE

2021 3rd PROMOTION • Expires September 30, 2021

EN: 1.800.667.8811 • FR: 1.800.211.1200

## Everyday Dentistry Made Better

Products Designed to Improve Patient Care and Comfort



Please visit [www.biscocanada.com](http://www.biscocanada.com)  
to take advantage of more products on sale!  
Receive **FREE SHIPPING** on all online orders.



# VacuLUX

## HVE ISOLATION SYSTEM

NEW  
PRODUCT

Introducing the **VacuLUX HVE** Isolation System.  
Better patient comfort, field and aerosol management.



### Introductory Offer

**BUY STARTER KIT  
GET MOUTHPIECE 4/PK FREE**

**\$575.00**

#### VacuLUX Isolation Starter Kit

- 1 HVE Adapter
- 1 Whisper Lite Kit
- 4 assorted autoclavable mouthpieces (Pedo, Sm, Med, and Lg)

**10% OFF**

**\$133.20 ~~\$148.00~~**

**VacuLUX Mouthpiece 4 pack**  
AVLUXM (Med), AVLUXS (Sm)  
AVLUXP (Pedo), AVLUXL (Lg)  
AVLUXA (Asst)

**10% OFF**

**\$292.50 ~~\$325.00~~**

**VacuLUX Adapter**  
VLUXD (Dryshield)  
VLUXZ (Zyris)

Autoclavable mouthpieces are now available and are universally compatible with Ascentcare's **VacuLUX** adapters as well as **Zyris Isolite** or **Dryshield** systems.

### Discover the features worth making a switch

In addition to a bite block, cheek retractor, tongue guard and hands-free, continuous high-volume suction, enjoy the added benefits of features such as:



#### Universal Fit

Designed to fit Zyris or Dryshield systems



#### Autoclavable

Save money with a more affordable option



#### Patient Comfort Flap

Soft silicone flap protects for a more comfortable experience



#### Liquid Rollover Edge

Prevent escaping fluids for increased spatter reduction



#### Evacuation Chamber

Consistent, predictable suction for a dry field you can count on



#### Integrated Light Guide

Guided, targeted illumination for VacuLUX or Isolite systems





# HYPERLUX



**Introductory Offer**  
**\$549.00 ~~\$589.00~~**

**HyperLUX (HLSTM)**

**- Portable Fiber Optic Light Engine:**

- HyperLUX light engine
- Light engine stand with docking ports
- Charging Cable and Adapter
- 5 foot Magnetically Coupled Light Guide Cable

The HyperLUX™ Fiber Optic light engine provides intense, cool illumination. Interfacing with a fiber optic light guide to direct high-powered light to Ascentcare's innovative dental instruments.



## Why the VacuLUX™ for illumination



### Plug & Play

Magnetically couple the fiber optic light guide for instant illumination without installation.



### Universal

Choose an adapter to illuminate your Dryshield™ or Zyris™ Isolite® & VacuLUX™ mouthpieces.



### Upgradable

Upgrade to illumination at any time as each VacuLUX adapter includes an integrated light guide.



### Portable

Share illumination at every chair with the portable, rechargeable light engine.



### Powerful

Choose from three intensity levels to send illumination where it's needed most with the HyperLUX™.

# VacuVUE

HVE DENTAL MIRROR



Introductory Offer

10% OFF

**\$943.20** ~~\$1048.00~~

VacuVUE HVE  
Mirror 3 Pack Starter Kit  
4VVUES3 (Size 4) or 5VVUES3 (Size 5)

Reduce aerosols and spatter by 90% or more with the **VacuVue HVE Mirror**.

## PureHD

Superior Mirror Technology:

- 99% color accuracy
- 40% higher reflective value
- front mirroring: 75 layers
- **replaceable mirror faces**




## Liquid Control

Reduces **spatter** and pooling fast with dual-sided liquid evacuation ports.



## Aerosol Control

Engineered with a large-bore port to help you reduce aerosols using continuous suction.



## ALUMINUM CONSTRUCTION

Constructed with aerospace-grade aluminium, this **autoclaveable** mirror is durable and lightweight.

## Ergonomic Handle

Triangular shape offers smooth, **agile** maneuverability.

## Simple Installation

Easily connect to your current valve or use the lightweight, **kink-resistant** Whisper Lite kit.



## 360° rotation

Designed to allow for **unhindered** rotation without losing the suction seal.






# QuietOn

## Active Noise Cancelling Earplugs

### Speech-Friendly Hearing Protection, Optimized For Dental Professionals



### How does it work?

- QuietOn earplugs combine active and passive noise cancellation technology
- The inner-device **microphone** "listens" to the noise
- The inner-device **speaker** generates opposite-phased soundwaves, cancelling the noise
- Attenuation varies by different frequencies
- Significant reduction of low and high frequencies
- Mid-frequencies like human speech reduced only slightly - important for staff-patient communication

QuietOn performance			
	LOW FREQUENCIES	MID FREQUENCIES	HIGH FREQUENCIES
QuietOn's noise reduction with foam tip	38.4 dB	27.4 dB	36.7 dB
QuietOn's noise reduction with silicone tip	39.1 dB	24.4 dB	31.1 dB
<b>NOISE TYPE</b>			
Speech			
Sirona Cerec MC XL			
Autoclave			
Assistina			
Compressor			
Vacuum pump			
Turbine (alone)			
Micromotor (alone)			
Ultrasonic scaler (alone)			
Light volume suction			
Saliva ejector			
3 way syringe (water & air)			
Turbine & High volume suction			
Micromotor & High volume suction			
Ultrasonic scaler & High volume suction			

### QuietOn's performance across the dental noise spectrum

In actual clinical situations, simultaneous noise sources can occur at different frequencies.

**The total noise level** is higher and spread across a wider spectrum, compared to noise from a single source.

**\$259.00** ~~**\$309.00**~~

**\$50 OFF**

### QuietOn Dental package

- 1 pair of QuietOn earplugs in a charging case
- Retention cord
- Micro USB cable
- User guide
- 1 pair of QuietOn Comply foam tips (Medium) attached to the earplugs
- 8 pairs of silicone tips in 2 sizes: 4 Blue (Large), 4 Green (Small)



# ALL-BOND UNIVERSAL™

Light-Cured Dental Adhesive

**A TRULY UNIVERSAL ADHESIVE**  
offering optimized performance for all of  
your daily restorative needs.



## One bottle. Countless Advantages.

### 1 COMPATIBLE

Compatible with all light-, self- and dual-cured resin composite and cement materials for all direct and indirect procedures

### 2 VERSATILE

Its chemistry offers the flexibility to either total-, self- or selective-etch

### 3 DURABLE

Offers high bond strength to all indirect substrates including: metal, glass ceramics, zirconia, alumina, porcelain (silica-based) and lithium disilicate

### 4 ONE BOTTLE

No additional activator is required

### 5 OPTIMAL

Low film thickness (less than 10 microns) and virtually no post-operative sensitivity, without sacrificing clinical performance



**BUY 4, Get 1 FREE\***

\*Free item must be the same item

**\$195.00** | Standard Kit (B-7202B0K-EU)

**\$184.00** | Refill (6ml) (B-7202P)

**\$219.00** | 50 Unit-Dose (0.20ml) (B-73050K)

**\$409.00** | 100 Unit-Dose (0.20ml) (B-73100K)

# ALL-BOND 3®

## The Gold Standard Just Got Better

ALL-BOND 3 Resin is **radiopaque allowing visibility on radiographs** in areas where adhesives commonly pool, reducing the risk of misdiagnosing caries



<b>\$299.00</b>	<b>ALL-BOND 3 KIT</b> (B-36200K)
<b>\$230.00</b>	<b>PART A &amp; B SET</b> (B-36010P)
<b>\$74.00</b>	<b>ALL-BOND 3 RESIN BOTTLE</b> (B-3608P)


**BUY 4, Get 1 FREE\***

\*Free item must be the same item

- ✓ Compatible with light-, self- and dual-cured materials (without the use of additional activators)
- ✓ Only 1-2 coats of Parts A & B are required, ensuring a more uniform and ideal film thickness while saving application time
- ✓ Ethanol-based formula allows less technique sensitivity to moisture resulting in more consistent and predictable results
- ✓ ALL-BOND 3 Resin is HEMA-free for more hydrophobic bonding to prevent degradation over time


The gold standard in total-etch bonding just got better with BISCO's ALL-BOND 3 Universal Dental Adhesive System.

ALL-BOND 3 utilizes patented technology that takes total-etch bonding to a whole new level of performance, assuring dentists of a simplified bonding process, as well as stronger, longer lasting bonds than ever before.



▶ ALL-BOND 3® is a truly universal adhesive. As shown to the above, ALL-BOND 3® has exceptional bond performance to all common dental substrates.

\* Licensed under 1 or more of the following U.S. Patents: 5,789,610, 5,270,351 & 5,401,783. See also U.S. Patent 5,348,988. Patent Pending.





# TheraBase™

Self-Adhesive Calcium Releasing Base/Liner

## Every great restoration starts with a great foundation.

TheraBase is a dual-cure, calcium and fluoride releasing, self-adhesive base/liner. Utilizing the THERA technology, TheraBase chemically bonds to tooth structure, and releases and recharges calcium and fluoride ions.<sup>1</sup> TheraBase's calcium release generates an alkaline pH which promotes pulp vitality.<sup>2</sup> It is a dual-cured material that will polymerize even in deep restorations where light cannot reach.

TheraBase is stronger and more durable than other base materials, glass ionomers and resin-modified glass ionomers.\* Additionally, it is radiopaque allowing for easy identification on radiographs, providing a quick and effective diagnosis.



**\$169.00**

TheraBase (H35001P)

**BUY 1, Get 1 FREE**



### Releases Calcium and Fluoride

Continuous release of Calcium and Fluoride ions<sup>1</sup>



### Self-adhesive

No bonding agents required - Save time and money



### High flexural strength

Stronger and more fracture resistant



### High compressive strength

Absorbs shock and stress from occlusal forces without fracturing



### Contains MDP

Contains the adhesion promoting monomer MDP, ensuring reliable and optimal bond to dentin<sup>3</sup>



### Easy to use

Auto-mix, dual-syringe provides a consistent mix for immediate delivery with zero to minimal waste of material.



### Dual-Cured

Material will fully cure even in deep restorations where light cannot reach.



### Alkaline pH

Generates an Alkaline pH (pH=11\*) in minutes, which promotes pulp vitality<sup>2</sup>



### High degree of conversion

Ensures enhanced physical properties



### Radiopaque


TheraBase is radiopaque allowing for identification on radiographs and effective diagnosis.



### THERA HYDROPHILIC MATRIX



BISCO has developed a hydrophilic matrix that allows for ion exchange. Previous or traditional resin matrices have been hydrophobic, but BISCO's matrix allows for ion exchange as water goes into the matrix, reacts, and calcium hydroxide ions and fluoride ions are released. At the calcium-rich interface between TheraBase and the dentin, an alkaline pH environment is created.



In deep restorations and when pulp exposure occurs, BISCO recommends the use of TheraCal LC for direct and indirect pulp capping and as a liner. Both TheraBase and TheraCal LC can be used together in a sandwich technique case scenario for optimum calcium and fluoride release benefits and ultimate pulp protection.

\* Data on file.

1. Gleave CM, Chen L, Suh BI. Calcium & fluoride recharge of resin cements. Dent Mater. 2016 (325):e26

2. T. Okabe, M. Sakamoto, H. Takeuchi, K. Matsushima. Effects of pH on Mineralization Ability of Human Dental Pulp Cells. Journal of Endodontics. Volume 32, Number 3, March 2006.

3. Hydrolytic stability of self-etch adhesives bonded to dentin, S Inoue 1, K Koshiro, Y Yoshida, J De Munck, K Nagakane, K Suzuki, H Sano, B Van Meerbeek, Journal of Dental Dentistry, December 2005



# TheraCal LC®

## 5 REASONS to Try TheraCal LC® Today!

- 1 **CALCIUM RELEASE¹**  
Stimulates hydroxyapatite and secondary dentin bridge formation¹,²
- 2 **ALKALINE pH**  
Promotes healing and apatite formation²,⁴
- 3 **PROMOTES HEALING¹**  
Liner insulates the pulp³, providing virtually no post-operative sensitivity
- 4 **INSULATION**  
Forms a protective barrier that insulates the pulp⁵,⁶
- 5 **MOISTURE TOLERANT¹ & RADIOPAQUE**  
Can be placed under restorative materials and cements



**BUY A Kit, Get 1 Syringe FREE**

Value  
**\$55**

**\$177.00** | TheraCal LC - Kit (H-33014P)  
4 Syringes TheraCal LC (1g each)

**\$55.00** | TheraCal LC (H-3301P)  
1 Syringe TheraCal LC (1g)

## TheraCal LC® CLINICAL TIPS

**APPLY TO  
VISIBLY MOIST  
DENTIN**



**APPLY IN  
INCREMENTS OF  
LESS THAN 1mm**



**LIGHT CURE  
FOR 20  
SECONDS**



1. Gandolfi MG, Siboni F, Prati C. Chemical-physical properties of TheraCal, a novel light-curable MTA-like material for pulp capping. International Endodontic Journal. 2012 Jun;45(6):571-9.  
2. ADA definitions for direct and indirect pulp capping at <http://www.ada.org/en/publications/cdt/glossary-of-dentalclinical-and-administrative-ter>  
3. Apatite-forming Ability of TheraCal Pulp-Capping Material, M.G. GANDOLFI, F. SIBONI, P. TADDEI, E. MODENA, and C. PRATI J Dent Res 90 (Spec Iss A): abstract number 2520, 2011 ([www.dentalresearch.org](http://www.dentalresearch.org))

4. Okabe T, Sakamoto M, Takeuchi H, Matsushima K (2006) Effects of pH on mineralization ability of human dental pulp cells. Journal of Endodontics 32, 198-201.  
5. Sangwan P, Sangwan A, Duhan J, Rohilla A. Tertiary dentinogenesis with calcium hydroxide: a review of proposed mechanisms. Int Endod J. 2013; 46(1):3-19  
6. Selcuk SAVAS, Murat S. BOTSALI, Ebru KUCUKYILMAZ, Tugrul SARI. Evaluation of temperature changes in the pulp chamber during polymerization of light-cured pulp-capping materials by using a VALO LED light curing unit at different curing distances. Dent Mater J. 2014;33(6):764-9.

# TheraCal PT™

Dual-Cured Resin Modified Calcium Silicate *Pulpotomy Treatment*

## New to the “Thera” family

TheraCal PT is a biocompatible, dual-cured, resin-modified calcium silicate that is used to treat exposed dentin. TheraCal PT maintains tooth vitality by performing as a barrier and protectant of the dental pulpal complex.


TheraCal PT's precise and direct placement allows its use in deep cavity preparations. The dual-cured set permits immediate placement of the restorative material.

### TheraCal PT is primarily indicated for:

- Pulpotomies: TheraCal PT may be placed after the partial or full removal of the coronal pulp.

### And can also be used for:

- Pulp exposures (direct pulp capping): TheraCal PT may be placed directly on pulp exposures after hemostasis is obtained.
- Protective liner (indirect pulp capping) or base under a variety of substrates.



20% OFF

**\$100.00** ~~**\$125.00**~~

**TheraCal PT**

Single Syringe Pack  
(H-34110P)



TheraCal PT syringe should remain in the pouch until ready to use.

**Once pouch is opened, use TheraCal PT within 90 days.**



Contents of one (1) syringe are good for **approximately 12-15 pulpotomy applications.**



The chemical formulation of TheraCal PT consists of synthetic Portland Cement calcium silicate particles in a hydrophilic matrix which facilitates calcium release.



TheraCal PT's physical properties resist breakdown and degradation leading to a durable seal.



Alkaline pH  
pH = 11.5 @ 7 days.




Moisture tolerant  
Low water solubility.



Radiopacity  
TheraCal PT is radiopaque (2.45mm) allowing for easy identification and differentiation from recurrent decay and other restorative materials.

## Cumulative Calcium Release





# Recommended Clinical Procedure for Zirconia Bonding



**TheraCem™**  
Self-Adhesive Resin Cement



## 5 REASONS to Try TheraCem™ Today!

- ✓ Continuous calcium & fluoride release
- ✓ Transitions from acidic to alkaline pH in minutes
- ✓ Contains MDP, allowing for strong bonds to zirconia, metal and alumina substrates without the use of a primer
- ✓ Quick and easy clean-up
- ✓ Highly radiopaque for quick and effective diagnosis



**BUY 4, Get 1 FREE**

**\$169.00**

1 Dual-Syringe (8g), Natural Shade (D-46311P)  
15 Mixing Tips and Intraoral Tips, Instructions

# UNI-ETCH® 32% w/BAC / ETCH® 37 w/BAC

**The only etchants on the market with BAC and thickened with polymers**

The addition of **Benzalkonium Chloride** (BAC) provides enhanced anti-microbial action and bond durability by inhibiting Matrix Metalloproteinases (MMPs) activated by acid-etching dentin. Thickened with polymers, Bisco's etchants rinse away easily and leave tooth surfaces free of residue.

The degradation of acid-etched dentin was shown to be caused by the presence of endogenous Matrix Metalloproteinases (MMPs) in dentin. BAC is effective at inhibiting both soluble recombinant MMPs and matrixbound dentine MMPs in the absence of resins.



**\$78.00**

Bulk Syringe (30ml)  
E56621P

**\$85.00**

4 Syringes (5g ea.)  
E5502EBM

**BUY 4, Get 1 FREE\***

\*Free item must be the same item

**Without BAC**

Etching without BAC activates MMP's that begin to **degrade** your adhesive layer.

**IMPROVE**

Your Bond's Durability with BAC

**With BAC**

Bisco Etches with BAC inhibit the MMPs, **creating a more durable bond.**



## SELECT HV ETCH

**Ideal pinpoint accurate viscosity for self-etch techniques**

SELECT HV ETCH delivers a smooth, non-stringy material which can be used in everyday restorative procedures.

- Thixotropic
- Coloured for easy visualization and contrast
- Washes off easily without leaving residue
- **35%** High Viscosity Phosphoric Acid Dentin/ Enamel Etchant

**\$78.00**

Bulk Syringe (30ml)  
(E59060P)

**\$85.00**

4 Syringes (5g ea.)  
(E59010P)

**BUY 4, Get 1 FREE\***

\*Free item must be the same item

# CAVITY CLEANSER™

**2% Chlorhexidine Digluconate**

**Easily applied with brush or cotton pellet.**

BISCO'S CAVITY CLEANSER is a 2% solution of Chlorhexidine Digluconate (CHX) intended for cleansing and moistening/rewetting cavity preparations. It is recommended for use upon completion of tooth preparation or etching prior to sealing the dentinal tubules. Cleansing of cavity preparations to remove debris and bacteria will decrease post-operative sensitivity.

**\$84.15** ~~**\$99.00**~~  
Cavity Cleanser (135ml)  
(B2816A)

**15% OFF**



# SUPER SEAL®

**Dental Desensitizer**

Super Seal® is a unique formula that (in one step) removes the smear layer, seals the tubules and desensitizes. No pumice cleaning or rinsing is required. Because of its acidic nature, Super Seal® demineralizes the smear layer (both organic and mineral debris) and the peritubular dentin (outermost ring of hard mineralized dentin of each tubule complex).

**\$58.50** ~~**\$65.00**~~  
Super Seal (4ml)  
(100294)

**10% OFF**



## DOCTOR USES:

- Periodontal Surgery
- Any Exposed Dentin
- Re-Wetting Agent
- During Temporization
- For Restorations
  - Amalgams
  - Composites
  - Crowns and Veneers
  - Inlays and Onlays



## HYGIENIST & ASSISTANT USES:

- Thermal Sensitivity
- Before and After Prophylaxis
- Before and After Whitening Treatments
- Exposed Root Surfaces
- Perio Scaling and Root Planing



## GENERAL FACTS

- **NO GLUTARALDEHYDE**
- **Soft Tissue Friendly**
- Water-Based Potassium Oxalate
- Chelation Starts Within 5 Seconds
- Removes Smear Layer
- Re-Wetting Agent
- Bubble-Smear Free and Accurate Impression
- Do Not Refrigerate
- No Light Curing
- Simple One-Step Process



# Start treatment with a smile

support of the cervicothoracic junction

support of the cervical spine

stable position of the head

**Handmade, extensively tested and developed into the perfect shape.** Happynecks® provide optimal support from the cervical spine to the cervicothoracic transition for maximum patient comfort and relaxation during dental treatment



HAPPYNECK® ANTHRACITE  
(21-HN004)



HAPPYNECK® SKY  
(21-HN003)



HAPPYNECK® AQUA  
(21-HN006)


## Quality & Durability

- Developed by a Physiotherapist
- Handmade in the Netherlands
- Durable exterior made from skai® (faux leather)
- High-quality Pantera foam

**\$189.00** ~~**\$209.00**~~ **Introductory Offer**

## Comfort & Ergonomics

The Happyneck® is designed in such a way that the shape of the pillow follows the curvature of the cervical spine. This even support of the cervical spine stimulates the extension of the neck, making it easier for the head to lean backwards. The end positions flexion/extension, lateral flexion and rotation give the practitioner an optimal view of the working area while maintaining an optimal ergonomic working posture.



# On the cusp of perfection.

## LMArte™ Cusp Misura EVALUATE EXPLICITLY



Unique instrument for estimating cusp thickness and height for evaluating the need for cusp preparation when doing adhesive restorations. It's ideal for evaluation of the strength of the natural tooth structure and decide which restorative technique to be used. If the forks shaped instrument tip penetrates up to its marks, the tooth structure is too weak to withstand occlusal forces with a direct restoration, and other clinical options should be considered.



NEW  
PRODUCT



**\$58.50** ~~**\$65.00**~~

LM468-469ES

10% OFF

## Your **New** Favourite Sharpen Free Instrument



NEW  
PRODUCT

SICKLE LM204S SD

- For the removal of supragingival calculus.
- Sharply angulated shank and very fine working end.
- With LM Sharp Diamond™ sharpen free coating

**\$75.60** ~~**\$84.00**~~

10% OFF



# LM Dental

## EXTRACT WITH CONFIDENCE



**\$93.80** ~~**\$134.00**~~ / ea  
Extraction Instrument

**30% OFF**

### LM SlimLift™

Luxating instrument designed for extremely narrow spaces.

#### SlimLift

- Fit in extremely narrow spaces
- Ideal for implant preparations
- Available blade sizes - Straight: 3mm, 4mm, 5mm  
- Curved: 3mm, 4mm, 5mm



### LM TwistOut™

Elevators for forceful extraction.

#### TwistOut

- Sturdy extraction instrument - you can apply more force, torque and twist
- Available blade sizes - Straight: 3mm, 4mm, 5mm  
- Curved: 3mm, 4mm, 5mm



### LM LiftOut™

Luxating instrument for atraumatic extraction.

#### LiftOut

- Finer blades; should not be used like elevators
- Available blade sizes - Straight: 2.5mm, 3mm, 4mm, 5mm  
- Curved: 3mm, 4mm, 5mm



All LM extraction instruments are supplied in a protective cassette.



LM extraction instruments uniquely combine ergonomics, Scandinavian design and functionality. They feature comfortable, non-slip ErgoTouch™ handles and are well-balanced and lightweight.



# LM Titanium Implant Instruments


## Maintain with confidence




**\$104.40** ~~**\$116.00**~~ / ea  
Implant instrument

**10% OFF**

**Replaceable tips offer many benefits:**



- Always sharp instruments
- Saves time from sharpening
- Minimizes waste
- Enables personalized tip combinations



**Lock Grip** - sealed and secure construction  
Use LM-ErgoMix in the normal way, clean and sterilize it like any other instrument. The sealed construction guarantees high level of hygiene.

To retip turn the lock grip.  
No tools needed!



Implant Universal Curette, Posterior/Anterior  
(LM283-284MTiEM)



Implant Mini Gracey 1-2, Anterior  
(LM201-202MTiEM)



Implant Mini Gracey 11-12, Posterior - mesial  
(LM211-212MTiEM)



Implant Mini Gracey 13-14, Posterior - distal  
(LM213-214MTiEM)

- Softer, Titanium Alloy Tips
  - Safe and compatible material
- Clinically Relevant Design
  - Mini-sized, slim tips fit under the tissue
- User-Friendly Replaceable Tips
  - Economical & ecological
  - No tools needed to replace tips
  - Tips do not come loose (not cone-sockets)
- Great Ergonomics
  - Large diameter (12mm) handle
  - Grippy medical-grade silicone surface

# LUCIDA

COMPOSITE  
GLOSS  
SYSTEM



## IT'S ALL ABOUT THE FINISH.

A simple and elegant chairside polishing system for composites, that provides a natural, high lustre, "enamel-like" finish.

Lucida™ powered by Style Italiano™ and DiaShine®, Lucida™ is your complete solution for the chairside polishing of composites, provisionals, porcelain veneers & denture-based acrylics. The Lucida system combines DiaShine's proprietary, water soluble, chairside polishing compound and disposable Lucida Felt Stars to give you a system that is both versatile and effective.

Developed by Style Italiano, the Lucida star-shaped felts, have rounded points to allow you to polish facially, interproximally and occlusally. The felts offer a variable thickness, depending on the pressure and speed applied, are disposable, and come with autoclavable mandrels.

**\$139.30** ~~**\$199.00**~~

**30% OFF**


### INTRO KIT (LUC-KIT)

- 1 Lucida Paste, 5 Lucida Mandrels, 70 Lucida Stars

Lucida Star Felts (100 pack)	LUC-FELT100	<b>\$35.00</b>
Lucida Diamond Paste (2g syringe)	LUC-PASTE	<b>\$89.00</b>
Lucida Mandrel (5 pack)	LUC-MAN5	<b>\$35.00</b>


DIASHINE.

STYLE ITALIANO  
SHADE GUIDES



Clinical cases by Dr. Marcin Krupinski – Style Italiano

## LUCIDA PASTE



**Lucida Paste** provides a luxurious gloss with its patented submicron hybrid water-soluble polishing compound. Its unique formula allows the paste to stay and spread perfectly on the tooth surface while polishing, yet goes away immediately when water is applied.

## LUCIDA STAR




The unique **Lucida Star** felt polishing wheel solves many issues of conventional felt polishing wheels on several fronts. Its unique thickness and shape allows access to hard-to-reach areas, provide intermittent touch and variable thickness, depending on the force and speed applied.

# K-LITE

by Katherine Losada



## MAKE THE INVISIBLE VISIBLE



The aim of a minimally invasive dentistry is to preserve as much of the sound tooth tissue as possible. Modern dental composites and adhesives contain fluorophores, giving them a fluorescent property. Under certain wavelengths of illumination (**UV-A 395 to 400nm**), these materials can be distinguished from tooth tissue. The K-Lite can help reduce the risk of iatrogenic tooth damage when these materials need to be removed.

The K-Lite is a ground-breaking dental device that aids in the:

- Identification of supra and subgingival calculus, dentine-enamel junction (DEJ).
- Detection of caries, demineralization areas of enamel, root canals by structure differentiation.
- Removal of excess or unwanted dental materials such as composites, adhesives and cements.
- Visualization of crack lines and fractures on natural teeth and crowns.



The K-Lite is a dual portable flashlight guided by a glass fibre optic, which consists of:

- a White LED (6000 °K)
- a UV-A LED (violet light)

**\$625**

**Buy Kit, Get Free Box Disposable sleeves**


- K-Lite: Refill, Disposable Protective Sleeves (200pcs.) **\$55**
- K-Lite, Light Guide, Glass Fiber **\$138**



Presence of pigments due to bacteria degradation in a cervical lesion.



Control after removal of infected area.



K-Lite shows the restoration margins and integrity of the composite filling.



Dual use of K-Lite : UV light shows fluorescent filling and white light shows cracks of tooth.






## NiTInol Sectional Matrix System

By Garrison Dental. A Bisco Canada Exclusive



The same **PEEK®** (polyether ether ketone) superpolymer reinforcement technology used to reduce ring stretch as with the award-winning **Composi-Tight®** systems.


Nitinol drawn-wire nickel titanium is superior to stainless steel or laser-cut, solid NiTi.



Precisely engineered tip geometry allows NiTin™ to grip where others pop-off, all while reducing flash.

**Your contacts are guaranteed - because NiTin retains an industry-best 93.7% of its original tooth-separating force after 300 cycles!\***

\*Here's how the torture test cycle works: The NiTin ring is placed on a typodont for 5 minutes, then removed and exposed to a curing light. It gets a full ultrasonic cleaning bath and finished up with a cycle through a steam autoclave. The separating force is measured at regular intervals.



### NiTIn Intro Kit



1 Ring Placement Forceps  
1 NiTin Standard Ring - White  
1 NiTin Long Ring - Black  
40 Assorted Matrices  
40 Assorted Wedges

20% Off

**\$358.40** ~~**\$448.00**~~  
10-NT-KMN-40

### Rings



NiTIn™ Standard  
(NT500)

NiTIn™ Long  
(NT400)

**\$236.00**  
2/Pack

**BUY 4, Get 1 FREE\***

\*Free item must be the same item

### Full Contour Bands



NTM100  
Premolar, 4.0 mm

NTM150  
Premolar, 4.1 mm

NTM175  
Molar, 5.0 mm



NTM200  
Tall Molar, 5.4 mm



NTM300  
Tall Molar, 6.0 mm

**\$335.00**  
Kit

**\$75.00**  
50/Pack

**BUY 4, Get 1 FREE\***


\*Free item must be the same item



# Fusion™ Anterior Matrix System



The Fusion Anterior Matrix System helps you easily express your artistic nature restoring beautiful smiles for your patients.





Fusion Anterior Matrix bands help you recreate ideal interproximal anatomy quickly and confidently:

- Only 0.0015" – significantly thinner than strips to facilitate better interproximal contact.
- Firm stainless steel resists deformation both during placement and while compressing composite into the preparation.
- Easy to slide into the sulcus for deeper restorations.
- Ideal anatomic curvature in the gingival-incisal direction AND facial-lingual.
- Far superior to strips when restoring large or deep carious lesions.

The clever features of the Fusion Anterior Wedge make it a 'must have' for all anterior restorations:

- Radical curvature ensures a firm seal at the cervical margin from facial to lingual.
- Helps maintain ideal anatomy by hugging the band to the tooth.
- Unique 'T' shape sits deeper in the interproximal space to minimize the dreaded 'black triangle'.
- Tapered separator strut allows you to adjust the amount of pressure exerted to control the interproximal contact.
- The cervical interface is half above and half below the sulcus for maximum adaptation with deeper restorations.



## Fusion Anterior Kit



50 short bands  
50 tall bands  
25 extra small wedges  
25 small wedges  
25 medium wedges  
25 large wedges

Introductory Price - 10% Off

**\$194.40** ~~\$216.00~~  
ANK01

## Fusion Anterior Bands



Tall

Short

**\$93.00**  
50/pk

**\$163.00**  
100/pk

**BUY 4, Get 1 FREE\***

\*Free item must be the same item

## Fusion Anterior Wedges



**\$123.00**  
Kit (200)

**\$34.00**  
50/pk

**\$62.00**  
100/pk

**BUY 4, Get 1 FREE\***

\*Free item must be the same item



# ONLINE ONLY PROMOTIONS

Please visit [www.biscocanada.com](http://www.biscocanada.com) to take advantage of more products on sale!

Receive **FREE SHIPPING** on all online orders.



## Quick-Stat™ FS

15.5% Ferric Sulfate Hemostatic Gel



20% OFF

## Stat-Flo™ Tip



20% OFF

### ONLINE ONLY PROMOTIONS

## Best Etch

37% Phosphoric Acid Etch

- Designed to etch and prepare dentin and enamel.



20% OFF

## SmearOFF

The only TRUE two-step, 2-in-1 irrigant



20% OFF

### ONLINE ONLY PROMOTIONS

## LM SharpDiamond™

100% Sharpen-Free



10% OFF

## LM DuraGradeMax™

Stainless Steel Instruments



10% OFF



# ONLINE ONLY PROMOTIONS

Please visit [www.biscocanada.com](http://www.biscocanada.com) to take advantage of more products on sale!

Receive **FREE SHIPPING** on all online orders.



**FitStrips** are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts.

**\$266.00**

**Universal Kit (FPSK04)**

- 2 handles
- 1 IPR gauge
- 1 single-sided of each grit (4 pcs)
- 1 double-sided of each grit (6 pcs)

**\$266.00**

**IPR Kit (FPSK05)**

- 2 handles
- 1 IPR gauge
- 2 double-sided of each grit (10 pcs)

**\$237.00**

**Subgingival Finishing Kit (FPSK07)**

- 2 handles
- 1 serrated saw
- 3 subg super fine, 3 subg fine, 3 subg medium

**NEW**



## Sanctuary™



**\$23.00**

6"x6" Sanctuary  
Latex Dental Dams (36 Sheets)

**BUY 4 GET 1 FREE**



**\$32.00**

6"x6" Sanctuary  
Non-Latex Dental Dams (15 Sheets)

**BUY 4 GET 1 FREE**



**\$34.00**

6"x6" Sanctuary YIN YANG  
Non-Latex Dental Dams (15 Sheets)

**BUY 4 GET 1 FREE**



**143.10 \$159.00**

Dental Dam Punch (157mm)

**10% OFF**

## ONLINE ONLY PROMOTIONS



Clean, deodorize & disinfect without toxic

**PRO Kit (FON-BUNDLE80-001)**

- Electrolyzer
- Power Cord
- 85 Activator Capsules
- 3 Spray Bottles
- Quick Start Guide

**\$151.30 \$178.00**

**15% OFF**



**85**

Activator Capsules

ON SALE



# AWARD-WINNING PRODUCTS

with unique benefits



TheraCal LC (2021 Product Award), Pulpal Protectant  
TheraCal PT (2021 Product Award), Pediatric Product  
TheraBase NEW! Self-Adhesive Base/Liner

TheraCem (2021 Product Award), Cement: Self-Adhesive  
All-Bond Universal (Editor's Choice)

Visit our **Virtual Booth** for product presentations & exclusive offers.

2021 **VIRTUAL**  
SASKATCHEWAN ORAL  
HEALTH CONFERENCE

SEPTEMBER 25-26, 2021

**EDDS** est. 1913  
EDMONTON & DISTRICT DENTAL SOCIETY

**Northwest Dental Exposition 2021**

**VIRTUAL EVENT** on OCTOBER 2, 2021

