



## BISCO Inc.

Backed by four decades of scientific research and innovation, BISCO products continue to change the game of adhesive dentistry and beyond.

Curion is proud to be the only authorized and exclusive distributor of BISCO products in Canada



To view all Bisco products, scan QR code.



### All-Bond Universal™

Universal, light-cured dental adhesive



\$237 Kit  
\$223 Bottle (6ml)  
\$264 50 Unit-Doses  
\$498 100 Unit-Doses

**BUY 3 + 1\***

### One-Step® One Step Plus®

Universal dental adhesives



\$212 One-Step  
\$212 One-Step Plus

**BUY 4 + 1**

### Select HV® Etch

35% high viscosity phosphoric acid



\$101 4 Syringes  
\$95 Bulk Syringe (30ml)

**BUY 3 + 1\***

### 32% Uni-Etch w/BAC Etch 37% w/BAC

Phosphoric acid etchant with BAC



\$101 4 Syringes  
\$95 Bulk Syringe (30ml)

**BUY 3 + 1\***

### TheraCal PT®

Pulpotomy Treatment



\$123.20 Reg. \$154 1 Syringe

**SAVE 20%**

### TheraBase®

Calcium & Fluoride releasing self-adhesive base/liner



\$144.90 Reg. \$207 1 Syringe

**SAVE 30%**

### Duo-Link Universal

Resin cement



\$442 Kit  
\$180 Refill Syringe

**BUY 3 + 1\***

### TheraCem®

Self-Adhesive Resin Cement



\$155.25 Reg. \$207 1 Syringe

**SAVE 25%**

### Z-Prime Plus

Zirconia, Alumina and Metal Primer



\$143.65 Reg. \$169 1 Bottle

**SAVE 15%**

### ZirClean

Restorative cleaner



\$48 1 Syringe

**BUY 3 + 1**

### Bis-Silane™ Porcelain Primer



**SAVE 10%**

\$114.30 Reg. \$127 Bis-Silane (Part A&B)  
\$107.10 Reg. \$119 Porcelain Primer (1 Bottle)

### Core-Flo™ DC Core-Flo™ DC Lite



\$127 Refill Syringe  
\$381 Core-Flo DC Lite Kit

**BUY 3 + 1\***

\*Free item must be of equal or lesser value.

**Slow-Speed Handpieces**

Designed for maximum comfort and optimal balance.

Curion is proud to be the only authorized and exclusive distributor of Curion Essence products in Canada.



To view all Curion essence products, scan QR code.



**E-Type 1:1 Ball Bearing Auto-Latch Contra Angle**

Slow speed, One Piece ISO E-Type attachment with a 1:1 gear ratio.



**\$289/ea**  
Reg. \$320

**BUY 1, 2 or 3**

**\$1,280**  
Net. \$256/ea

**BUY 4 + 1**

**\$2,880**  
Net. \$240/ea

**BUY 9 + 3**

**NSK Type 20K Air Motor**

The NSK Type 20K Air Motor features four drive air holes ideally used for slow-speed handpieces.



**\$349/ea**  
Reg. \$400

**BUY 1, 2 or 3**

**\$1,600**  
Net. \$320/ea

**BUY 4 + 1**

**\$3,600**  
Net. \$300/ea

**BUY 9 + 3**

**E-Type 4:1 Screw-in Prophy Contra Angle**



Slow speed, one piece ISO E-Type attachment with a 4:1 gear ratio.



**1 YEAR**  
Warranty

**\$289/ea**  
Reg. \$320

**BUY 1, 2 or 3**

**\$1,280**  
Net. \$256/ea

**BUY 4 + 1**

**\$2,880**  
Net. \$240/ea

**BUY 9 + 3**

**E-Type 4:1 Snap-on Prophy Contra Angle**



**E-Type 4:1 BB Push Button Latch Contra Angle**



**HAPPYNECKS**

Superior quality, handmade in the Netherlands.

Curion is proud to be the only authorized and exclusive distributor of Happynecks in Canada.



To view all Happynecks products, scan QR code.



**Happynecks**

Optimal neck support pillow for a better patient experience.



\$160 Regular Size

**BUY 4 + 1**

**Happynecks - Kids**

The children's version of the Happynecks.



\$140 Kids Size

**BUY 4 + 1**

**Happynecks - Lumbar**

Optimal neck lumbar support pillow for a better patient experience.



\$175 Lumbar

**BUY 4 + 1**

\*Free item must be of equal or lesser value.



# ON SALE NOW!

2024 2nd Edition | Prices expire June. 30, 2024

EN: 1.800.667.8811

FR: 1.800.211.1200

curion.ca

## Smile Line

made in switzerland



### SMILE LINE

Something about Smile Line products conjure a feeling of uniqueness, exclusivity and a dreamy quality that is often not present in dental products.

Curion is proud to be the only authorized and exclusive distributor of Smile Line products in Canada.



To view all Smile Line products, scan QR code.



#### Smile Lite MDP-2

Taking smartphone dental photography to a new level.

NEW! Bluetooth Connected Handle



**\$859** Reg. \$949  
Smile Lite MDP-2 Kit

SAVE \$90

#### Optishade StyleItaliano™

Digital shade matching in the palm of your hands!




**\$2,569.50** Reg. \$2,855  
Optishade StyleItaliano Full Set

SAVE 10%

#### Compo-Vibes

A magical tool designed for easier, faster and better composite modeling!



**\$402.30** Reg. \$447 Master Set  
**\$358.20** Reg. \$398 Basic Set  
**\$76.50** Reg. \$85 2 Refill Tips  
**\$127.80** Reg. \$142 4 Refill Tips

SAVE 10%

#### K-Lite

Make the invisible, visible!



**\$591.30** Reg. \$657 K-Lite Kit  
**\$127.80** Reg. \$142 Refill Light Guide  
**\$52.20** Reg. \$58 Disposable Sleeves

SAVE 10%



## DMG America

A smile ahead together

Curion is proud to be an authorized distributor of DMG products in Canada.



To view all DMG products, scan QR code.



#### Icon - Caries Infiltration

The solution to treating caries while preserving healthy tooth structure.



BEST SELLER

SAVE 10%

**\$293.40** Reg. \$326 Mini Kit - 2 Patient Packs (Smooth Surface or Proximal)  
**\$924.30** Reg. \$1024 Cube - 7 Patient Packs (Smooth Surface or Proximal)  
**\$189.00** Icon Etch Refill

#### Luxatemp Ultra

The best just got better!  
Rated #1 in flexural strength



TOP PICK

**\$388** Refill (A1, A2, A3, B1, Bleach Light)

BUY 3 + 1

#### LuxaFlow Ultra

LuxaFlow Ultra for life-like repairs



**\$84** Refill (A1, A2, A3, B1)

BUY 3 + 1

#### TempoCem ID

This aesthetic, temporary cement is invisible when you want it to be, yet detectable when you need it to be.



TOP PICK

**\$133** Smartmix

BUY 3 + 1

\*Free item must be of equal or lesser value.

## TRI HAWK

50 years of designing, building, and providing burs with an unmatched combination of cutting speed, strength, and safety to bring our customers the best burs in the dental industry.

Curion is proud to be the only authorized and exclusive distributor of Tri Hawk burs in Canada.



To view all Tri Hawk products, scan QR code.



### TALON

Unmatched combination of cutting speed, strength and safety.

\$120 25 Roll  
\$439 100 Roll

BUY 4 + 1\*

### Metal & Crown Cutters



**Talon 8**      **Talon 10**      **Talon 12**

	Talon 8	Talon 10	Talon 12
Head Diameter (mm)	0.8	1.0	1.2
Head Length (mm)	4.2	4.2	4.2
Overall Length (mm)	19.5	19.5	19.5

### Crown Cutters Surgical

\$135 25 Roll  
\$490 100 Roll

BUY 4 + 1\*



**Talon 10SU**      **Talon 12SU**      **Talon 14SU**

	Talon 10SU	Talon 12SU	Talon 14SU
Head Diameter (mm)	1.0	1.2	1.4
Head Length (mm)	4.2	4.2	4.2
Overall Length (mm)	25.5	25.0	25.0

### AMAZILIA

Solid carbide one piece construction to reduce breakage.

Durable blade design for long bur life.

\$197 25 Roll  
\$720 100 Roll

BUY 4 + 1\*



**FG151**      **FG151 L**      **FG152**

	FG151	FG151 L	FG152
Head Diameter (mm)	1.6	1.6	1.6
Head Length (mm)	10.0	10.0	10.0
Overall Length (mm)	25.0	28.0	25.0

### Finishing Carbide

\$135 25 Roll  
\$490 100 Roll

BUY 4 + 1\*



**FG7901**      **FG7902**      **FG7406**

	FG7901	FG7902	FG7406
Head Diameter (mm)	0.9	1.0	1.6
Head Length (mm)	3.5	3.5	3.5
Overall Length (mm)	19.5	19.5	19.5

### Standard Carbides

\$65 25 Roll  
\$228 100 Roll

### Surgical Carbides

\$135 25 Roll  
\$490 100 Roll

BUY 4 + 1\*



Shapes	Amalgam Cutter	Crosscut Straight Fissure	Straight Fissure	Crosscut Tapered Fissure	Tapered Fissure	Pear	Inverted Cone	Round
<b>Standard Carbides</b>	FG245	FG556, FG557, FG558	FG56, FG57, FG58	FG700, FG701	FG170, FG171	FG330, FG331, FG332	FG33 1/2, FG35, FG37	FG 1/2, FG2, FG4, FG6
<b>Surgical Carbides</b>		FG557SU, FG558SU	FG57SU, FG58SU	FG701SU, FG702SU	FG170SU, FG171SU			FG2SU, FG4SU, FG6SU

\*Free item must be of equal or lesser value.



# ON SALE NOW!

2024 2nd Edition | Prices expire June. 30, 2024

EN: 1.800.667.8811

FR: 1.800.211.1200

curion.ca



## EVE Rotary

Because the tool defines the result.

Curion is proud to be an authorized distributor of EVE products in Canada.




To view all EVE products, scan QR code.



**MOST POPULAR**

### DIACOMP® PLUS RA

Polishing System for Composite, Hybrid Ceramic, 3D Printed Hybrid Material



**2** PROCESSING STEPS  
ÉTAPES DE TRAITEMENT

DIAMOND IMPREGNATED POLISHING SYSTEM

#### Starter Kit

- \$95 Diacomp Plus - Twist 14mm
- \$85 Diacomp Plus - Mix
- \$124 Diacomp Plus - Occluflex

#### Refill

- \$102 - \$155 Medium | Pre-polishing
- Fine | High-shine Polishing

**BUY 4 + 1\***


**MOST POPULAR**

### DIAPOL® RA

Polishing System for Feldspathic Ceramic, Lithium Disilicate, Zirconia

Now available in HP

Recommended by GC Dental for polishing of GC LISI Blocks



**3** PROCESSING STEPS  
ÉTAPES DE TRAITEMENT

DIAMOND IMPREGNATED POLISHING SYSTEM

#### Starter Kit

- \$95 Diapol Twist
- \$112 Diapol Occluflex
- \$112 Diapol Set RA 307

#### Refill


- \$102 - \$155 Coarse | Smoothing
- Medium | Pre-polishing
- Fine | High-shine Polishing

**BUY 4 + 1\***

### DIACERA RA

Polishing System for Zirconia, Especially suitable for full contour zirconia restorations

Now available in HP



**2** PROCESSING STEPS  
ÉTAPES DE TRAITEMENT

DIAMOND IMPREGNATED POLISHING SYSTEM

#### Starter Kit

- \$95 Twist Mix 14mm

#### Refill

- \$102 - \$155 Medium | Pre-polishing
- Fine | High-shine Polishing

**BUY 4 + 1\***

### FLEXI-D RA

Polishing Discs for Composite, Amalgam, and Gold.

FLEXI-D 1/2" (10mm)



FLEXI-D 5/8" (14mm)



Coarse Medium Fine Extra Fine Mandrel

**4** PROCESSING STEPS  
ÉTAPES DE TRAITEMENT

#### Starter Kit

- \$58 Mix Set - 1/2" (10mm)
- \$58 Mix Set - 5/8" (14mm)

#### Refill


- \$52 1/2" (10mm)
- \$52 5/8" (14mm)
- \$15 3x Mandrel

**BUY 4 + 1\***

### FLEXISTRIPS

Composite, Amalgam, Gold

**NEW!**



- \$24 Set FS-50
- \$24 FS-60/40
- \$24 FS-20/7


- Flexible polishing strips for interproximal application
- For polishing composite and metal
- With gap between abrasives for easier application
- Coated on one side with aluminium oxide

**BUY 4 + 1**

### DIABRUSH RA

Composite, Feldspathic Porcelain

**NEW!**



- \$85 WTDI
- \$85 WZDI
- \$85 WSDI

- For high-shine polishing on composite and feldspar ceramic
- For polishing occlusal fillings
- Fibre with diamond grit
- Application without polishing paste
- Heat resistant and autoclavable polyamide fibre

**BUY 4 + 1**



**RE-GEN™ Bioactive EndoSealer**

**NEW!**

Bioceramic Root Canal Sealer

- **Contains Bioglass:** to promote bio-regenerative healing while attracting and exchanging ions that promote remineralization of hydroxyapatite
- **Ready to use:** injectable, single component paste
- **Biocompatible and Osteogenic:** bioactive root repair material that is effective on non-vital cases
- **Anti-bacterial:** creates a high alkaline environment to prevent bacteria formation
- **Dimensionally stable:** low solubility, high wash out resistance
- **Quick set time:** under 1 hour



**\$133** RE-GEN Bioactive EndoSealer Refill (2g)  
(20) MST Visco-Tips  
**\$51** MST Visco-Tips (20)  
**\$108** MST Visco-Tips (50)

**BUY 2 + 1\***

**Voyager Irrigating Tips**

**NEW!**

Dual side ports for safety and efficacy



27GA

29GA

30GA

31GA



**Dual Side Port**  
**\$92** 27GA, 29GA, 30GA, or 31GA  
Size (17mm or 27mm)



**Blunt End**  
**\$77** 29GA  
Size (17mm or 27mm)

**BUY 2 + 1**

**Endo Ultra**

Cordless Ultrasonic Activator



**BONUS OFFER**

**\$1048** Endo Ultra Kit  
**\$231** Activator Tips

**Buy 1 Kit, Get a Pack of  
Activator Tips for Free**

**CHX-Plus™**

Enhanced 2% Chlorhexidine Gluconate Solution



**\$104** 480ML Bottle

**BUY 3 + 1**

**SmearOFF™ 2-in-1**

The only TRUE two-step 2-in-1 irrigant



**\$152** 480ML Bottle

**BUY 3 + 1**

**Chlor-XTRA™**

6% sodium hypochlorite solution



**\$59** 480ML Bottle

**BUY 3 + 1**



# ON SALE NOW!

2024 2nd Edition | Prices expire June. 30, 2024

EN: 1.800.667.8811

FR: 1.800.211.1200

curion.ca

## Garrison

Dental Solutions

## GARRISON Dental

Humble beginnings. Future focused.

Curion is proud to be an authorized distributor of Garrison products in Canada.




To view all Garrison products, scan QR code.



### Strata-G Matrix System

270° Seal for ultimate flash control.



**\$827.10** Reg. \$919 Complete Kit  
**\$790.20** Reg. \$878 Basic Intro Kit

**SAVE 10%**

### FitStrip

Interproximal Finishing & Contouring

FitStrip make your finishing and contouring procedure easier!



**\$264** Starter Kit, Single-Sided Kit, or Subgingival Finishing Kit  
**\$297** Universal Kit or IPR Kit  
**\$264** Refill (10 Fitstrips + 2 Handles)  
**\$124** Refill (4 Fitstrips)  
**\$62** IPR Gauge

**BUY 4 + 1\***


### Strata-G Rings



**\$545** Kit 3/pk (1 ring ea)  
**\$373** Short Ring 2/pk  
**\$373** Tall Ring 2/pk  
**\$373** Wide Ring 2/pk

**BUY 4 + 1\***

### Strata-G Matrix Bands



**\$665** 400 Kit (100 ea)  
**\$352** 200 Kit (50 ea)  
**\$194** 100 Refill  
**\$113** 50 Refill

**BUY 4 + 1\***

### Strata-G Wedges



**\$298** 400 Kit (100 ea)  
**\$158** 200 Kit (50 ea)  
**\$87** 100 Refill  
**\$47** 50 Refill

**BUY 4 + 1\***

## Ascentcare

Dental Products

## ASCENTCARE Dental

VacuLUX Isolation systems are designed for simplicity and crafted with premium materials.

Curion is proud to be the only authorized distributor of Ascentcare Dental products in Canada.



To view all Ascentcare Dental products, scan QR code.



### VacuLUX™ Air Chairside Isolation System



**BONUS OFFER** Buy 2 Chairside Kits, Get 4 Boxes of 4-Mouthpiece FREE

**\$635** Chairside Isolation System - VacuLUX  
**\$694** Chairside Isolation System - DRYSHIELD or ZYRIS/ISOLITE

### VacuLUX™ Air Portable Isolation System



**BONUS OFFER** Buy 2 Portable Kits, Get 4 Boxes of 4-Mouthpiece FREE

**\$575** Portable Isolation System - VacuLUX  
**\$633** Portable Isolation System - DRYSHIELD or ZYRIS/ISOLITE

### VacuLUX™ Air Adapter



**BONUS OFFER** Buy 2 Air Adapters, Get 2 Boxes of 4-Mouthpiece FREE

**\$329** VacuLUX Adapter  
**\$399** DRYSHIELD or ZYRIS/ISOLITE Adapter

\*Free item must be of equal or lesser value.



**ON SALE NOW!**

2024 2nd Edition | Prices expire June. 30, 2024

EN: 1.800.667.8811

FR: 1.800.211.1200

curion.ca

**LM** feel the difference

**LM DENTAL**  
Feel the difference

Curion is proud to be the only authorized and exclusive distributor of LM Dental instruments in Canada.



To view all LM products, scan QR code.



**SharpDiamond™**  
100% sharpen-free instruments

**\$75 - \$87**  
Reg. \$100 - \$116

**SAVE 25%**

**New Designs**

**NV1 Sickle** (Compares to Nevi 1)  
• Popular sickle and disc combination design for efficient anterior instrumentation.

**NV2 Sickle** (Compares to Nevi 2)  
• Popular posterior sickle scaler design with heavily angled, long shank and thin, long blades for easy access to interproximals.

**LM23 Slim Sickle**  
• Same clinical application as the beloved 2045, but easier to adapt – now even slimmer!

**Mini McCall 13S-14S**  
• Universal curette for removal of supra and subgingival calculus from all tooth surfaces.

**Most Popular Items**

**SharpJack Sickle** (Compares to Nevi 3)  
• Universal Scaler with long, thin, dual-edged blades for easy and efficient supragingival scaling.

**Mini Dual-Gracey (Syntette) Curette**  
• Unique combination of universal curette convenience with the placement accuracy of 2 different Graceys.

**DuraGradeMax™**  
Stainless steel instruments

**\$68.40 - \$76.50**  
Reg. \$76 - \$85

**SAVE 10%**

**Most Popular Items**

204S Sickle Scaler

Standard Dual-Gracey (Syntette) Curette

Furcator KS

Standard Gracey 1/2 Curette

Standard Gracey 13/14 Curette

Standard Gracey 15/16 Curette

**LM Dark Diamond™**  
Non-stick coating instruments

**\$79.20**  
Reg. \$88

**SAVE 10%**

Plastic Instrument - Mini Universal

Plastic Instrument - Anterior

Plastic Instrument - Posterior

Burnisher

Burnisher-Pluggler

Packer-Modeller IO

Modeller MH

Modeller MH Mini

**Find us here!**

**ODA/ASM**  
April 18-19, 2024 | Toronto  
Booth #1219

**JDIQ**  
May 27-28, 2024 | Montreal  
Booth #1509

**RDH Evolve**  
May 3-4, 2024 | Edmonton

**Live Webinar**

**Predictable Bonding of Indirect Restorations**

**BISCO**

**Instructor: Joseph Kim, DDS, JD**

**Date and time:**  
Wednesday, May 8, 2024  
5:00 PM - 6:30 PM PST / 8:00 PM - 9:30 PM EST

**1.5 Hours**  
**Earn 1.5 CE Credits**  
**Tuition: Free**

Scan to register