



AELITEFLO™ AELITEFLO™ LV

Call Today to Place Your Order

1-800-667-8811

BISCO DENTAL PRODUCTS CANADA



LIGHT-CURED FLOWABLE MICROHYBRID COMPOSITE

AELITEFLO is a light-cured composite with flow characteristics which make it ideal for Class V restorations, as well as Class III's and small Class IV's. Other uses of this material are as a liner for Class I and II restorations, pit and fissure sealants, the repair of marginal defects, small core build-ups, blocking out undercuts, porcelain repair, porcelain veneer cementation or as a flowable composite for bonding fibrous splints into place. Even though AELITEFLO is flowable, its properties prevent it from slumping for better control in hard-to-access areas. Using the tips provided, AELITEFLO is easily injected directly into the cavity preparation and usually can be shaped with the provided tips.

AELITEFLO LV is a light-cured, low viscosity microhybrid composite with flow characteristics which make it ideal for lining Class I and Class II restorations and suited as a pit and fissure sealant.

- **AELITEFLO Intro Kit: 1.5g Syringes:** A1, A2, A3.5, C2, and D3, disposable syringe tips, instructions / MSDS
- **Refill Syringe (1.5g):** A1, A2, A2-Opaq, A3, A3.5, A5, B1, B3, C2, C3, D3, Opaque White, Translucent / Incisal
- **AELITEFLO LV Extra-Flowable (1.5g Syr:** A2, A3.5)

UNIQUE BENEFITS

- Flowable composite can be easily dispensed and injected into the cavity preparation without dripping or slumping
- Low modulus of elasticity microscopically bends with the natural flexure of the tooth to prevent de-bonding and marginal breakdown, especially in Class V restorations
- Uniform particle size is easily polished and stain resistant
- Glass-filled composite resin provides excellent wear resistance and strength

CLINICAL SIGNIFICANCE

- Due to its elastic properties, AELITEFLO "bends with the tooth", leading to better retention. Excellent for abfraction lesions and Class V defects
- AELITEFLO LV exhibits low viscosity and is excellent as the first increment in the gingival floor of a Class II prep, where the material will flow against the matrix band, filling in the irregularities
- Due to its higher filler content, AELITEFLO is an excellent pit and fissure sealant

+ RELATED MATERIALS



**ONE-STEP®
ONE-STEP® PLUS**



BISCO Etchants

Uni-Etch (32%)
Etch37 (37%)
Regular or w/BAC



**BisBlock™
Dentin Desensitizer**