



BISCO Adhesives

All-Bond Universal®
Light-Cured Dental Adhesive

Universal Primer™
Dual-Cured Adhesive

All-Bond SE®
Self-Etching Adhesive

All-Bond 3®
Universal Dental Adhesive System

All-Bond 2®
Universal Dental Adhesive

**One-Step®/
One-Step® Plus**
Universal Dental Adhesive/
Filled Universal Dental Adhesive

Light-Cured



Dual-Cured



Compatible with
Light-, Dual-, and
Self-Cured Materials



Total-Etch Technique



Self-Etch Technique



Selective-Etch Technique



No Additional Activator
Required for Compatibility
with Dual-Cured Products



Single Bottle



Compatible with
Both Direct & Indirect
Restorations



Storage at
Room Temperature



MDP Containing



Best Uses



Consult
Instructions
For Use

A light-cured, single bottle with no additional activator for all direct and indirect restorations to be used with any bonding technique.

Ideal for those who want ONE adhesive in the office for all restorations.

A dual-cure primer/adhesive designed for those who prefer to not light-cure the adhesive layer under indirect restorations.

A light-cured adhesive for all direct and indirect restorations to be used in the self-etch bonding technique.

A dual-cured, hydrophobic adhesive for all direct and indirect restorations to be used in the total-etch bonding technique.

A dual-cure, hydrophilic, BPDM containing monomer adhesive system to be used in the total-etch bonding technique.

A light-cured, single bottle with no additional activator for all direct and indirect restorations to be used in the total-etch bonding technique.

*Use One-Step Plus for a "Filled" adhesive option to reduce post-operative sensitivity.

Exclusively distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200

