



BISCO Cements

TheraCem®

Dual-Cured, Self-Adhesive Calcium & Fluoride Releasing Cement

Duo-Link Universal™

Dual-Cured Resin Luting Cement

Choice™ 2

Light-Cured Veneer Cement

BisCem®

Dual-Cured, Self-Adhesive Resin Cement

Curing Mode

Dual Cured

Dual Cured

Light Cured

Dual Cured

Requires Adhesive

•

•

Substrates

metal, alumina, zirconia, glass ceramics, lithium disilicate, composite

metal, alumina, zirconia, glass ceramics, lithium disilicate, composite

metal, alumina, zirconia, glass ceramics, lithium disilicate, composite

metal, alumina, zirconia, glass ceramics, lithium disilicate, composite

Inlays/Onlays

•

•

•

Crowns/Bridges

•

•

•

Endodontic Posts

•

•

•

Maryland Bridges

•

Veneers

•

Primer required for Metal, Zirconia, Alumina, Composite and Post Restorations

•

•

•

Silane recommended for Glass Ceramics/Porcelain or Lithium Disilicate/Silicate Restorations

•

•

•

•

Best Uses

TheraCem is a self-adhesive cement, ideal for retentive preparations when "bonding" the restoration is NOT needed. Releases both calcium and fluoride. Contains MDP.

Ideal for less retentive preparations when "bonding" the restoration IS needed. Ideal for most restorations including CAD/CAM restorations.

BISCO's recommended cement for veneer restorations. Choice 2 is specifically formulated for color stability.

Dual-Cured, self-adhesive resin cement.

